



EPIGENETICS &
GENE REGULATION



Introduction

Most people grasp genetics as it applies to them as easily as a seahorse grasps things in his surroundings. It begins when a child hears she has her father's eyes and considers how that could be. The father wonders, in turn, if he will get the 'good genes' of his great-grandfather, still alive at 98. In the meantime, Dr. Oz is telling his viewing audience about the links between heart problems and genetics. But where is the discussion, the understanding, the headlines about epigenetics?

Consider this: Gregor Mendel's pea studies were published in the 1860s, when Darwin's *On Origin of the Species* was making the rounds. Conrad Hal Waddington had to understand genetics in order to even suggest an epigenetic landscape, "only" 70 years ago. A query in PubMed using 'epigenetic' as a search term reveals a bounce from 321 publications in 2000 to 1,210 in 2005, jumping again to 4,752 in 2012. A similar search for 'genetics' in PubMed indicates that the same bounce in that field happened in the early 1960s, 50 years ago. Clearly, epigenetics is in its early years.

Cayman Chemical's goal is to help make research possible. In anticipation of a need for better tools for epigenetics research, Cayman has synthesized and purified a collection of research reagents and recombinant proteins to study nucleosomal modification. Furthermore, Cayman has built a dedicated Epigenetic Screening Laboratory and hired scientists to develop assays to profile and screen compounds that modulate epigenetic signaling, described in this catalog. Moreover, the scientists in our Epigenetic Screening Laboratory can provide screening and profiling services according to your needs.

Visit www.caymanchem.com/episcreen for more information.

Abbreviations

ADP	Adenosine Diphosphate	IP	Immunoprecipitation
ChIP	Chromatin Immunoprecipitation	JmJc	Jumonji C
DNA	Deoxyribonucleic Acid	LSD	Lysine-Specific Demethylase
EC₅₀	Half maximal effective concentration	MLL	Mixed-Lineage Leukemia or Myeloid/Lymphoid Leukemia
ELISA	Enzyme-linked Immunosorbent Assay	NAD	Nicotinamide Adenine Dinucleotide
ER	Estrogen Receptor	NF-κB	Nuclear Factor κ-light-chain-enhancer of activated B cells
Gcn5	General control of amino acid synthesis protein	NSD	Nuclear Receptor SET domain-containing protein
GSK	Glycogen Synthase Kinase	PAD	Protein Arginine Deiminase
GST	Glutathione S-Transferase	PHD	Plant Homeodomain
HAT	Histone Acetyltransferase	PRMT	Protein Arginine Methyltransferase
HDAC	Histone Deacetylase	RNA	Ribonucleic Acid
His	Histidine	SCID	Severe Combined Immunodeficiency
IC₅₀	Half maximal inhibitory concentration	SIR	Silent Information Regulator
ICC	Immunocytochemistry	SIRT	Sirtuin
ID₅₀	Infectious dose	WB	Western Blot
IF	Immunofluorescence		
IHC	Immunohistochemistry		

Ordering Information

Orders are accepted by telephone, fax, mail, e-mail, or via the Cayman Chemical website. We will accept telephone orders Monday through Friday from 8 AM to 6 PM EST. All orders received by 1 PM EST will be shipped the same day if stock is available (Monday through Thursday only). Confirming purchase orders must be clearly marked as such to avoid possibility of duplication.

Domestic Shipments

In most instances we ship FedEx Standard Overnight Delivery (not available to all locations), with delivery by 3:30 PM of the next business day. Product availability may vary. Local delivery is available for the Ann Arbor area only. Other shipping options will be considered upon request, but can be granted only under conditions that will ensure the quality of the product. Freight is prepaid and added to the invoice. Please inquire at the time of order for an estimate of the freight charges. If you wish us to ship collect, please supply a valid account number when ordering. Please address all orders to:

Cayman Chemical Company
1180 E. Ellsworth Road
Ann Arbor, MI 48108 USA
Phone: (734) 971-3335

Toll-free Phone: (800) 364-9897
Fax: (734) 971-3640
E-mail: custserv@caymanchem.com
www.caymanchem.com

Include the following information with your order:

- Catalog number, description, size, and quantity desired.
- Complete shipping address. (Delivery is not available to post office box numbers.)
- A complete billing address.
- A purchase order number or major credit card (Visa, MasterCard, or American Express), account number, and expiration date.
- Name of the end user.

Terms

- U.S. funds only, drawn on a U.S. bank.
- Net 30 days.
- F.O.B. Ann Arbor, Michigan, U.S.A.
- Bank fees and wire transfer fees are not to be deducted from the invoice amount.

Returns

Products cannot be returned without prior authorization from Cayman Chemical Company. Please contact our Customer Service Department for return shipping instructions. Custom orders and radioactive material cannot be accepted for return credit if due to a customer's error.

Technical Assistance

Technical assistance is available from 8 AM to 5:30 PM EST. If inquiring about a purchased product, please provide the catalog number, lot number, and date of purchase to our technical staff so they may answer your questions quickly. Technical assistance may be reached toll free at 888-526-5351, via e-mail at techserv@caymanchem.com, or on the web at www.caymanchem.com/techserv.

Use of Research Products

The products in this catalog are not for human or veterinary disease diagnosis or therapeutic drug use. They should be used only by technically qualified individuals or those under their direct supervision. Any individual working directly with these products should have free access to the applicable Material Safety Data Sheet (MSDS) and should read and understand it completely prior to use. Please contact our Customer Service Department or visit the specific product page on our website if you require additional copies of any MSDS.

Patent Disclaimer

The end-user assumes full responsibility for appropriate licensing and/or non-infringement for any proprietary claim or patent.

Warranty and Limitation of Remedy

Cayman Chemical Company makes no warranty or guarantee of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular purpose, non-infringement, suitability, and merchantability, which extends beyond the description of the chemicals hereof. Cayman warrants only to the original customer that the material will meet our specifications at the time of delivery. Cayman will carry out its delivery obligations with due care and skill. Thus, in no event will Cayman have any obligation or liability, whether in tort (including negligence) or in contract, for any direct, indirect, incidental, or consequential damages, even if Cayman is informed about their possible existence. This limitation of liability does not apply in the case of intentional acts or negligence of Cayman, its directors or its employees. Buyer's exclusive remedy and Cayman's sole liability hereunder shall be limited to a refund of the purchase price, or at Cayman's option, the replacement, at no cost to Buyer, of all material that does not meet our specifications. Said refund or replacement is conditioned on Buyer giving written notice to Cayman within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material. Neither party shall be liable to the other in any manner for the failure of or delay in the fulfillment of all or part of their obligations, resulting from causes or circumstances beyond their reasonable control including, but not limited to, floods, fires, hurricanes, tornadoes, earthquakes, other natural calamities and extraordinary weather, insurrections, wars, acts of terrorism, riots, embargoes, governmental refusals to issue approval for export, other governmental orders or restrictions, shortages of shipping vehicles, delays in transportation, inability to obtain supplies and materials, strikes, and lockouts.

Epigenetics & Gene Regulation

4 DNA Methylation Establishing & Editing the Methylome By Olivia L. May, Ph.D.

[Products]
6 Histones & Chromatin Regulators

[Products]
9 Acetyltransferases

[Products]
12 Deacetylases

16 Histone Acetylation, Deacetylation, and Genomic Bistability By Thomas G. Brock, Ph.D.

[Products]
24 Acetyl Readers

26 Bromodomains, ChromoHub, and STRING By Thomas G. Brock, Ph.D.

[Products]
36 Methyl Readers

[Products]
38 Methyltransferases

48 Histone Methylation and the Language of Epigenetics By Thomas G. Brock, Ph.D.

[Products]
53 Demethylases

58 Histone Demethylation By Daniel A. Bochar, Ph.D.

[Information]
61 Article References

[Information]
62 Indices

DNA Methylation

establishing & editing the methylome

by [Olivia L. May, Ph.D.]

Throughout our lifespan, chemical changes subtly occur in our DNA. A comparison of the DNA of a newborn baby with that of a centenarian shows that the scope of these changes can be dramatic, potentially having an influence on disease.¹ This may help explain why our risk of cancer and other diseases increases as we get older. Such exogenous influences can also be inherited, impacting the genetics of an individual's offspring. Epigenetic regulation involves genetic control by factors other than an individual's DNA sequence. These heritable changes alter DNA accessibility and chromatin structure, thereby regulating patterns of gene expression. Through epigenetic programming, genes are switched on or off as needed to determine which proteins are transcribed and expressed. This process is crucial to normal development (controlling genomic imprinting and X-chromosome inactivation) as well as for the suppression of transposable elements and the maintenance of stable cellular identities. Unfortunately, changes in DNA methylation patterns also contribute to human diseases, like cancer, for which risk increases with age.

DNA methylation, the chemical process that adds a methyl group to DNA, is the principal manner in which genes are transcriptionally repressed or silenced. This modification occurs specifically at CpG sites, regions in which a cytosine nucleotide is located next to a guanine nucleotide that is linked by




Figure 1. DNA methylation and demethylation. *De novo* DNA methylation at the cytosine in CpG dinucleotides is initiated by DNMT3A and DNMT3B. After replication, DNMT1 maintains the methylation state in the daughter strands. Base excision repair mechanisms facilitate removal of methylation after deamination of methyl cytosine (5mC) creates a T-G mismatch.

a phosphate (Figure 1). Inserting methyl groups changes the appearance and structure of DNA, which may directly block DNA recognition and binding of transcription factors, or may attract other factors that preferentially bind to DNA to interfere with transcription factor accessibility. Three families of proteins which bind methylated DNA have been identified so far.^{2,3} These include MBD domain proteins, Kaiso and Kaiso-like proteins, and SRA domain proteins. By recruiting these proteins, DNA methylation marks can promote the persistence of certain histone states, such as deacetylation, thus enabling posttranslational histone modifications. As an example, methyl CpG-binding domain protein 2 (MeCP2), a member of the MBD family, binds to methyl CpG and recruits HDACs, which promote chromosome condensation and transcriptional repression.

Establishing methylation marks

CpG sites are methylated by one of three DNA methyltransferases (DNMTs). During embryogenesis, *de novo* methylation is typically performed by DNMT3A and DNMT3B. Both have similar structures, comprising a PWWP domain, PHD-like or ADD domain, and a carboxy-terminal catalytic domain. The PWWP domains are necessary for binding of DNMT3A and DNMT3B to chromatin *in vivo*. DNMT3L, which lacks the PWWP and a functional catalytic domain, forms a heterotetramer with DNMT3A or DNMT3B and stimulates the activity of its *de novo* partners, guiding the recognition of DNA

targets.⁴ The PHD domain of DNMT3L interacts with the amino terminal tail of histone H3 and is activated only when bound to unmethylated DNA.⁵

Appropriate combinations of histone modifications create either protective or permissive conditions for the docking of *de novo* methylation complexes. Specific methylation of the lysine at residue 4 (H3K4) by SETD1 inhibits DNMT3L binding. H3K4 methylation is a marker for active genes, and there is an inverse correlation between the presence of trimethylated H3K4 and DNA methylation at promoters.^{5,6} Indeed, chromatin at the majority of unmethylated CpGs is enriched in H3K4 di- and tri-methylation. Removal of H3K4 methylation by the histone demethylase LSD2 (KDM1B) is required for establishment of DNA methylation at the promoters of certain imprinted genes. Additionally, DNMT3A binding is promoted by SETD2 trimethylation of lysine 36 of histone 3 (H3K36me3), whereas the H3K36me2 demethylase KDM2A binds to unmethylated CpGs resulting in depletion of H3K36me2.⁶ Many questions still remain however as to where and how patterns of methylation are established to target appropriate, region-specific histone modifications.⁷

Maintaining methylation marks for long term stability

The DNA methylation patterns established during embryonic development are faithfully copied through somatic cell divisions in order to maintain a gene's transcriptionally active or inactive state. The ubiquitously expressed DNMT1 is predominantly responsible for maintaining cellular levels of CpG methylation. With guidance from the UHRF1 domain it recognizes hemimethylated DNA and methylates appropriate cytosines in newly synthesized daughter strands formed during replication (Figure 1). The base pairing of CpG allows for the reciprocal preservation of methylation during subsequent replication cycles. CpG islands keep their overall unmethylated state (or possibly methylated state) extremely stably through multiple cell generations, and DNMT1 is partly responsible for this stability.


In early embryogenesis, following sex-determination of the embryo, methylation is erased throughout the genome (in order to reset germ-cell specification) and then reestablished in all but CpG islands.⁸ This guarantees renewal of totipotency at each generation allowing new methylation marks to be established.⁹ As developmental differentiation proceeds, these marks accumulate especially in promoter or other gene regulatory regions to repress transcription of certain key pluripotency genes. This cycle of early embryonic demethylation followed by *de novo* methylation is critical in determining somatic DNA methylation patterns. Once established, somatic DNA methylation (a nongenetic trait) is passed through daughter cell generations, and with it, the contextual effects on gene expression. Methylation is considered a long-term, relatively stable epigenetic trait that contributes to maintenance of the cellular phenotype. While a significant fraction of CpG islands are prone to progressive methylation in certain tissues to maintain permanent cell lines, *de novo* methylation also occurs in mature somatic cells, especially in abnormal cells such as cancers.

Fine-tuning methylation marks


The removal of CpG methylation predictably occurs in cellular reprogramming during gametogenesis, after zygote fertilization, or to preserve induced pluripotency. This demethylation process requires the action of both cytidine deaminases and DNA repair mechanisms (Figure 1). Enzymatic deamination of 5-methylcytosine (5mC) leads to formation of thymine and T:G base-pair mismatches. Base excision repair mechanisms subsequently delete thymine and restore C:G base pairing during epigenetic reprogramming. Spontaneous deamination of 5mC also requires base excision repair mechanisms to repair base-pairing mismatch. This process is highly inefficient in most differentiated cells, however, as spontaneous deamination of 5mC typically results in an overall depletion of CpG dinucleotide sequences. Until very recently, CpG methylation was considered long-lasting and difficult to eliminate post-differentiation, save for inactivation of DNMTs

“With the recent discovery that cytosines can be hydroxymethylated to 5-hydroxymethylcytosine (5hmC), an active mechanism for demethylation has since been eagerly endorsed.”

Passive DNA Demethylation




DNMTs




deamination


Active DNA Demethylation




TDG/SMUG1



DNMTs



TET



AID

Figure 2. Potential mechanism of passive and active DNA demethylation. Passive DNA demethylation was thought to occur by a reduction in activity or absence of DNMTs. Active DNA demethylation involves 3 different enzyme families to enable DNA repair: (1) 5-methylcytosine (5mC) is hydroxylated by TET to form 5-hydroxymethylcytosine (5hmC). (2) 5mC (or 5hmC) is deaminated by AID to form 5-methyluracil (5mU) or 5-hydroxymethyluracil (5hmU). (3) TDG and SMUG1 base excision repair (BER) glycosylases replace these intermediates (5hmU), culminating in cytosine replacement and DNA demethylation.

or after spontaneous deamination and mismatch repair. DNA methylation marks were generally thought to be lost ‘passively’ by lack of DNA methylation maintenance at replication (Figure 1), resulting in progressive loss of methylation at each cell division.¹⁰

With the recent discovery that cytosines can be hydroxymethylated to 5-hydroxymethylcytosine (5hmC), an active mechanism for demethylation has since been eagerly endorsed (Figure 2).¹¹ 5hmC, which impairs remethylation by DNMTs until it is replaced by DNA repair, has been hypothesized to serve as an intermediate in the removal of methylated cytosines. Three different enzyme families are thought to drive active demethylation: the ten-eleven translocation (TET) family, the AID (Activation-Induced cytidine Deaminase) family, and a family of base excision repair (BER) glycosylases (Figure 2). TET has been shown to hydroxylate and then further oxidize methylated cytosines to catalyze their conversion to 5hmC. Targeting specific hydroxymethylated loci, AID

deaminates 5hmC to uracil (5hmU), which is eventually replaced with an unmethylated cytosine by thymine-DNA glycosylase (TDG) and single-strand-selective monofunctional uracil-DNA glycosylase 1 (SMUG1).

Thus, DNA methylation is proving to be a dynamic process, requiring continuous regulation and potentially having an important editing role for cellular signaling or tissue-specific differentiation. Further understanding of the nuances of addition and removal of methylation marks on DNA will continue to inspire therapeutic strategies for targeting DNA methylation in the prevention of cancer and other human diseases.¹²⁻¹⁴ In contribution to this effort, Cayman offers several potent DNMT inhibitors including the nucleoside analogs: 5-Azacytidine (Item No. 11164), Decitabine (Item No. 11166), and Zebularine (Item No. 10975) as well as the less cytotoxic, non-nucleoside analog: RG-108 (Item No. 13302).

Histones & Chromatin Regulators

Antibodies

- 6 13503 Chromosome Associated Protein-C Polyclonal Antibody (aa 47-61)
- 6 13501 Chromosome Associated Protein-C Polyclonal Antibody (aa 281-297)
- 7 13535 Histone H2A Polyclonal Antibody
- 7 13538 Histone H2B (C-Term) Polyclonal Antibody
- 7 13539 Histone H2B (N-Term) Polyclonal Antibody
- 7 13540 Histone H3 (Phospho-Ser²⁸) Monoclonal Antibody (Clone 117C826)
- 8 13784 Histone H3.3 Polyclonal Antibody
- 8 13543 Histone H4 Polyclonal Antibody
- 8 13776 Mammalian STE-20-Like Kinase 1 Polyclonal Antibody

Proteins

- 6 11010 Core Histones (human)
- 6 10261 Histone H2A (Xenopus recombinant)
- 7 10262 Histone H2B (Xenopus recombinant)
- 7 10263 Histone H3 (human recombinant)
- 8 10877 Histone H3 Peptide Substrate (1-21)
- 8 10530 Histone H3 Trimethyl Lys9 Peptide
- 8 10264 Histone H4 (human recombinant)
- 8 10854 Histone H4 Peptide Substrate (1-21)
- 8 10380 Histone H4 Peptide Substrate (15-24)


Chromosome Associated Protein-C Polyclonal Antibody (aa 47-61)

13503

CAP-C

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human CAP-C amino acids 47-61 • Host: rabbit • Cross Reactivity: (+) human CAP-C • Application(s): WB • CAP-C plays a critical role in the structural maintenance of chromosomes, including proper condensation and segregation.




Chromosome Associated Protein-C Polyclonal Antibody (aa 281-297)

13501

CAP-C

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human CAP-C amino acids 281-297 • Host: rabbit • Cross Reactivity: (+) human CAP-C • Application(s): WB • CAP-C plays a critical role in the structural maintenance of chromosomes, including proper condensation and segregation.




Core Histones (human)

11010

HeLa Core Histones

Purity: ≥95% **Stability:** ≥6 months at -80°C

Source: Highly purified mixture of human core histones (H2A, H2B, H3, and H4) isolated *via* hydroxyapatite chromatography from HeLa S3 (human cervical adenocarcinoma) nuclear pellet • A histone octamer consists of two copies of each of the core histones, H2A, H2B, H3, and H4 which dimerize to create the nucleosome core.




Histone H2A (Xenopus recombinant)

10261

M_r: 13.9 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C

Lyophilized powder Source: Recombinant protein consisting of amino acids 1-129 expressed in *E. coli* • Histone 2A is one of the four histones that comprise a nucleosome protein core.




Histone H2A Polyclonal Antibody

13535

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human histone H2A amino acids 1-15 and 81-96 • Host: rabbit • Cross Reactivity: (+) human and mouse histone H2A • Application(s): ELISA and WB • Histone H2A is one of the five main histone proteins involved in the structure of chromatin in eukaryotic cells. It is considered a core histone and forms a dimer with H2B; the core molecule is complete when H3-H4 also attaches to form a tetramer.




Histone H2B (Xenopus recombinant)

10262

M_r: 13.7 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C

Lyophilized Powder Source: Recombinant protein consisting of amino acids 1-123 expressed in *E. coli* • Histone H2B is one of the core nucleosomal histones. It undergoes many modifications which include acetylation, methylation, and phosphorylation that are important for regulation of gene transcription.




Histone H2B (C-Term) Polyclonal Antibody

13538

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human histone H2B amino acids 111-125 • Host: rabbit • Cross Reactivity: (+) chicken, canine, *Drosophila*, human, mouse, rat, most mammals, and zebrafish histone H2B • Application(s): WB • Histone H2B is one of the 5 main histone proteins involved in the structure of chromatin in eukaryotic cells. It features a main globular domain and a long N-terminal tail and forms a dimer with H2A to compose part of the nucleosome core.




Histone H2B (N-Term) Polyclonal Antibody

13539

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human histone H2B • Host: rabbit • Cross Reactivity: (+) human histone H2B • Application(s): WB • Histone H2B is one of the 5 main histone proteins involved in the structure of chromatin in eukaryotic cells. It features a main globular domain and a long N-terminal tail and forms a dimer with H2A to compose part of the nucleosome core.




Histone H3 (human recombinant)

10263

M_r: 15.5 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C

Lyophilized powder Source: Recombinant protein consisting of amino acids 1-136 expressed in *E. coli* • Histone H3 is one of the core nucleosomal histones. It undergoes many modifications which include acetylation, methylation, and phosphorylation that are important for regulation of gene transcription.




Histone H3 (Phospho-Ser²⁸) Monoclonal Antibody (Clone 117C826)

13540

PHH3

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human histone H3 • Host: mouse, clone 117C826 • Cross Reactivity: (+) human histone H3 • Application(s): WB • H3 phosphorylation at serine 28 is coupled with mitotic chromosome condensation in diverse mammalian cell lines.



Histone H3 Peptide Substrate (1-21) 10877*H3 Peptide***FW:** 2,280.7 **Purity:** ≥95% by HPLCA lyophilized peptide **Stability:** ≥1 year at -20°C**Summary:** A target substrate for several of the histone modifying enzymes including lysine methyltransferases, arginine methyltransferases, acetyltransferases, and others.


1 mg

Histone H3 Trimethyl Lys9 Peptide 10530*H3K9me3 Peptide***FW:** 1,601 **Peptide Sequence:** ARTKQTARK(Me)₃-STGGKAA lyophilized peptide **Stability:** ≥1 year at -20°C**Summary:** Peptide contains a trimethylated lysine residue at position nine which is a substrate for several of the JmjC domain-containing class of histone demethylases; JMJD2A and JMJD2D exhibit higher affinity with H3K9me3 compared to H3K9me2


1 mg

Histone H3.3 Polyclonal Antibody 13784*H3.3A, H3.3B, H3F3A, H3F3B*Protein A-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: human histone H3.3 amino acids 100-136 • Host: rabbit • Cross Reactivity: (+) chicken, ovine, *Drosophila*, equine, human, mouse, and opossum histone H3.3 • Application(s): IHC and WB • Histone H3.3 constitutes the predominant form of histone H3 in non-dividing cells and is incorporated into chromatin independently of DNA synthesis. It has a central role in transcription regulation, DNA repair, DNA replication, and chromosomal stability.

1 ea



Lane 1: 293 cell lysate (Absence of immunizing peptide)
Lane 2: 293 cell lysate (Presence of immunizing peptide)
Lane 3: NIH 3T3 cell lysate

Histone H4 (human recombinant) 10264**M_r:** 11.5 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°CLyophilized powder **Source:** Recombinant protein consisting of amino acids 1-103 expressed in *E. coli* • Histone H4 is one of the core nucleosomal histones. The N-terminal tail of histone H4 undergoes many modifications which include acetylation, methylation, and phosphorylation that are important for regulation of gene transcription.50 µg
100 µg
250 µg
1 mg

Lane 1: MW Markers
Lane 2: Purified H4 (0.5 µg)
Lane 3: Purified H4 (1 µg)
Lane 4: Purified H4 (2 µg)

Histone H4 Peptide Substrate (1-21) 10854*H4 Peptide***FW:** 2,091.5 **Purity:** ≥95% by HPLCA lyophilized peptide **Stability:** ≥1 year at -20°C**Summary:** A target substrate for several of the histone modifying enzymes including lysine methyltransferases, arginine methyltransferases, acetyltransferases, and others.

1 mg


Histone H4 Peptide Substrate (15-24) 10380*Histone H4K20 Peptide, SET8 Methyltransferase Acceptor Peptide***FW:** 1,278 **Peptide Sequence:** AKRHRKVLRLD-NH₂Peptide lyophilized from ammonium bicarbonate buffer **Stability:** ≥1 year at -20°C**Summary:** Peptide contains a lysine at position 20 which is a substrate or acceptor peptide for the lysine methyltransferases KMT5A (SET8) and KMT5B (SUV4-20H1)

1 mg

5 mg

Histone H4 Polyclonal Antibody 13543Protein G-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: human histone H4 amino acids 15-30 • Host: rabbit • Cross Reactivity: (+) human histone H4 • Application(s): WB • Histone H4 is a structural component of the nucleosome and is subject to covalent modification including acetylation and methylation, which may alter expression of genes located on DNA associated with its parent histone octamer.


1 ea



Lane 1: Histone H4

Mammalian STE-20-Like Kinase 1 Polyclonal Antibody 13776*KRS2, MST-1, STK4*Protein G-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: human MST-1 amino acids 372-390 • Host: rabbit • Cross Reactivity: (+) human MST-1 • Application(s): WB • MST-1 is a serine/threonine kinase that has been implicated in the promotion of chromatin condensation.

1 ea



Lane 1: Jurkat cell lysate

Acetyltransferases

Antibodies

- 9 10010567 **Acetyl Lysine Monoclonal Antibody (Clone 7F8)**
9 13725 **Acetyl Lysine Polyclonal Antibody-biotin**
9 13726 **Acetyl Lysine Polyclonal Antibody HRP Conjugate**
11 13789 **TIP60 Polyclonal Antibody**

Biochemicals

- 9 13144 **Anacardic Acid**
9 12095 **Butyrolactone 3**
9 10547 **4-pentynoyl-Coenzyme A (trifluoroacetate salt)**
9 12086 **CPH2**
10 11012 **Delphinidin chloride**
10 10566 **Garcinol**

Kits

- 10 10006515 **HAT Inhibitor Screening Assay Kit**

Proteins

- 10 10782 **Gcn5 (human recombinant)**
10 10009115 **pCAF Histone Acetyltransferase**
10 10783 **TIP60 (human recombinant)**

Acetyl Lysine Monoclonal Antibody (Clone 7F8) 10010567Purified IgG₁ lyophilized **Stability:** ≥1 year at -20°C**Summary:** Antigen: Acetylated KLH • Host: mouse • Isotype: IgG₁ • Cross Reactivity: (+) acetylated lysine residues; (-) non-acetylated lysine residues • Application(s): ELISA, ICC, and WB • This antibody is useful for monitoring levels of acetylation on various proteins (e.g., histones and p53).





1 ea

Acetyl Lysine Polyclonal Antibody-biotin 13725Rabbit immunoglobulin in PBS **Stability:** ≥1 year at -20°C**Summary:** Antigen: acetylated KLH • Host: rabbit • Cross Reactivity: (+) acetylated lysine residues; (-) non-acetylated proteins • Application(s): ELISA, IF, IP, and WB • This biotin-tagged antibody is useful for monitoring levels of acetylation on various proteins.

400 µl

Acetyl Lysine Polyclonal Antibody HRP Conjugate 13726Rabbit immunoglobulin in PBS **Stability:** ≥1 year at -20°C**Summary:** Antigen: acetylated KLH • Host: rabbit • Cross Reactivity: (+) multi-species • Application(s): ELISA, IF, IP, and WB • This HRP-conjugated antibody is useful for monitoring levels of acetylation on various proteins.


400 µl

Anacardic Acid 13144*[16611-84-0] 6-pentadecyl Salicylic Acid***MF:** C₂₂H₃₆O₃ **FW:** 348.5 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** An alkyl salicylic acid isolated from cashew shells; inhibits the HAT activity of p300 and pCAF (IC₅₀ = 8.5 and 5 µM, respectively); suppresses NF-κB activation, inhibits IκB-α phosphorylation, and prohibits p65 nuclear translocation1 mg
5 mg
10 mg
25 mg**Butyrolactone 3** 12095*[778649-18-6]***MF:** C₉H₁₂O₄ **FW:** 184.2 **Purity:** ≥95%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** Specifically inhibits the histone acetyltransferase Gcn5 (IC₅₀ = 100 µM) and can inhibit pre-RNA splicing with an IC₅₀ value of 0.5 mM1 mg
5 mg**4-pentynoyl-Coenzyme A (trifluoroacetate salt)** 10547*Click Tag™ 4-pentynoyl-CoA***MF:** C₂₆H₄₀N₇O₁₇P₃S • CF₃COOH **FW:** 961.6 **Purity:** ≥95%A lyophilized powder **Stability:** ≥2 years at -20°C**Summary:** An acyl-CoA donor that can be metabolically transferred onto lysine residues of proteins by lysine acetyltransferases; an azide-alkyne bioconjugation reaction, known as 'click chemistry', can then be used to tag the acetylated proteins with fluorescent or biotinylated labels for subsequent analysis500 µg
1 mg
5 mg
10 mg**CPH2** 12086*[357649-93-5]***MF:** C₁₄H₁₄ClN₃S **FW:** 291.8 **Purity:** ≥95%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** Specifically inhibits Gcn5-dependent acetylation of histone H3K14 at a concentration of 0.8 mM both *in vitro* and *in vivo*5 mg
10 mg

Delphinidin chloride

11012


[528-53-0] Ephedrine

MF: C₁₅H₁₁ClO₇ **FW:** 338.7 **Purity:** ≥98%
A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A natural plant pigment which induces the release of nitric oxide by vascular endothelium, causing vasorelaxation; inhibits signaling through EGFRs, suppressing the expression of ERα and inducing both apoptosis and autophagy at a dose of 1-40 μM; inhibits the HAT activities of p300/CBP (IC₅₀ = ~30 μM)1 mg
5 mg
10 mg

Garcinol


10566

[78824-30-3] Camboginol

MF: C₃₈H₅₀O₆ **FW:** 502.8 **Purity:** ≥95%
A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** An inhibitor of the HATs p300 and pCAF (IC₅₀ = 7 and 5 μM, respectively) that also inhibits the HAT Gcn5; promotes neurogenesis and *ex vivo* expansion of human hematopoietic stem cells; induces apoptosis in several types of cancer cells and has anti-inflammatory actions1 mg
5 mg
10 g
25 g

Gcn5 (human recombinant)

10782

*General control of amino acid synthesis protein 5-like 2, KAT2A, Lysine acetyltransferase 2A, STAF97***M:** 96.3 kDa **Purity:** ≥80% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal His-tagged protein consisting of amino acids 2-837 expressed in Sf21 cells • Recombinant Gcn5 preferentially acetylates lysine 14 on histone H3 *in vitro*; however, alone it is unable to acetylate nucleosomal core histone substrates. Acetylation of the nucleosomal histones requires that Gcn5 be a part of either the multisubunit SAGA or ATAC protein complexes, which have a broad substrate specificity, including H3K9, H3K18, H4K8, and H4K16, as well as additional sites on histone H2B.25 μg
50 μg
100 μgLane 1: MW Markers
Lane 2: Gcn5 (1 μg)
Lane 3: Gcn5 (2 μg)
Lane 4: Gcn5 (4 μg)

HAT Inhibitor Screening Assay Kit

10006515

*Histone Acetyltransferase***Stability:** ≥1 year at -20°C**Summary:** Cayman's HAT Inhibitor Screening Assay Kit provides a fast, fluorescence-based method for evaluating pCAF HAT inhibitors. The procedure requires only three easy steps, all performed in the same microwell plate. In the first step of the protocol, HAT is incubated with acetyl-CoA and the histone H3 peptide. During this time, HAT catalyzes the enzymatic transfer of acetyl groups from acetyl-CoA to the H3 peptide producing an acetylated peptide and CoASH. Following addition of isopropanol to stop the enzymatic reaction, CPM is added to the wells of the plate. CPM reacts with the free thiol groups present on CoASH forming a highly fluorescent product that is detected using excitation and emission wavelengths of 360-390 and 450-470 nm, respectively.

96 wells


pCAF Histone Acetyltransferase

10009115

*HAT, p300/(CREB binding protein) Associated Factor***M:** ~40 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C**Source:** Active recombinant GST-tagged protein consisting of amino acids 2-837 expressed in *E. coli* • pCAF belongs to the GCN5/pCAF family of nuclear HATs. Cayman's pCAF preparation contains 165 amino acids from the HAT activity domain of human pCAF fused to GST at the N-terminus. Enzyme activity was determined using a fluorescent HAT assay and is comparable to that found in the literature.25 μg
50 μg
100 μg

TIP60 (human recombinant)

10783


*cPLA(2) Interacting Protein, Esa1, HIV-1 Tat Interacting Protein 60 kDa, Hs.6364, HTATIP, KAT5, PLIP***M:** 60.3 kDa **Purity:** ≥80% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal His-tagged protein consisting of amino acids 2-513 expressed in Sf21 cells • TIP60 is a member of the MYST family of lysine acetyl transferases. It has been shown to acetylate histones, p53, and the Ataxia Telangiectasia Mutant protein kinase.25 μg
50 μg
100 μgLane 1: TIP60 (5 μg)
Lane 2: TIP60 (10 μg)
Lane 3: TIP60 (20 μg)
Lane 4: MW Markers

TIP60 Polyclonal Antibody

13789

*KAT5, Lysine Acetyltransferase 5*Protein A-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: peptide within the region of human TIP60 amino acids 480-530 • Host: rabbit • Cross Reactivity: (+) human (isoform CRA_b), chimpanzee, orangutan, equine, canine, mouse, ovine, rat, opossum, zebrafish, and *Xenopus* TIP60 • Application(s): IHC and WB • TIP60 belongs to the MYST family of HATs. It is a catalytic subunit of the NuA4 HAT complex which is involved in transcriptional activation of select genes principally by acetylation of nucleosomal histones H4 and H2A.

1 ea

Lane 1: Human brain lysate (Absence of immunizing peptide)
Lane 2: Human brain lysate (Presence of immunizing peptide)
Lane 3: Mouse brain lysate
Lane 4: Rat brain lysate

Deacetylases

Antibodies

15	13493	HDAC3 Polyclonal Antibody
15	13494	HDAC4 Polyclonal Antibody
15	13499	HDAC6 Polyclonal Antibody
15	13500	HDAC7 (Phospho-Ser ¹⁵⁵) Polyclonal Antibody
18	13504	HDAC11 Polyclonal Antibody
18	13778	Metastasis Associated 1 Family Member 2 Polyclonal Antibody
19	13785	p66α Polyclonal Antibody
22	13477	SIRT7 Polyclonal Antibody

Biochemicals

12	13145	AGK2
12	14004	AK-7
12	10575	Apicidin
13	89740	CAY10398
13	10005019	CAY10433
13	10009797	CAY10591
13	13146	CAY10603
13	13172	CBHA
13	13686	Chidamide
13	12084	CI-994
13	10009798	EX-527
13	10576	HC Toxin
13	13277	(S)-HDAC-42
18	13295	HNHA
18	10641	JGB1741
18	13174	M 344
18	13284	MS-275
19	13176	Oxamflatin
19	10444	PCI 34051
19	13212	Pimelic Diphenylamide 106
19	13870	Pyroxamide
19	10009929	SAHA
19	10675	SAHA-BPyne
19	10671	coumarin-SAHA
20	10495	4-iodo-SAHA
20	13178	Salermide
20	10443	SB 939
20	10572	Scriptaid
23	10523	Sirtinol
23	13121	Sodium Butyrate
23	13168	Splitomicin
23	10574	Suberohydroxamic Acid
23	13085	Tenovin-1
23	13086	Tenovin-6
23	89730	Trichostatin A
23	10559	Tubastatin A (trifluoroacetate salt)
23	13033	Valproic Acid (sodium salt)

Kits

14	10011563	HDAC Activity Assay Kit
14	600150	HDAC Cell-Based Activity Assay Kit
14	10011564	HDAC1 Inhibitor Screening Assay Kit
18	700230	HDAC8 Inhibitor Screening Assay Kit
20	10010401	SIRT1 Direct Fluorescent Screening Assay Kit
21	10010991	SIRT1 FRET-Based Screening Assay Kit
20	700280	SIRT2 Direct Fluorescent Screening Assay Kit
21	10011566	SIRT3 Direct Fluorescent Screening Assay Kit
21	700290	SIRT6 Direct Fluorescent Screening Assay Kit


Proteins

14	10009231	HDAC1 (human recombinant)
14	10009377	HDAC2 (human recombinant)
14	10009232	HDAC3/NCOR2 (human recombinant)
15	10009652	HDAC4 (human recombinant)
15	10009379	HDAC5 (human recombinant)
15	10009465	HDAC6 (human recombinant)
15	19380	HDAC8 (human recombinant)
18	10009466	HDAC9 (human recombinant)
19	11633	NCOR2/SMRT (human recombinant)
21	10011190	SIRT1 (human recombinant)
21	10011191	SIRT2 (human recombinant)
22	10011194	SIRT3 (human recombinant)
22	10317	SIRT4 (human recombinant)
22	10318	SIRT5 (human recombinant)
22	10315	SIRT6 (human recombinant)
22	10316	SIRT7 (human recombinant)

AGK2

13145


[304896-28-4]

MF: C₂₃H₁₃Cl₂N₃O₂ FW: 434.3 Purity: ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A cell-permeable, selective inhibitor of SIRT2 (IC₅₀ = 3.5 μM) that minimally affects either SIRT1 or SIRT3; rescues dopamine neurons from α-synuclein toxicity in both *in vitro* and *in vivo* Parkinson's disease models1 mg
5 mg
10 mg
25 mg

AK-7

14004


[420831-40-9]

MF: C₁₉H₂₁BrN₂O₃S FW: 437.4 Purity: ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A cell- and brain-permeable inhibitor of SIRT2 (IC₅₀ = 15.5 μM); diminishes neuronal cell death induced by mutant huntingtin fragment in culture; down-regulates cholesterol biosynthetic gene expression and reduces total cholesterol levels in neurons *in vivo*5 mg
25 mg

Apicidin

10575


[183506-66-3] OSI 2040

MF: C₃₄H₄₉N₅O₆ FW: 623.8 Purity: ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A fungal toxin that demonstrates selective inhibition of HDAC3/NCOR over HDAC6 (IC₅₀s = 15.8 and 665.1 nM, respectively); has broad spectrum activity against Apicomplexan parasites and exhibits antiproliferative activity against various cancer cell lines (IC₅₀s = 0.13-2.36 μM)1 mg
5 mg
10 mg

CAY10398

89740


[193551-00-7] MD 85, PX 089274

MF: C₁₅H₂₃N₃O₃ FW: 293.4 Purity: ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** An inhibitor of HDAC (IC₅₀ = 10 μM)1 mg
5 mg
10 mg
25 mg

CAY10433

10005019


[537034-17-6] BML-210, N-phenyl-N'-(2-Aminophenyl)hexamethylenediamide

MF: C₂₀H₂₅N₃O₂ FW: 339.4 Purity: ≥95%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** An HDAC inhibitor with an IC₅₀ value of 30 μM when tested in HeLa cell nuclear extracts using 200 μM acetylated fluorometric substrate1 mg
5 mg
10 mg
25 mg

CAY10591

10009797


[839699-72-8] SIRT1 Activator 3, Sirtuin 1 Activator 3

MF: C₂₀H₂₅N₃O₂ FW: 367.5 Purity: ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** An activator of SIRT1 that decreases TNF-α levels from 325 pg/ml (control) to 104 and 53 pg/ml at 20 and 60 μM, respectively; exhibits a significant dose-dependent effect on fat mobilization in differentiated adipocytes1 mg
5 mg
10 mg
25 mg

CAY10603

13146


[1045792-66-2]

MF: C₂₂H₃₀N₄O₆ FW: 446.5 Purity: ≥95%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A potent and selective inhibitor of HDAC6 (IC₅₀ = 0.002 nM, as compared with 271, 252, 0.42, 6851, and 90.7 nM for HDAC1, 2, 3, 8, and 10, respectively); prevents the growth of several pancreatic cancer cell lines (IC₅₀ = 0.1-1 μM)500 μg
1 mg
5 mg
10 mg

CBHA

13172


[174664-65-4] m-Carboxycinnamic Acid bis-Hydroxamine, HDAC Inhibitor II

MF: C₁₀H₁₀N₂O₄ FW: 222.2 Purity: ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** HDAC1 and HDAC3 inhibitor (ID₅₀ = 0.01 and 0.07 μM, respectively, *in vitro*); induces apoptosis in nine different neuroblastoma cell lines in culture (0.5-4.0 μM) and completely suppresses neuroblastoma tumor growth in SCID mice at 200 mg/kg5 mg
10 mg
25 mg
50 mg

Chidamide

13686


[743420-02-2]

MF: C₂₂H₁₉FN₄O₂ FW: 390.4 Purity: ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** An HDAC inhibitor that increases histone H3 acetylation levels in LoVo and HT29 colon cancer cells at concentrations as low as 4 μM; dose-dependently decreases the activation of several oncogenic signaling kinases and induces cell cycle arrest in colon cancer cells1 mg
5 mg
10 mg
25 mg

CI-994

12084


[112522-64-2] N-Acetyldinaline, Goe 5549, PD 123654, Tacedinaline

MF: C₁₃H₁₅N₃O₂ FW: 269.3 Purity: ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** An inhibitor of class I HDACs with IC₅₀ values of 0.9, 0.9, 1.2, and >20 μM for recombinant human HDAC 1, 2, 3, and 8, respectively; displays a wide spectrum of antitumor activity, particularly in tumors normally refractory to conventional anticancer agents5 mg
10 mg
50 mg

EX-527

10009798


[49843-98-3]

MF: C₁₃H₁₃ClN₂O FW: 248.7 Purity: ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A cell-permeable, selective inhibitor of SIRT1 (IC₅₀ = 98 nM); inhibits other SIRTs only at much higher concentrations and has no effect on other HDACs1 mg
5 mg
10 mg
25 mg

HC Toxin

10576


[83209-65-8] Toxin I (Helminthosporium carbonum)

MF: C₂₂H₃₄N₄O₆ FW: 450.5 Purity: ≥95%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A cell-permeable, reversible inhibitor of HDACs (IC₅₀ = 30 nM)500 μg
1 mg

(S)-HDAC-42

13277

[935881-37-1] AR42

MF: C₁₈H₂₀N₂O₃ FW: 312.4 Purity: ≥95%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A potent inhibitor of HDACs (IC₅₀ = 16 nM *in vitro*); decreases the viability of prostate cancer cell lines (IC₅₀ = 0.40 μM); strongly suppresses the growth of PC-3 tumor xenografts1 mg
5 mg
10 mg
25 mg


HDAC Activity Assay Kit

10011563

Stability: ≥1 year at -80°C

Summary: Cayman's HDAC Activity Assay Kit provides a fast, fluorescence-based method for measuring Class I and II HDAC activity that eliminates radioactivity, extraction, or chromatography. The procedure requires only two easy steps, both performed in the same microplate. The fluorescent reaction product is analyzed using a plate reader with excitation wavelengths of 340-360 nm and emission wavelengths of 440-465 nm.

96 wells




HDAC Cell-Based Activity Assay Kit

600150

Stability: ≥6 months at -80°C

Summary: Cayman's HDAC Cell-Based Assay Kit provides an easy tool for studying HDAC activity modulators in whole cells. By using a cell-permeable HDAC substrate, the activity of various protein lysine-specific deacetylases including HDAC1-containing complexes can be measured in intact cells in a simple and homogenous manner. The fluorescence of the deacetylated reaction product can be analyzed using a plate reader or a fluorometer with excitation wavelengths of 340-360 nm and emission wavelengths of 440-465 nm. An HDAC inhibitor, trichostatin A, is included for checking specificity of the HDAC reaction. This assay parallels Cayman's HDAC Activity Assay Kit (Item No. 10011563), which uses a nuclear extract rather than whole cells for the assay. Together, both assays will help to identify whether an inhibitor/activator has a direct effect on the enzyme.

96 wells



HDAC1 (human recombinant)

10009231

M_r: ~79.9 kDa **Purity:** >10% by SDS-PAGE **Stability:** ≥6 months at -80°C

Source: Active recombinant protein containing a C-terminal GST-tag expressed in Sf21 cells • HDAC1 is a Class I HDAC that catalyzes the deacetylation of core histones and other non-histone proteins to control complex biological events, including cell development, differentiation, programmed cell death, angiogenesis, and inflammation.

1 ea


HDAC1 Inhibitor Screening Assay Kit

10011564

Stability: ≥1 year at -80°C

Summary: Cayman's HDAC1 Inhibitor Screening Assay Kit provides a fast, fluorescence-based method for screening HDAC1 inhibitors. The procedure requires only two easy steps, both performed in the same microplate. The fluorescent reaction product is analyzed using a fluorometer with excitation wavelengths of 340-360 nm and emission wavelengths of 440-465 nm. Sufficient purified HDAC1 is provided for 96 assays.

96 wells



HDAC2 (human recombinant)

10009377

M_r: ~60 kDa **Purity:** ≥70% by SDS-PAGE **Stability:** ≥6 months at -80°C

Source: Active full-length recombinant protein containing a C-terminal His-tag expressed in Sf9 cells • HDAC2 is a Class I HDAC that catalyzes the deacetylation of core histones, resulting in tightening of nucleosomal integrity, restricting access to transcription factors, and suppression of transcription.

1 ea

HDAC3/NCOR2 (human recombinant)

10009232

M_r: ~49.7 kDa **Purity:** ≥50% **Stability:** ≥6 months at -80°C

Source: Active recombinant protein containing a complex of human HDAC3 with a C-terminal His-tag and human NCOR2 amino acids 395-489 with an N-terminal GST-tag • HDAC3 is a Class I HDAC that is inactive alone and requires binding with the deacetylase activation domain of NcoR2 for activation.

1 ea


HDAC3 Polyclonal Antibody

13493

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human HDAC3 amino acids 2-17 • Host: rabbit • Cross Reactivity:(+) human HDAC3 • Application(s): ChIP, IP, and WB • HDAC3 is a class I HDAC that plays an important role in cell development, differentiation, programmed cell death, angiogenesis, and inflammation.

1 ea



HDAC4 (human recombinant)

10009652

M_r: 75.2 kDa **Purity:** ≥50% **Stability:** ≥6 months at -80°C

Source: Active N-terminal GST-tagged protein consisting of amino acids 627-1,085 expressed using a baculovirus expression system • HDAC4 is a Class IIa HDAC that can shuttle between the nucleus and cytoplasm, suggesting potential extranuclear functions by regulating the acetylation status of nonhistone substrates.

1 ea


HDAC4 Polyclonal Antibody

13494

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human HDAC4 amino acids 194-209 • Host: rabbit • Cross Reactivity: (+) human and mouse HDAC4 • Application(s): ChIP, IP, and WB • HDAC4 is a class II HDAC that can shuttle between the nucleus and cytoplasm, suggesting potential extranuclear functions by regulating the acetylation status of non-histone substrates.

1 ea



HDAC5 (human recombinant)

10009379

M_r: 51 kDa **Purity:** ≥90% **Stability:** ≥6 months at -80°C

Source: Active recombinant protein consisting of amino acids 657-1,123 with a C-terminal His-tag expressed in Sf9 cells • HDAC5 is a Class IIa HDAC that can shuttle between the nucleus and cytoplasm, suggesting potential extranuclear functions by regulating the acetylation status of non-histone substrates.

1 ea

HDAC6 (human recombinant)

10009465

M_r: ~159 kDa **Purity:** ≥80% **Stability:** ≥6 months at -80°C

Source: Active recombinant protein with an N-terminal GST-tag expressed in Sf9 cells • HDAC6 is a Class II HDAC that can shuttle between the nucleus and cytoplasm, suggesting potential extranuclear functions by regulating the acetylation status of non-histone substrates.

1 ea


HDAC6 Polyclonal Antibody

13499

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human HDAC6 amino acids 1-16 • Host: rabbit • Cross Reactivity: (+) human and mouse HDAC6 • Application(s): ChIP, IP, and WB • HDAC6 is a class II HDAC that can shuttle between the nucleus and cytoplasm, suggesting potential extranuclear functions by regulating the acetylation status of non-histone substrates.

1 ea


HDAC7 (Phospho-Ser¹⁵⁵) Polyclonal Antibody

13500

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: peptide from human HDAC7 containing phospho-Ser¹⁵⁵ • Host: rabbit • Cross Reactivity: (+) chimpanzee, bovine, canine, human, monkey, mouse, and rat HDAC7 • Application(s): WB • HDAC7 is a class IIa HDAC that plays a specific role in maintaining vascular integrity by repressing the expression of matrix metalloproteinase 10. It promotes repression mediated by transcriptional corepressor NCOR2 and is an efficient corepressor of the androgen receptor. It is also responsible for the deacetylation of lysine residues on the N-terminal part of the core histones.

1 ea




HDAC8 (human recombinant)

19380

M_r: 45.3 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C

Source: Active recombinant protein containing a C-terminal His-tag expressed in *E. coli* • HDAC8 is a Class I HDAC that catalyzes the deacetylation of core histones, resulting in tightening of nucleosomal integrity, restricting access to transcription factors, and suppression of transcription. It can also play an important role in mediating nuclear receptor functions by forming co-repressor complexes with nuclear receptors in the absence of ligands.

25 µg
50 µg
100 µg

HISTONE

Acetylation, Deacetylation, & Genomic Bistability

by [Thomas G. Brock, Ph.D.]

If you search PubMed using the keyword ‘epigenetics’, you’ll find over 5,500 papers. Closer inspection will reveal that almost half were published in the last two years! This rapid growth of the field reflects, at least in part, a broadening definition of the term. The original concept of epigenetics, developed in the late 20th century, focused on changes at the nuclear level that are heritable but do not involve changes in the DNA sequence. In that vision, such changes should be passed on for at least one generation, whether that is a mitotic division (at the cellular level) or an offspring (in complex organisms). In the new millennium, epigenetics has come to embrace an ever-increasing number of nuclear events which alter gene expression without changing DNA sequence. The absolute requirement for heritability has been lessened, replaced by an interest in histone and nucleosomal modification.

The study of histone acetylation and deacetylation has benefitted from these changes. The acetylation of histones has been described for decades, and early studies showed that acetylation status changed during cell division or development. Interest grew markedly when histone deacetylases were found to alter chromatin structure and regulate gene transcription. However, concerns that acetylation marks are not maintained as cells divide has brought into question whether these enzymes are truly involved in epigenetic processes.¹ This concern considers de/acetylation by itself, whereas models embracing the complexity of histone modification suggest roles for acetyl marks in heritable expression states. This article touches on these issues.

Histone Acetylases and Deacetylases

Acetylation refers to the addition of an acetyl group (CH₃CO) to organic compounds. Histone acetyltransferases (HATs) catalyze the transfer of an acetyl group from acetyl-CoA to the terminal amine on the side chain of lysine residues; since many HATs also acetylate lysine (denoted ‘K’) on other proteins, there is also a KAT nomenclature for most HATs (Table 1). While the twenty or so human HATs have a common general enzymatic capacity, they actually diverge in function, with some pairs of HATs showing overlap in histone targets. Importantly, most HATs contain distinct binding domains in addition to the acetyltransferase domain. For example, the MYST enzymes contain a C2HC-type zinc finger domain for binding other molecules, while CBP and P300 have bromodomains for binding acetylated lysine residues. This means that each HAT shows certain specificity in both where they bind and where they acetylate. The conversion of the positively charged lysine to non-charged acetyl-lysine, like the addition of negative phosphates to uncharged amino acids during phosphorylation, alters protein structure and interactions with other biomolecules. For example, acetylation of histones typically promotes the recruitment of effector proteins, relaxation of chromatin conformation, and an increase in transcription.

Like phosphorylation, acetylation is reversible. Histone deacetylases (HDACs, aka KDACs) are a smaller group of evolutionarily conserved enzymes. The human class I HDACs are homologous to the yeast enzyme Rpd3 and include

Coding Gene	KAT Name	Site of Histone Modification
HAT 1	KAT1	H2AK5, H4K5, H4K12
GCN5	KAT2A	H3K9, H3K14, H3K56, H4K5, H4K8, H4K12, H4K16, H4K91
PCAF	KAT2B	H3K9, H3K14
CBP	KAT3A	H3K14, H3K18, H3K27, H3K56, H4K5, H4K8, H4K12, H4K16
P300	KAT3B	H3K14, H3K18, H3K27, H3K56, H4K5, H4K8, H4K12, H4K16
TAF1	KAT4	H3K14
TIP60	KAT5	H2AK5, H4K5, H4K8, H4K12, H4K16
MYST3	KAT6A	H3K9, H3K14
MYST4	KAT6B	
MYST2	KAT7	H3K14, H4K5, H4K8, H4K12
MYST1	KAT8	H4K16
ELP3	KAT9	H3K9, H3K18
GTF3C4	KAT12	H3K14
NCOA1	KAT13A	H3K14
NCOA3	KAT13B	H3K14
CLOCK	KAT13D	H3K14
CDY1		
CDY2		
CDYL		
MGEA5		H4K8, H3K14
NAT10		

Table 1. Some human lysine acetyltransferases

HDAC1, 2, 3, and 8. Class II HDACs are homologous to yeast Hda1 and are divided into class IIa (HDAC4, 5, 7, 9) and class IIb (HDAC6 and 10) based on structure. The human class III HDACs, homologous to the yeast Sir2 protein, include the sirtuin family of NAD⁺-dependent protein deacetylases (SIRT1-3, 5, 6). The novel HDAC11 has a distinct structure and is a class IV HDAC. The HDACs often participate in the formation of transcriptional repressor complexes, inducing chromatin compaction through histone deacetylation, and silencing gene expression.

The Devil in the Details

Given an assortment of HATs and HDACs, what can we say about how they really work? Let’s start with a look at the histones. Nucleosomes consist of

DNA wrapped twice around an octamer composed of two sets of the core histones H2A, H2B, H3, and H4, with N-terminal tails of each histone projecting from the structure, as displayed in Figure 1. HATs acetylate lysine residues, which, intriguingly, are regularly spaced along the tails, commonly with 2 or 3 intervening residues. The spacing is apparent whether you look at the primary sequence given in letters or at the secondary structure of the nucleosome (lysines denoted in red). Why might they be spaced in this way? Moreover, most (if not all) lysine residues may be methylated or acetylated, but not both. Thus, most lysines have at least three alternative states: unmarked, acetylated, or methylated. Up to three methyl groups can accumulate on each lysine residue, with one methyltransferase (KMT) usually adding the first group and a second KMT executing the further additions. This suggests that the monomethylated lysine state is distinct from the polymethylated state, since the former is more readily transitioned back to the unmarked residue. If a given lysine state is viewed as a factor in affecting, say, protein binding, then the multiple states of the numerous lysines represents an abundance of information. This is evident without considering the impact of lysine ubiquitination or sumoylation, arginine methylation or citrullination, or phosphorylation of serines and threonines.

Unlike some enzymes, HATs and HDACs rarely work alone. For example, the HAT nuclear receptor coactivator 1 (NCOA1) directly binds dimerized nuclear receptors in a hormone-dependent fashion and, in this context, recruits additional HATs, including P300 and CBP.² HDAC1 and HDAC2 bind several partners to form the nucleosome remodeling and histone deacetylation (NuRD) complex that directs both histone deacetylation and ATP-dependent chromatin remodeling. Complexes containing HATs or HDACs bind to DNA or proteins, including modified histones. For example, the HAT-containing NuA3 complex binds methylated histone H3, localizing its HAT closer to the target site, H3K14.³ Importantly, and perhaps counterintuitively, HAT or HDAC complexes that bind to one nucleosome are not targeting residues on the same nucleosome. Instead, these complexes, typically much bigger than the histone octamers

which surround them, may act on multiple histone tails. This is not a trivial point, as it may link acetylation with heritability.

Genomic Bistability

An iconic image from the earliest days of epigenetics is that of the ‘epigenetic landscape’, presented by Conrad Hal Waddington in 1932 to visualize the external manifestation of genetic activity in an era when genes were considered discrete heritable units but their structures and functions were yet undiscovered (Figure 2). In this model, meant to pertain to cell differentiation, marbles (cells) move varying ways down a landscape whose contour is affected by genes. Details within the contours are further defined by factors above (‘epi-’) the fixed genetic level, and these details determine the final resting state of differentiation for each cell type.

Histone modifications may turn over rapidly, suggesting such marks might be too labile to carry information across cell divisions.⁴ How might

“Histone modifications may turn over rapidly, suggesting such marks might be too labile to carry information across cell divisions.”

acetylation contribute to heritable changes? This question relates to another topic: bistable systems. Some biological processes are analog and graded, where the output is a varying function of the input and the output returns to baseline when the input is removed. Other processes may be digital, yes-or-no by nature. For example, cell division, fertilization, and apoptosis are directional, digital decisions, as





Figure 2. Epigenetic valleys (left) are like bistable conditions (right). In epigenetics, differentiation is like contours, affected by underlying genes and ‘epi-’genetic factors, shaping cell fate. In bistable systems, two stable states (balls 1 and 3) may be attained from an unstable state (ball 2), depending on the stimulus.




Figure 1. Histone tails can have many marks. The amino termini project from the octamer core. N-terminal sequences for each of the core histones are given, with known sites of acetylation (*), methylation (**), either acetylation or methylation (!), or phosphorylation (p) indicated. Lysine residues within the nucleosome are colored in red.

changes between states are hard to reverse. Some biochemical reactions are reversibly bistable, moving between active and inactive states. A simple model of bistability, taken from physics, is reminiscent of Waddington’s landscape: particles may have three critical states, two of which represent minimum free energies and the third a maximum free energy (Figure 2). Mathematically, the unstable maximum must lie between the two stable minima and thus represents a barrier.

In an early application of bistability theory to epigenetic cell memory, Ian B. Dodd and colleagues start with the understanding that stability depends on positive feedback, where modified nucleosomes recruit enzymes that similarly modify nearby nucleosomes.⁵ There are several established examples of histone modifying enzymes that can both recognize and create the same modification.⁶⁻⁹ For example, CBP both binds acetylated proteins using a bromodomain and acetylates histones and other proteins. Robust stability requires cooperativity of modified nucleosomes in the modification reaction as well as modification beyond nearest neighbors.⁵ Current theoretical and empirical research is helping to refine our understanding of bistability in biological processes, including epigenetic signaling. For example, some systems employ multiple positive feedback loops, or combinations of positive and negative feedback loops, as well as spatial or temporal staggering of loops, in order to generate or enhance bistability.¹⁰⁻¹² The assortment of histone marks, in concert with complexes which can read, write, and erase these marks, clearly represents machinery capable of using bistability to provide heritability. The challenge will be to gain insight into how the signals work, how they go awry, and how they might be mended.


HDAC8 Inhibitor Screening Assay Kit

700230

Stability: ≥1 year at -80°C

Summary: Human HDAC8 is a class I HDAC and has been identified in a variety of human cancer tissues. Cayman's HDAC8 Inhibitor Screening Assay provides a convenient fluorescence-based method for screening HDAC8 inhibitors. The procedure requires only two easy steps, both performed in the same microplate. The fluorescent reaction product is analyzed with an excitation wavelength between 350-360 nm and an emission wavelength between 450-465 nm. Sufficient HDAC8 is provided for 96 assays.

96 wells



HDAC9 (human recombinant)

10009466

M_r: ~50.7 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C

Source: Active recombinant protein consisting of amino acids 604-1,066 with a C-terminal His-tag expressed in Sf9 cells • HDAC9 is a Class IIa HDAC that can shuttle between the nucleus and cytoplasm, suggesting possible extranuclear functions including regulating the acetylation status of non-histone substrates.

1 ea


HDAC11 Polyclonal Antibody

13504

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: peptide from human HDAC11 amino acids 182-199 • Host: rabbit • Cross Reactivity: (+) human, mouse, and rat HDAC11 • Application(s): WB • HDAC11 is a class IV HDAC that is expressed in kidney, heart, brain, skeletal muscle, and testis.


1 ea



HNHA

13295

[926908-04-5] Histone Deacetylase Inhibitor VI

MF: C₁₇H₂₁NO₂S **FW:** 303.4 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A cell-permeable inhibitor of HDAC activity (IC₅₀ = 100 nM)5 mg
10 mg
25 mg
50 mg


JGB1741

10641

[1256375-38-8] ILS-JGB-1741

MF: C₂₇H₂₄N₂O₂S **FW:** 440.6 **Purity:** ≥98%A crystalline solid **Stability:** ≥1 year at -20°C

Summary: A SIRT1-specific inhibitor (IC₅₀ = 15 μM); inhibits metastatic breast cancer MDA-MB 231 cell proliferation (IC₅₀ = 512 nM), dose-dependently increasing p53 acetylation and p53-mediated apoptosis in these cells

1 mg
5 mg
10 mg
25 mg


M 344

13174

[251456-60-7] D237, Histone Deacetylase Inhibitor III, MS 344

MF: C₁₆H₂₅N₃O₃ **FW:** 307.4 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C

Summary: An inhibitor of HDACs, inhibiting maize HDAC (IC₅₀ = 100 nM) as well as human HDAC1 (IC₅₀ = 46 nM); shows a 3-fold selectivity for HDAC6 over HDAC1

5 mg
10 mg
25 mg
50 mg

Metastasis Associated 1 Family Member 2


Polyclonal Antibody

13778

*MTA-L1 Protein, MTA2, p53 Target Protein in Deacetylase Complex, PID*Affinity-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: peptide from a portion of human MTA2 amino acids 650-700 • Host: rabbit • Cross Reactivity: (+) chimpanzee, human, and Rhesus monkey MTA2 • Application(s): IHC and WB • MTA2 is a nuclear protein that interacts with HDAC1 and HDAC2 and has a functional role in chromatin remodeling and deacetylase activity. It interacts with p53 and represses p53-dependent transcriptional activation, thereby regulating p53-mediated cell growth arrest and apoptosis.

1 ea




MS-275

13284

[209783-80-2] Entinostat, SNDX 275

MF: C₂₁H₂₀N₄O₃ **FW:** 376.4 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C

Summary: An inhibitor of HDACs that preferentially inhibits HDAC1 (IC₅₀ = 300 nM) over HDAC3 (IC₅₀ = 8 μM); does not inhibit HDAC8; induces p21/C1P1/WAF1, slowing cell growth, differentiation, and tumor development *in vivo*

1 mg
5 mg
10 mg
25 mg

NCOR2/SMRT (human recombinant)

11633

CTG Repeat Protein 26, Nuclear Receptor Corepressor 2, Silencing Mediator of Retinoic Acid and Thyroid Hormone Receptor, SMAP270, SMRT, T3 Receptor-Associating Factor, Thyroid-Retinoic-Acid-Receptor-Associated Corepressor, TRAC

M_r: 39 kDa **Purity:** ≥60% **Stability:** ≥6 months at -80°C

Source: Recombinant N-terminal GST-tagged protein consisting of amino acids 395-489 expressed in *E. coli* • NcoR2 is a transcriptional corepressor that plays an essential role in the regulation of development and metabolism. HDACs can mediate nuclear receptor functions by forming co-repressor complexes with nuclear receptors in the absence of ligands.

50 μg


Oxamflatin

13176

[151720-43-3] Metacept 3

MF: C₁₇H₁₄N₂O₄S **FW:** 342.4 **Purity:** ≥95%A crystalline solid **Stability:** ≥2 years at -20°C

Summary: A potent inhibitor of HDACs (IC₅₀ = 15.7 nM); has been shown to alter the expression of several genes whose products are involved in cell morphology, motility, apoptosis, and cell cycle control, reducing the proliferation of cancer cells

1 mg
5 mg
10 mg


p66α Polyclonal Antibody

13785

*GATA Zinc Finger Domain-Containing Protein 2A, Transcriptional Repressor p66α*IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: mouse p66α amino acids 572-585 • Host: rabbit • Cross Reactivity: (+) human, mouse, and rat p66α • Application(s): WB • p66 is one of the components of the MeCP1 complex, an HDAC core complex involved in methylated DNA silencing.

1 ea




PCI 34051

10444

[950762-95-5]

MF: C₁₇H₁₆N₂O₃ **FW:** 296.3 **Purity:** ≥95%A crystalline solid **Stability:** ≥1 year at -20°C

Summary: A potent HDAC8 inhibitor (IC₅₀ = 0.01 μM) with >200-fold selectivity over HDAC isoforms 1, 2, 3, 6, and 10 (IC₅₀s = 4, >50, >50, 2.9, and 13 μM, respectively); induces caspase-dependent apoptosis in cell lines derived from T-cell lymphomas or leukemias (GI₅₀s = 2.4 - 4 μM)

5 mg
10 mg
50 mg
100 mg


Pimelic Diphenylamide 106

13212

[937039-45-7] TC-H 106, Histone Deacetylase Inhibitor VII

MF: C₂₀H₂₅N₃O₂ **FW:** 339.4 **Purity:** ≥98%A crystalline solid **Stability:** ≥1 year at -20°C

Summary: A slow, tight-binding inhibitor of class I HDACs, progressively binding HDACs and remaining bound after wash-out; inhibits class I HDACs (IC₅₀ = 150, 760, 370, and 5,000 nM for HDAC1, 2, 3, and 8, respectively) but not class II HDACs (IC₅₀ >180 μM for HDAC4, 5, and 7)

1 mg
5 mg
10 mg
25 mg


Pyroxamide

13870

[382180-17-8]

MF: C₁₃H₁₉N₃O₃ **FW:** 265.3 **Purity:** ≥95%A crystalline solid **Stability:** ≥2 years at -20°C

Summary: An inhibitor of HDACs, including HDAC1 (IC₅₀ = 0.1-0.2 μM); induces growth suppression and cell death of certain types of cancer cells in culture

5 mg
10 mg
25 mg
50 mg


SAHA

10009929

[149647-78-9] Suberoylanilide Hydroxamic Acid, Vorinostat, Zolinza™

MF: C₁₄H₂₀N₂O₃ **FW:** 264.3 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C

Summary: An HDAC inhibitor of class I and class II HDACs at around 50 nM; arrests cell growth in a wide variety of transformed cells in culture at 2.5-5.0 μM

100 mg
250 mg
500 mg
1 g


SAHA-BPpyne

10675

[930772-88-6] Suberoylanilide Hydroxamic Acid-BPpyne

MF: C₂₇H₃₁N₃O₃ **FW:** 477.6 **Purity:** ≥98%A solution in methanol **Stability:** ≥1 year at -20°C

Summary: A SAHA derivative with a benzophenone crosslinker and an alkyne tag to be used for profiling HDAC activities in proteomes and live cells; labels HDAC complex proteins both in proteomes at 100 nM and in live cells at 500 nM; IC₅₀ = ~3 μM for inhibition of HDAC activity in HeLa cell nuclear lysates in an HDAC activity assay

50 μg
100 μg
500 μg


coumarin-SAHA

10671

[1260635-77-5] coumarin-Suberoylanilide Hydroxamic Acid

MF: C₁₈H₂₂N₂O₅ **FW:** 346.4 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C

Summary: A fluorescent probe that competitively binds HDAC; demonstrates fluorescence excitation and emission maxima of 325 and 400 nm, respectively, which is quenched by 50% when bound to HDAC

1 mg
5 mg
10 mg
25 mg

4-iodo-SAHA

10495


[1219807-87-0] 4-iodo-Suberoylanilide Hydroxamate

MF: C₁₄H₁₉IN₂O₃ FW: 390.2 Purity: ≥98%

A crystalline solid Stability: ≥2 years at -20°C

Summary: A hydrophobic derivative of the class I and class II HDAC inhibitor SAHA that demonstrates >60% inhibition of HDAC1 and HDAC6 activity in a deacetylase activity assay; inhibits proliferation of SKBR3 breast-derived, HT29 colon-derived, and U937 leukemia cell lines with EC₅₀ values of 1.1, 0.95, and 0.12 μM, respectively

50 mg
100 mg
250 mg
500 mg



Salermide

13178


[1105698-15-4]

MF: C₂₆H₂₂N₂O₂ FW: 394.5 Purity: ≥98%

A crystalline solid Stability: ≥2 years at -20°C

Summary: An inhibitor of SIRT1 and SIRT2, causing tumor-specific apoptotic cell death; causes 90% apoptosis within 72 hours (IC₅₀ = 20 μM) by reactivating proapoptotic genes that are repressed by SIRT1 in MOLT4 leukemia cells

5 mg
10 mg
50 mg
100 mg



SB 939

10443


[929016-96-6] Pracinostat

MF: C₂₀H₃₀N₄O₂ FW: 358.5 Purity: ≥98%

A crystalline solid Stability: ≥2 years at -20°C

Summary: A pan-HDAC inhibitor (IC₅₀ = 77 nM in an *in vitro* HDAC1 activity assay) that prevents proliferation of ovarian (A2780), colon (COLO 205 and HCT-116), and prostate cancer (PC-3) cell lines at IC₅₀ values of 0.48, 0.56, 0.48, and 0.34 μM, respectively; binds all HDAC isozymes with similar affinity (K_s = 16-28 nM) with the exception of HDAC6 and 7 (K_s = 247 and 104 nM, respectively)

1 mg
5 mg
10 mg
25 mg



Scriptaid

10572


[287383-59-9] GCK 1026

MF: C₁₈H₁₈N₂O₄ FW: 326.4 Purity: ≥98%

A crystalline solid Stability: ≥2 years at -20°C

Summary: An HDAC inhibitor that has an optimal effective concentration of 6-8 μM in a cell-based assay, is less toxic than trichostatin A, and works in a wide variety of biological systems; induces cell cycle arrest in cancer cells *in vitro* and *in vivo*; facilitates the cloning of inbred mouse strains produced by somatic cell nuclear transfer

1 mg
5 mg
10 mg
25 mg




SIRT Direct Fluorescent Screening Assay Kits

The sirtuins represent a distinct class of trichostatin A-insensitive lysyl-deacetylases (class III HDACs) that catalyze a reaction coupling lysine deacetylation to the formation of nicotinamide and O-acetyl-ADP-ribose. Cayman's Direct Fluorescent Screening Assay Kits provide a convenient fluorescence-based method for screening SIRT inhibitors or activators. The procedure requires only two easy steps, both performed in the same microplate. In the first step, the substrate is incubated with human recombinant SIRT along with its cosubstrate NAD⁺. Deacetylation sensitizes the substrate such that treatment with the developer in the second step releases a fluorescent product. The fluorophore can be analyzed with an excitation wavelength of 350-360 nm and an emission wavelength of 450-465 nm.

SIRT1 Direct Fluorescent Screening Assay Kit 10010401

Stability: ≥1 year at -80°C


96 wells



SIRT2 Direct Fluorescent Screening Assay Kit 700280

Stability: ≥1 year at -80°C


96 wells



SIRT3 Direct Fluorescent Screening Assay Kit 10011566

Stability: ≥1 year at -80°C


96 wells



SIRT6 Direct Fluorescent Screening Assay Kit 700290

Stability: ≥1 year at -80°C

96 wells




SIRT1 FRET-Based Screening Assay Kit 10010991

Stability: ≥1 year at -80°C

Summary: Human SIRT1 is the homolog of yeast Sir2 and has been shown to regulate the activity of the p53 tumor suppressor and inhibit apoptosis. Small molecule activators of SIRT1, such as resveratrol, extend lifespan in yeast and *C. elegans* in a manner that resembles caloric restriction. Cayman's SIRT1 FRET-based Screening Assay provides a convenient fluorescence-based method for screening SIRT1 inhibitors or activators. The procedure requires only two easy steps, both performed in the same microplate. In the first step, the substrate, which is coupled to the fluorophore and quencher, is incubated with human recombinant SIRT1 along with its cosubstrate NAD⁺. Deacetylation sensitizes the substrate such that treatment with the developer in the second step results in the separation of the quencher and fluorophore. The resulting fluorescence is analyzed using an excitation wavelength of 335-345 nm and emission wavelength of 440-465 nm.

96 wells



SIRT1 (human recombinant) 10011190

NAD-dependent Deacetylase 1, Silent Information Regulator 2, SIR2L1, SIR2-like Protein 1, Sirtuin 1

M_r: 89.2 kDa Purity: ≥60% Stability: ≥9 months at -80°C

Source: Active recombinant N-terminal GST-tagged enzyme amino acids 193-747 expressed in *E. coli* • SIRT1 is the human sirtuin with the greatest homology to yeast Sir2 and has been shown to regulate the activity of the p53 tumor suppressor and inhibit apoptosis.

25 units
50 units
100 units

SIRT2 (human recombinant) 10011191

NAD-dependent Deacetylase 2, Silent Information Regulator 2, SIR2L2, SIR2-like Protein 2, Sirtuin 2

M_r: 44.2 kDa Purity: ≥90% Stability: ≥9 months at -80°C

Source: Active recombinant N-terminal His-tagged enzyme amino acids 2-389 expressed in *E. coli* • SIRT2 is a cytoplasmic protein responsible for the deacetylation of histone H4 and α-tubulin, a modification important for controlling the cell cycle. SIRT2 co-localizes with HDAC6 and microtubules and functions as a mitotic checkpoint in preventing chromosomal instability that can lead to hyperploid cells.


25 μg
50 μg
100 μg

SIRT3 (human recombinant)

10011194

Mitochondrial Nicotinamide Adenine Dinucleotide-dependent Deacetylase, NAD-dependent Deacetylase 3, Silent Information Regulator 3, SIR2L3, SIR2-like Protein 3, Sirtuin 3
M_r: 37 kDa **Purity**: ≥90% **Stability**: ≥9 months at -80°C

Source: Active recombinant N-terminal His-tagged enzyme amino acids 101-399 expressed in *E. coli* • SIRT3, is a mitochondrial protein that is synthesized as an enzymatically inactive protein. Human SIRT3 is activated by a matrix-processing peptidase. The constitutive expression of SIRT3 promotes the expression of PGC-1α, UCP1, and other genes involved in mitochondrial functions, indicating that SIRT3 modulates adaptive thermogenesis in BAT.

25 µg
50 µg
100 µg

Lane 1: MW Markers
Lane 2: SIRT3 (5 µg)
Lane 3: SIRT3 (2.5 µg)
Lane 4: SIRT3 (1.25 µg)


SIRT4 (human recombinant)

10317

NAD-dependent ADP-ribosyltransferase Sirtuin 4, Silent Information Regulator 4, SIR2L4, SIR2-like Protein 4, Sirtuin 4

M_r: 61.9 kDa **Purity**: ≥95% **Stability**: ≥9 months at -80°C

Source: Recombinant N-terminal GST-tagged enzyme expressed in *E. coli* • SIRT4 is a mitochondrial ADP-ribosyltransferase responsible for the transfer of ADP-ribose from NAD to specific substrates such as glutamate dehydrogenase.

25 µg
50 µg
100 µg

Lane 1: MW Markers
Lane 2: SIRT4 (4 µg)
Lane 3: SIRT4 (2 µg)
Lane 4: SIRT4 (1 µg)

SIRT5 (human recombinant)


10318

NAD-dependent Deacetylase 5, Silent Information Regulator 5, SIR2L5, SIR2-like Protein 5, Sirtuin 5

M_r: 60.6 kDa (GST-tag); 26 kDa (native) **Purity**: ≥90%

Stability: ≥9 months at -80°C

Source: Recombinant N-terminal GST-tagged enzyme expressed in *E. coli* • SIRT5 is located in the mitochondrial matrix and its functions are largely still being elucidated, however a few promising substrates have been studied. SIRT5 has been shown to deacetylate carbamoyl phosphate synthetase I, activating the enzyme to catalyze the first step of the urea cycle.

25 µg
50 µg
100 µg

Lane 1: MW Markers
Lane 2: SIRT5 (4 µg)
Lane 3: SIRT5 (2 µg)
Lane 4: SIRT5 (1 µg)


SIRT6 (human recombinant)

10315

NAD-dependent Deacetylase 6, Silent Information Regulator 6, SIR2L6, SIR2-like Protein 6, Sirtuin 6

M_r: 43.7 kDa **Purity**: ≥95% **Stability**: ≥6 months at -80°C

Source: Active recombinant N-terminal His-tagged enzyme amino acids 1-355 expressed in *E. coli* • SIRT6 associates specifically with telomeres and functions at chromatin to decrease NF-κB signaling. Mammalian cells depleted of SIRT6 display abnormal telomere structures similar to defects found in Werner syndrome, a premature aging disorder associated with a shortened life span.

25 µg
50 µg
100 µg

Lane 1: MW Markers
Lane 2: SIRT6 (2 µg)
Lane 3: SIRT6 (1 µg)
Lane 4: SIRT6 (0.5 µg)


SIRT7 (human recombinant)

10316

NAD-dependent deacetylase 7, Silent Information Regulator 7, SIR2L7, SIR2-like protein 7, Sirtuin 7

M_r: 49.3 kDa **Purity**: ≥85% **Stability**: ≥6 months at -80°C

Source: Recombinant N-terminal His-tagged enzyme amino acids 2-400 expressed in *E. coli* • SIRT7 activates transcription by RNA polymerase I and deacetylates p53. It prevents progressive deterioration of the heart, and is suggested to play an important role in regulation of stress responses and cell death in the heart.

25 µg
50 µg
100 µg

Lane 1: MW Markers
Lane 2: SIRT7 (1 µg)
Lane 3: SIRT7 (2 µg)
Lane 4: SIRT7 (5 µg)

SIRT7 Polyclonal Antibody


13477

SIR2, Sirtuin 7

Protein G-purified IgG **Stability**: ≥1 year at -20°C

Summary: Antigen: human SIRT7 amino acids 35-51 and 361-377 • Host: rabbit • Cross Reactivity: (+) human SIRT7 • Application(s): WB • SIRT7 is a member of the sirtuin family of proteins, which are able to metabolize NAD⁺. Reports of histone-activated SIR2-mediated NAD⁺ metabolism and NAD⁺-activated SIR2-mediated histone deacetylation suggest a coupled reciprocal activation mechanism involving interactions of SIR2 with NAD⁺ and the N-ε-acetyl-lysine groups of acetylated histone.

1 ea



Lane 1: Human liver homogenate
Lane 2: PBMC lysate

Sirtinol


10523

[410536-97-9] *Sir Two Inhibitor Naphthol*

MF: C₂₆H₂₂N₂O₂ **FW**: 394.5 **Purity**: ≥98%

A crystalline solid **Stability**: ≥2 years at -20°C

Summary: A cell-permeable inhibitor of sirtuin NAD⁺-dependent deacetylases, inhibiting the yeast sirtuin Sir2p with an IC₅₀ value of 68 µM and the human sirtuins SIRT1 and SIRT2 with IC₅₀ values of 131 and 38 µM, respectively; does not alter HDAC1 activity

1 mg
5 mg
25 mg**Sodium Butyrate**


13121

[156-54-7] *Butyric Acid*

MF: C₄H₈O₂ • Na **FW**: 111.1 **Purity**: ≥95%

A crystalline solid **Stability**: ≥2 years at room temperature

Summary: A short chain fatty acid that inhibits HDACs, induces growth arrest, differentiation and apoptosis in cancer cells, and suppresses inflammation by reducing the expression of pro-inflammatory cytokines

50 g
100 g
250 g
500 g**Splitomicin**


13168

[5690-03-9] *1-Naphthalenepropanoic Acid*

MF: C₁₃H₁₀O₂ **FW**: 198.2 **Purity**: ≥98%

A crystalline solid **Stability**: ≥2 years at -20°C

Summary: A small molecule inhibitor of Sir2p HDAC activity, displaying higher activity *in vivo* (minimal inhibitory concentration = 0.49 µM) than *in vitro* (IC₅₀ = 60 µM); has diverse effects on mammalian cells

5 mg
10 mg
25 mg
50 mg**Suberohydroxamic Acid**


10574

[38937-66-5] *SBHA, Suberic bis-Hydroxamic Acid*

MF: C₈H₁₆N₂O₄ **FW**: 204.2 **Purity**: ≥98%

A crystalline solid **Stability**: ≥2 years at -20°C

Summary: A competitive HDAC inhibitor that inhibits HDAC1 (IC₅₀ = 0.25 µM) and HDAC3 (IC₅₀ = 0.30 µM); causes cell differentiation, cell cycle arrest, or apoptosis

100 mg
250 mg
500 mg
1 g**Tenovin-1**


13085

[380315-80-0]

MF: C₂₀H₂₃N₃O₂S **FW**: 369.5 **Purity**: ≥98%

A crystalline solid **Stability**: ≥2 years at -20°C

Summary: A small molecule activator of p53 that decreases the growth of BL2 Burkitt's lymphoma and ARN8 melanoma cells; inhibits the deacetylase activity of purified human SIRT1 and 2

5 mg
10 mg
50 mg
100 mg**Tenovin-6**


13086

[1011557-82-6]

MF: C₂₅H₃₄N₄O₂S **FW**: 454.6 **Purity**: ≥95%

A crystalline solid **Stability**: ≥2 years at -20°C

Summary: An analog of tenovin-1; elevates p53 activity in MCF-7 cells at 10 µM and reduces growth of ARN8 melanoma xenograft tumors in SCID mice at a dose of 50 mg/kg

1 mg
5 mg
10 mg
25 mg**Trichostatin A**


89730

[58880-19-6] *TSA*

MF: C₁₇H₂₂N₂O₃ **FW**: 302.4 **Purity**: ≥98%

A crystalline solid **Stability**: ≥1 year at -20°C

Summary: A potent, reversible inhibitor of HDAC (IC₅₀ = 70 nM)


500 µg
1 mg
5 mg
10 mg**Tubastatin A (trifluoroacetate salt)**

10559

MF: C₂₀H₂₁N₃O₂ • CF₃COOH **FW**: 449.2 **Purity**: ≥95%

A crystalline solid **Stability**: ≥2 years at -20°C

Summary: Potent HDAC6 inhibitor (IC₅₀ = 15 nM) with 1,000-fold selectivity against all other HDAC isoforms (IC₅₀'s >16 µM), excluding HDAC8 (IC₅₀ = 0.9 µM); induces α-tubulin hyperacetylation at 2.5 µM in primary cortical neuron cultures; displays dose-dependent neuronal protection of primary cortical neuron cultures at 5-10 µM

500 µg
1 mg
5 mg
10 mg**Valproic Acid (sodium salt)**


13033

[1069-66-5] *2-Propylvaleric Acid, Sodium Valproate*

MF: C₈H₁₅O₂ • Na **FW**: 166.2 **Purity**: ≥95%

A crystalline solid **Stability**: ≥2 years at -20°C

Summary: An analog of valeric acid, long used as an anti-convulsant; inhibits Class I HDACs with an IC₅₀ value of ~2 mM; also inhibits GSK3 and depletes cellular 1,4,5-IP₃

10 g
25 g
50 g
100 g

Acetyl Readers

Antibodies

30 13497 BRD4/HUNK1 Polyclonal Antibody

Biochemicals

32 11187 (+)-JQ1
32 11232 (-)-JQ1
32 11155 PFI-1


Kits

24 600710 BAZ2B bromodomain TR-FRET Assay Kit
25 600500 BRD2 bromodomain 1 TR-FRET Assay Kit
28 600510 BRD2 bromodomain 2 TR-FRET Assay Kit
28 600630 BRD3 bromodomain 1 TR-FRET Assay Kit
29 600520 BRD4 bromodomain 1 TR-FRET Assay Kit
30 600530 BRD4 bromodomain 2 TR-FRET Assay Kit
30 600650 BRD4 bromodomain 2 TR-FRET Assay Kit
31 600730 BRM bromodomain TR-FRET Assay Kit


Proteins

24 11917 BAZ2B bromodomain (human recombinant)
25 11918 BPTF bromodomain (human recombinant)
25 11507 BRD1 bromodomain (human recombinant)
25 11071 BRD2 bromodomain 1 (human recombinant; GST-tagged)
25 11069 BRD2 bromodomains 1 and 2 (human recombinant)
25 11070 BRD2 bromodomain 2 (human recombinant; GST-tagged)
28 11285 BRD3 bromodomain 1 (human recombinant)
28 11068 BRD4 bromodomain 1 (human recombinant; GST-tagged)
29 11720 BRD4 bromodomain 1 (human recombinant; His-tagged)
29 11052 BRD4 bromodomains 1 and 2 (human recombinant)
29 11066 BRD4 bromodomain 2 (human recombinant; GST-tagged)
29 11721 BRD4 bromodomain 2 (human recombinant; His-tagged)
30 11509 BRD9 bromodomain (human recombinant)
30 11548 BRDT bromodomain 1 (human recombinant)
31 11649 BRDT bromodomain 2 (human recombinant)
31 11284 BRG1 bromodomain (human recombinant)
31 11289 BRM bromodomain (human recombinant)
31 11650 BRPF3 bromodomain 1 (human recombinant)
32 14133 CECR2 bromodomain (human recombinant)
32 11288 CREB-binding protein bromodomain (human recombinant)
32 11920 PCAF bromodomain (human recombinant)
33 11652 Polybromo-1D bromodomain 1 (human recombinant)
33 11922 TAF1 bromodomain 1 (human recombinant)
33 11653 TRIM24 bromodomain (human recombinant)
33 11549 WDR9 bromodomain 2 (human recombinant)


BAZ2B bromodomain (human recombinant) 11917

*Bromodomain Adjacent to Zinc Finger Domain 2B, KIAA1476, WALp4***M_r**: 40.2 kDa **Purity**: ≥80% **Stability**: ≥1 year at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 2,064-2,168 expressed in *E. coli* • BAZ2B is expressed in heart, skeletal, testis, and pancreatic tissues. A rare allele of BAZ2B has been identified to be a predictor of Sudden Cardiac Death. The full-length protein contains several DNA-targeting domains, including a methyl-CpG binding domain, a DNA-binding DDT domain, and a tandem PHD-bromodomain.25 µg
50 µg
100 µg


BAZ2B bromodomain TR-FRET Assay Kit 600710

*Bromodomain Adjacent to Zinc Finger Domain 2B, KIAA1476, WALp4***Stability**: ≥6 months at -80°C **Z' Factor**: 0.6**Summary**: BAZ2B is a novel bromodomain-containing protein whose biological function, while not yet confirmed, is believed to function similar to ACF1, the *Drosophila* BAZ2B ortholog. ACF complexes play a role in establishing regular nucleosome spacing during chromatin assembly and influencing different remodeling outcomes at target loci. Cayman's BAZ2B bromodomain TR-FRET Assay Kit is a homogeneous, TR-FRET assay method amenable to rapid characterization of inhibitors of bromodomain/acetylated peptide interaction in a high throughput format.384 wells
1,920 wells
9,600 wells


BPTF bromodomain (human recombinant) 11918

*Bromodomain PHD Finger Transcription Factor***M_r**: 40.1 kDa **Purity**: ≥80% **Stability**: ≥1 year at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 2,796-2,907 expressed in *E. coli* • BPTF is the largest component of the NURF chromatin remodeling complex. It includes adjacent PHD and bromodomains which recognize trimethylation of H3K4 or acetylation of lysines in histone 4, respectively. BPTF is an essential regulator of gene expression in early mouse embryo.25 µg
50 µg
100 µg


BRD1 bromodomain (human recombinant) 11507

*BR140-like protein, BRL, Bromodomain and PHD finger-containing protein 2, Bromodomain containing 1, BRPF1, BRPF2, DKFZp686F0325***M_r**: 41.1 kDa **Purity**: ≥95% **Stability**: ≥1 year at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 556-680 expressed in *E. coli* • BRD1 is a bromodomain containing protein that has been identified as a susceptibility gene in neurological disorders, such as schizophrenia and bipolar affective disorder.25 µg
50 µg
100 µg


BRD2 bromodomain 1 (human recombinant; GST-tagged) 11071

*Bromodomain containing 2, RING3, RNF3***M_r**: 42.4 kDa **Purity**: ≥95% **Stability**: ≥6 months at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 65-187 expressed in *E. coli* • The isolated individual or tandem bromodomains of BRD2 have been shown to bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters.25 µg
50 µg
100 µg


BRD2 bromodomain 1 TR-FRET Assay Kit 600500

*Bromodomain containing 2, RING3, RNF3***Stability**: ≥6 months at -80°C **Z' Factor**: 0.71**Summary**: The isolated individual or tandem bromodomains of BRD2 bind acetylated histone tails, which couples histone acetylation marks to the transcriptional regulation of target promoters. Small molecule inhibitors of bromodomain interactions hold promise as useful therapeutics for human disease. Cayman's BRD2 Bromodomain 1 TR-FRET Assay Kit is a homogeneous, TR-FRET assay method amenable to rapid characterization of inhibitors of bromodomain/acetylated peptide interaction in a high throughput format.384 wells
1,920 wells
9,600 wells

BRD2 bromodomain 1 and 2 (human recombinant) 11069

*Bromodomain containing 2, RING3, RNF3***M_r**: 71.2 kDa **Purity**: ≥90% **Stability**: ≥6 months at -80°C**Source**: Recombinant GST-tagged protein consisting of amino acids 65-459 expressed in *E. coli* • The isolated individual or tandem bromodomains of BRD2 have been shown to bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters.25 µg
50 µg
100 µg

BRD2 bromodomain 2 (human recombinant; GST tagged) 11070

*Bromodomain containing 2, RING3, RNF3***M_r**: 42 kDa **Purity**: ≥95% **Stability**: ≥6 months at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 339-459 expressed in *E. coli* • The isolated individual or tandem bromodomains of BRD2 have been shown to bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters.25 µg
50 µg
100 µg

Bromodomains, ChromoHub, & STRING

by [Thomas G. Brock, Ph.D.]

Our Epigenetic Screening Laboratory is designed to be flexible and innovative in helping you meet your research needs.

More information on products is available at www.caymanchem.com/epigenetics for products and www.caymanchem.com/episcreen for epigenetic screening and profiling services.

In the everyday parlance of students of epigenetics, the enzymes that put marks on histones, like acetyltransferases, are 'writers', while those that remove the marks are 'erasers'. In keeping with the literary motif, the proteins that interact with the marks are called 'readers'. Bromodomains (BRDs) are the modules on certain proteins which act as the readers of ϵ -N-lysine acetylation marks placed on histones as well as other proteins. Proteins containing BRDs, and their co-regulators, are involved in chromatin remodeling, modulation of transcription, and cell signaling.¹ Dysfunction involving BRDs has been implicated in broad categories of diseases, including cancer, obesity, type 2 diabetes, and inflammation.^{2,3} Examples of more specific diseases associated with mutations or fusions of genes expressing proteins with BRDs are veno-occlusive disease with immunodeficiency syndrome (SP110), X-linked mental retardation (BRWD3), and infant pro-B acute lymphoblastic leukemia (MLL).^{4,6} Thus, these domains, proteins, and interacting complexes are of great interest. This brief article introduces BRDs, as well as some resources that are useful for their study.




Figure 1. The structures of the BRDs of CBP and BRPF1, the latter bound to acetylated lysine (Kac), reveal a common overall structure. From 3DWY1 and 2RS9 in RCSB Protein Data Bank

Bromodomains

In 1992, John Tamkun, a young scientist studying gene regulation in *Drosophila* with James Kennison and others, described a new gene involved in the transcriptional activation of homeotic genes, which they named brahma after the prominent Hindu god.⁷ They noted that brahma contained a 77 amino acid (aa) motif which is found in other *Drosophila*, yeast, and human regulatory proteins, suggesting that this pattern may be characteristic of a new family of regulatory proteins. Now, over 20 years later, these BRDs are indeed known to be distinguishing features of several transcriptional regulators. However, they are also found on proteins that are not directly involved in gene transcription, including acetyltransferases and chromatin remodeling factors.¹ Just as Brahma has four heads, the domains named after him typically consist of four helices (Figure 1). Two loops, one long and the other short, bind the acetylated lysine. The presentation of BRDs on this scale may be misleading. Proteins which contain BRDs range from 807 to 2,969 amino acids, averaging about 1,500 amino acids. Given that histones are ~ 130 amino acids, a protein binding a histone using a BRD is like a person, averaging 150 lbs but ranging up to 296 lbs, biting a 13 lb turkey at a Thanksgiving dinner.

Of course, the BRD/acetylated lysine association is more of a lock and key interaction than mouth on drumstick. Some histone acetylation sites are very specific in which proteins can bind to them, while others are more promiscuous.¹ Similarly, proteins vary in how specific they are for acetylation sites. In addition, recent evidence indicates that neighboring posttranslational modifications, such as acetylation and phosphorylation, affect binding. That is, BRDs recognize patterns of marks on histones. For example, the combination of acetylation on H3K9 and H4K16 with phosphorylation at H3S10 provides the requisites for BRD4 binding, recruitment of TEFb, and gene transcription.⁸ The recruitment of partner proteins is an important element in signaling through BRDs. This topic raises the questions: where can one find information about these interactions, how important are they, and what is their impact?

ChromoHub

ChromoHub (http://apps.thesgc.org/resources/phylogenetic_trees/index.php) is a data hub for navigators of chromatin-mediated signaling.⁹ Developed by researchers at the Structural Genomics Consortium at the University of Toronto, this site is a valuable resource for anyone trying to keep up with research related to chromatin-mediated regulation of epigenetic mechanisms. It initially segregates data according to 19 different domains involved in writing, reading, or erasing histone marks, generating phylogenetic trees which are either gene-based, produced using ClustalW, or domain-based, developed from seed sequence alignments using ICM. These trees are then overlaid with the information of interest to the user.

The first option listed in ChromoHub refers to Disease Associations and provides references for each protein indicating involvement in inflammation, cancer, viral infections, neurological diseases, metabolic disorders, immune disorders, or regenerative medicine. Perhaps not surprisingly, BRD-containing proteins have been predominantly linked with cancer. This may in part reflect low-hanging fruit, since associating chromatin remodeling proteins to cancer is easier than identifying more specific roles for each protein. Several additional options are available at ChromoHub to extend the information related to cancer and provide quality information. For example, evaluation of gene fusions in cancer shows that there is abundant evidence for fusions of MYST histone acetyltransferases with CBP in leukemias. PBRM1 shows the highest rate of somatic mutations in cancer, found in 16% of 288 patients with kidney renal clear cell carcinoma. Of course, these cancer links will be useful for all of the chromatin-modifying proteins.

Moving beyond the disease links, a user might choose options regarding funding and "jump in activity" to evaluate which BRD-containing proteins have been of greatest recent interest. ChromoHub shows that proteins like CECR2 and PBRM1 have shown upswings in interest during the past year (Figure 2). The funding option also provides publications

for each protein, segregated by year. According to abstracts provided by ChromoHub, ATAD2 has been identified as a predictor of poor prognosis in breast and lung cancers, helping to explain a jump in funding for ATAD2 research in 2012. ChromoHub also indicates that BAZ2B is drawing interest for its documented link with sudden cardiac death, while money to study CECR2 is going to a computational biologist to do an evolutionary analysis of protein superfamilies. By using options like these, a scientist can find recent information specifically about each protein and begin to understand where interests in the field of chromatin remodeling are changing.

STRING 9.0

Like a scientist who leverages teamwork to accomplish great things, proteins with BRDs work with other proteins to achieve their goals. In ChromoHub, the option called 'protein interactions' allows one to interrogate potential partners for each protein, using the application STRING.¹⁰ For example, the BRD and PHD finger-containing protein 1, BRPF1, links to several proteins by varying degrees (Figure 3). Different colored lines indicate different evidence supporting the relationship, with red lines indicating experimentation, blue lines for databases, and green lines for textmining;




Figure 3. Potential protein interactions with BRD and PHD finger-containing protein 1 (BRPF1) are given by STRING 9.0.¹⁰




Figure 2. A partial presentation of a ChromoHub BRD tree shows where information is available regarding disease associations (graphs), inhibitors (purple boxes), funding (dollar signs), and recent jump in activity (red triangle) for some proteins with BRDs. Hovering over a symbol opens a link to annotated information.


other criteria also exist. Correlation scores are also provided, as well as basic information about each of the proteins. A quick examination of this page provides the bottom line: BRPF1 is a component of a complex that also contains MYST3 (aka MOZ, KAT6A), MEAF6, MYST4 (aka MORE, KAT6B), and ING5. Additional information on the page says that BRPF1 is a "component of the MOZ/MORF complex which has a histone H3 acetyltransferase activity" and that it "positively regulates the transcription of RUNX1 and RUNX2." In short, STRING provides enough information for an inquisitive scientist to begin to understand the actions of any given chromatin-modifying protein.

Of course, STRING has limitations. The proteins linked to BRPF1 by textmining only (green lines) merely appear in the same abstract together. Thus, BRPF1 and intersectin 1 (ITSN1) are both substrates of granzyme B, while BRPF1 and synaptonemal complex protein 1 (SYCP1) genes are regulated by hydrogen peroxide. These kinds of associations redefine the meaning of 'protein interactions', and not in a favorable way. On the upside, the website is updated every two years. STRING 9.0 was generated in 2011, which means that STRING 9.1 is due out soon. In fact, a preview version of STRING 9.1 reveals that BRPF1, with the MOZ/MORF complex, associates with its target histone H3, and associations with ITSN1 and SYCP1 are no longer listed (although new tenuous associations with a collagen protein, COL4A3, and a demethylase, JMJD1C, are added). Liberally speaking, STRING errs on the side of giving extra information.


Cayman Chemical

Perhaps the best source for products and services that you can actually use to advance your research is Cayman, which has a dedicated team of scientists studying epigenetics with an emphasis on BRDs. They have engineered a large collection of human recombinant BRDs, synthesized inhibitors, and developed assays to screen for novel modulators of BRD/acetyllsine binding. Most of these products are listed in the Acetyl Readers section of this catalog.


BRD2 bromodomain 2 TR-FRET Assay Kit 600510

*Bromodomain containing 2, RING3, RNF3***Stability:** ≥6 months at -80°C **Z' Factor:** 0.82**Summary:** The isolated individual or tandem bromodomains of BRD2 bind acetylated histone tails, which couples histone acetylation marks to the transcriptional regulation of target promoters. Small molecule inhibitors of bromodomain interactions hold promise as useful therapeutics for human disease. Cayman's BRD2 Bromodomain 2 TR-FRET Assay Kit is a homogeneous, TR-FRET assay method amenable to rapid characterization of inhibitors of bromodomain/acetylated peptide interaction in a high throughput format.384 wells
1,920 wells
9,600 wells


BRD3 bromodomain 1 (human recombinant) 11285

*Bromodomain containing protein 3, ORFX, RING3L, RING3-like protein***M_r:** 41.2 kDa **Purity:** ≥95% **Stability:** ≥1 year at -80°C**Source:** Recombinant N-terminal GST-tagged protein consisting of amino acids 24-144 expressed in *E. coli* • The isolated individual or tandem bromodomains of BRD3 have been shown to bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters.25 µg
50 µg
100 µg


BRD3 bromodomain 1 TR-FRET Assay Kit 600630

*Bromodomain containing protein 3, ORFX, RING3L, RING3-like protein***Stability:** ≥6 months at -80°C **Z' Factor:** 0.86**Summary:** The isolated individual or tandem bromodomains of many BET family members, including BRD2, BRD3, BRD4, and BRDT, bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters. Small molecule inhibitors of bromodomain interactions hold promise as useful therapeutics for human disease. Cayman's BRD3 bromodomain 1 TR-FRET Assay Kit is a homogeneous, TR-FRET assay method amenable to rapid characterization of inhibitors of bromodomain/acetylated peptide interaction in a high throughput format.384 wells
1,920 wells
9,600 wells


BRD4 bromodomain 1 (human recombinant; GST tagged) 11068

*Bromodomain containing 4, HUNK1, MCAP***M_r:** 41.4 kDa **Purity:** ≥60% **Stability:** ≥6 months at -80°C**Source:** Recombinant N-terminal GST-tagged protein consisting of amino acids 49-170 expressed in *E. coli* • The isolated individual or tandem bromodomains of BRD4 have been shown to bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters.25 µg
50 µg
100 µg


BRD4 bromodomain 1 (human recombinant; His-tagged) 11720

*Bromodomain containing 4, HUNK1, MCAP***M_r:** 16.6 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C**Source:** Recombinant N-terminal His-tagged protein consisting of amino acids 49-170 expressed in *E. coli* • The isolated individual or tandem bromodomains of BRD4 have been shown to bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters.25 µg
50 µg
100 µg


BRD4 bromodomain 1 TR-FRET Assay Kit 600520

*Bromodomain containing 4, HUNK1, MCAP***Stability:** ≥6 months at -80°C **Z' Factor:** 0.73**Summary:** The isolated individual or tandem bromodomains of BRD4 bind acetylated histone tails, which couples histone acetylation marks to the transcriptional regulation of target promoters. Small molecule inhibitors of bromodomain interactions hold promise as useful therapeutics for human disease. Cayman's BRD4 bromodomain 1 TR-FRET Assay Kit is a homogeneous, TR-FRET assay method amenable to rapid characterization of inhibitors of bromodomain/acetylated peptide interaction in a high throughput format.384 wells
1,920 wells
9,600 wells


BRD4 bromodomain 1 and 2 (human recombinant) 11052


*Bromodomain containing 4, HUNK1, MCAP***M_r:** 73.4 kDa **Purity:** ≥90% **Stability:** ≥6 months at -80°C**Source:** Recombinant GST-tagged protein consisting of amino acids 49-460 expressed in *E. coli* • The isolated individual or tandem bromodomains of BRD4 have been shown to bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters.25 µg
50 µg
100 µg

BRD4 bromodomain 2 (human recombinant; GST-tagged) 11066










*Bromodomain containing 4, HUNK1, MCAP***M_r:** 40.6 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C**Source:** Recombinant N-terminal GST-tagged protein consisting of amino acids 342-460 expressed in *E. coli* • The isolated individual or tandem bromodomains of BRD4 have been shown to bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters.25 µg
50 µg
100 µg

BRD4 bromodomain 2 (human recombinant; His-tagged) 11721

*Bromodomain containing 4, HUNK1, MCAP***M_r:** 15.8 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C**Source:** Recombinant N-terminal His-tagged protein consisting of amino acids 342-460 expressed in *E. coli* • The isolated individual or tandem bromodomains of BRD4 have been shown to bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters.25 µg
50 µg
100 µg


BRD4 bromodomain 2 TR-FRET Assay Kit 600530*Bromodomain containing 4, HUNK1, MCAP***Stability:** ≥6 months at -80°C **Z' Factor:** 0.68**Summary:** The isolated individual or tandem bromodomains of BRD4 bind acetylated histone tails, which couples histone acetylation marks to the transcriptional regulation of target promoters. Small molecule inhibitors of bromodomain interactions hold promise as useful therapeutics for human disease. Cayman's BRD4 Bromodomain 2 TR-FRET Assay Kit is a homogeneous, TR-FRET assay method amenable to rapid characterization of inhibitors of bromodomain/acetylated peptide interaction in a high throughput format.384 wells
1,920 wells
9,600 wells**BRD4/HUNK1 Polyclonal Antibody** 13497*Bromodomain containing 4*Protein G-purified IgG **Stability:** ≥1 year at -80°C**Summary:** Antigen: peptide from human BRD4 within the region of amino acids 150-200 • Host: rabbit • Cross Reactivity:(+) chimpanzee, human, mouse, and rat BRD4 • Application(s): WB • BRD4 is a chromatin-binding protein whose expression is induced in response to growth stimuli. It acts at different stages of the cell cycle to regulate cell growth and chromosomal dynamics during mitosis.

1 ea

**BRD9 bromodomain (human recombinant)** 11509*Bromodomain containing 9***M_r:** 40.7 kDa **Purity:** ≥95% **Stability:** ≥1 year at -80°C**Source:** Recombinant N-terminal GST-tagged protein consisting of amino acids 21-137 expressed in *E. coli* • Human BRD9 contains a single bromodomain and has five isoforms that are produced by alternative splicing.25 µg
50 µg
100 µg**BRDT bromodomain 1 (human recombinant)** 11548*BRD6, BRD-containing protein testis specific, Cancer/testis antigen 9, CT9, RING3-like protein***M_r:** 40.6 kDa **Purity:** ≥95% **Stability:** ≥1 year at -80°C**Source:** Recombinant N-terminal GST-tagged protein consisting of amino acids 21-137 expressed in *E. coli* • BRDT is similar to the RING3 protein family and possesses 2 bromodomain motifs and a PEST sequence motif, which is a region rich in proline (P), glutamic acid (E), serine (S), and threonine (T) residues known to have a short intracellular half-life.25 µg
50 µg
100 µg**BRDT bromodomain 1 TR-FRET Assay Kit** 600650*BRD6, BRD-containing protein testis specific, Cancer/testis antigen 9, CT9, RING3-like protein***Stability:** ≥6 months at -80°C **Z' Factor:** 0.82**Summary:** The isolated individual or tandem bromodomains of many BET family members, including BRD2, BRD3, BRD4, and BRDT, bind acetylated histone tails, serving to couple histone acetylation marks to the transcriptional regulation of target promoters. Small molecule inhibitors of bromodomain interactions hold promise as useful therapeutics for human disease. Cayman's BRDT Bromodomain 1 TR-FRET Assay Kit is a homogeneous, TR-FRET assay method amenable to rapid characterization of inhibitors of bromodomain/acetylated peptide interaction in a high throughput format.384 wells
1,920 wells
9,600 wells**BRDT bromodomain 2 (human recombinant)** 11649*BRD6, Bromodomain testis-specific protein, Cancer/testis antigen 9, CT9, RING3-like protein***M_r:** 41.2 kDa **Purity:** ≥95% **Stability:** ≥1 year at -80°C**Source:** Recombinant N-terminal GST-tagged protein consisting of amino acids 259-379 expressed in *E. coli* • BRDT shares homology with the RING3 protein. The two bromodomains of BRDT recognize acetylated histone H4. Loss of BRDT leads to defects in spermatogenesis. In addition to testis specific expression, BRDT was found in approximately 20% of non-small cell lung cancers.25 µg
50 µg
100 µg**BRG1 bromodomain (human recombinant)** 11284*ATP-dependent helicase SMARCA4, BAF190A Mitotic growth and transcription activator, BRG1-associated factor 190A, Protein BRG-1***M_r:** 41.8 kDa **Purity:** ≥95% **Stability:** ≥1 year at -80°C**Source:** Recombinant N-terminal GST-tagged protein consisting of amino acids 1,448-1,575 expressed in *E. coli* • BRG1 is a member of the SWI/SNF protein family, which forms part of a large ATP-dependent chromatin remodeling complex. BRG1 is mutated in many cancer cell lines, such as breast, prostate, lung, pancreas and colon. Further, BRG1 has an important role as a tumor suppressor.25 µg
50 µg
100 µg**BRM bromodomain (human recombinant)** 11289*SMARCA2A/B, SWI/SNF ATPase***M_r:** 43.7 kDa **Purity:** ≥95% **Stability:** ≥1 year at -80°C**Source:** Recombinant N-terminal GST-tagged protein consisting of amino acids 1,367-1,511 expressed in *E. coli* • BRM is a member of the SWI/SNF protein family, which forms part of a large ATP-dependent chromatin remodeling complex. This complex is required for transcriptional activation of genes normally repressed by chromatin.25 µg
50 µg
100 µg**BRM bromodomain TR-FRET Assay Kit** 600730**Stability:** ≥6 months at -80°C **Z' Factor:** 0.62**Summary:** The SWI/SNF family of ATP-dependent chromatin remodeling enzymes catalyzes structural changes that allow proteins to access the nucleosomal DNA, alter the position of the nucleosomes on the DNA, or eject the histone octamer from the template. Mammalian SWI/SNF remodelers contain either BRM (SMARCA2) or BRG1 (SMARCA4) as their catalytic subunit. These subunits also contain a C-terminal bromodomain that binds to acetylated residues on histone tails. Small molecule inhibitors of BRM interactions with peptide binding partners hold promise as useful therapeutics for human disease. Cayman's BRM bromodomain TR-FRET Assay Kit is a homogeneous, TR-FRET assay method amenable to rapid characterization of inhibitors of bromodomain/acetylated peptide interaction in a high throughput format.384 wells
1,920 wells
9,600 wells**BRPF3 bromodomain 1 (human recombinant)** 11650*Bromodomain and PHD Finger Containing 3, KIAA1286, MGC 58603***M_r:** 41.2 kDa **Purity:** ≥95% **Stability:** ≥1 year at -80°C**Source:** Recombinant N-terminal GST-tagged protein consisting of amino acids 576-701 expressed in *E. coli* • BRPF3 is a component of the MOZ/MORF HAT complex. The addition of BRPF proteins to MOZ/MORF increases its HAT activity. Consequently, BRPF3 is likely to play a role in regulation of transcriptional activation by MOZ/MORF.25 µg
50 µg
100 µg


CECR2 bromodomain (human recombinant)

14133

*Cat Eye Syndrome Critical Region Protein 2***M_r**: 40.6 kDa **Purity**: ≥95% **Stability**: ≥1 year at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 424-538 expressed in *E. coli* • CECR2 is a transcription factor that forms a heterodimeric complex with CECR2-containing-remodeling factor (CERF). The CECR2/CERF forms a complex with the ATP-dependent chromatin remodeler SNF2L playing a critical role in neurulation. More recently, the bromodomain of CECR2 was shown to have strong γ-H2AX inhibition activity suggesting that CECR2 may play a role in DNA damage response.

CREB-binding protein bromodomain (human recombinant)


11288

*cAMP-responsive element-binding protein 1 CREB-1, CBP, CREBBP***M_r**: 40.8 kDa **Purity**: ≥95% **Stability**: ≥1 year at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 1,081-1,197 expressed in *E. coli* • CREBBP bromodomain has been shown to modulate the stability and function of the tumor suppressor protein p53. CREBBP bromodomain recognizes the acetylated lysine residue 382 on p53.

(+) -JQ1

11187

[1268524-70-4]


MF: C₂₃H₂₅ClN₄O₂S **FW**: 457.0 **Purity**: ≥98%A crystalline solid **Stability**: ≥2 years at -20°C**Summary**: Displaces BET proteins from chromatin by competitively binding to the acetyl-lysine recognition pocket of BET bromodomains; binds BRD4 bromodomains 1 and 2 with K_d values of ~50 and 90 nM, respectively1 mg
5 mg
10 mg

NOTE: Manufactured, marketed, and sold with authorization from Tensha Therapeutics, Inc. Patent Pending relating to PCT Publ. No. WO/2011/143669, and any related U.S. and foreign patents and patent applications.

(-)-JQ1

11232


[1268524-71-5]

MF: C₂₃H₂₅ClN₄O₂S **FW**: 457.0 **Purity**: ≥98%A crystalline solid **Stability**: ≥2 years at -20°C**Summary**: The inactive stereoisomer of a selective BET bromodomain inhibitor1 mg
5 mg
10 mg

NOTE: Manufactured, marketed, and sold with authorization from Tensha Therapeutics, Inc. Patent Pending relating to PCT Publ. No. WO/2011/143669, and any related U.S. and foreign patents and patent applications.


PCAF bromodomain (human recombinant)

11920

*KAT2B, p300/CBP-associated factor***M_r**: 40.9 kDa **Purity**: ≥95% **Stability**: ≥6 months at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 714-831 expressed in *E. coli* • PCAF is a transcriptional coactivator that works both as a histone lysine acetyltransferase, through its HAT domain, and as an acetyl-lysine reader through its conserved bromodomain located directly C-terminal to the HAT domain. The PCAF bromodomain binds acetylated histone H3 and H4 as well as non-histone targets. Bromodomain binding is dictated by the position of the acetylated lysine as well as interactions with specific residues flanking the acetyl-lysine.25 µg
50 µg
100 µg


PFI-1

11155

MF: C₁₆H₁₇N₃O₄S **FW**: 347.4 **Purity**: ≥98%A crystalline solid **Stability**: ≥2 years at -20°C**Summary**: A BET bromodomain inhibitor that exhibits inhibitory activity at BRD2 bromodomain 2 and BRD4 bromodomain 1 with IC₅₀ values of 98 nM and 0.22 µM, respectively1 mg
5 mg
10 mg
25 mg


Polybromo-1D bromodomain 1 (human recombinant)

11652

*BAF180, BRG1-associated factor 180, hPBI, MGC165155, MGC165156, PB1, PBRM1, Protein Polybromo-1***M_r**: 42.8 kDa **Purity**: ≥95% **Stability**: ≥1 year at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 23-156 expressed in *E. coli* • PBRM1 contains six bromodomains and is a component of the SWI/SNF complex. PBAF. PBAF is targeted to acetylated sites in chromatin by the PBRM1 bromodomains, where it plays a role in cell cycle regulation and tumor suppression.


TAF1 bromodomain 1 (human recombinant)

11922

*BA2R, CCG1, CCGS, DYT3, KAT4, NSCL2, TAF1 RNA Polymerase II, TAFII250, TAF2A, TATA Box Binding Protein (TBP)-Associated Factor, XDP***M_r**: 41.7 kDa **Purity**: ≥95% **Stability**: ≥1 year at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 1,371-1,496 expressed in *E. coli* • TAF1 is a component of transcription factor IID, and binds to core promoter sequences at the transcription start site. TAF1 helps control transcription by both its kinase and histone acetyltransferase enzymatic activities. It interacts with transcriptional activators such as the androgen receptor to promote transcription.25 µg
50 µg
100 µg


TRIM24 bromodomain (human recombinant)

11653

*E3 Ubiquitin-protein Ligase TRIM24, bTIF1, PTC6, RING finger protein 82, RNF82, Transcriptional Intermediary Factor 1-α, Tripartite Motif Containing 24***M_r**: 40.9 kDa **Purity**: ≥95% **Stability**: ≥1 year at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 896-1,014 expressed in *E. coli* • TRIM24 is a transcriptional cofactor, whose inactivation leads to hepatocellular carcinoma in mice. The N-terminal TRIM domain of TRIM24 binds ligand-bound nuclear receptors, while its tandem C-terminal plant homeo-domain and bromodomain target TRIM24 to acetylated histones in chromatin.25 µg
50 µg
100 µg

WDR9 bromodomain 2 (human recombinant)

11549

*BRD and WD repeat-containing protein 1, WD repeat-containing protein 9***M_r**: 40.3 kDa **Purity**: ≥85% **Stability**: ≥1 year at -80°C**Source**: Recombinant N-terminal GST-tagged protein consisting of amino acids 1,310-1,426 expressed in *E. coli* • WDR9 possesses two bromodomain motifs and eight WD repeats. WDR9 is also known to interact with BRG1 (SMARCA4). This product contains the bromodomain 2 region of WDR9.25 µg
50 µg
100 µg

NOTE: SUMO Affinity tag used under license from LifeSensio

Cayman Chemical is helping make research possible with

Epigenetic Screening and Lead Discovery **SERVICES**

Learn more

For more information or to receive a quote for contract services please contact our sales department at

1-800-364-9897

or email us at

sales@caymanchem.com



Epigenetics Screening Library (96-Well) 11076

10 mM solutions in DMSO Stability: ≥ 2 years at -20°C

Summary: The Epigenetics Screening Library contains various molecules that are known to modulate the activity of a variety of epigenetic 'writers and erasers' and 'reader' proteins in a 96-well Matrix tube rack format as 10 mM stocks in DMSO.

It may include compounds that modulate the activity of methyltransferases, demethylases, histone acetyltransferases, histone deacetylases, and acetylated histone binding proteins. The composition of this screening library will always vary somewhat depending upon our inventory. Stability data is not available for compounds as supplied in the screening library.

50 μl • 100 μl • 200 μl



Screening and Profiling Services

Cayman offers a dedicated Epigenetic Screening Laboratory designed to be flexible and innovative. High-throughput capabilities allow us to work with you to screen a chemical library against specific epigenetic targets. Alternatively, our broad collection of epigenetic enzymes, substrates, and assays enable profiling the activity of a few compounds against several targets. Our experienced staff and expanding suite of assays are designed to get the results you need in a timely manner.

- **Comprehensive screening laboratory focused on epigenetics**
- **Backed by Cayman's core strengths in protein production, assay development, chemical synthesis, and medicinal chemistry**
- **Capacity to test compounds at a rate of up to 100,000 compounds per week**
- **Ability to profile compounds against a broad epigenetic enzyme panel**

Includes assays for:

methyltransferases

demethylases

acetyltransferases

**deacetylases
(HDACs and SIRT5)**

**histone code readers
(e.g. bromodomain proteins)**

Gene-to-Structure Services

Cayman Chemical has been making recombinant proteins and developing assays for many years. As an extension of these core competencies, we now offer full gene-to-structure contract services. Our experienced lab staff will clone, express, purify, and determine the structure of your protein in the presence or absence of ligands. We are fully equipped to aid you in the drug discovery process.

- **Custom and large scale protein production**
- **Macromolecular X-ray crystallography**
- **Biophysical characterization**
- **Structure Based Drug Design (SBDD)**



Methyl Readers

Antibodies

- 36 13771 MBD1 Monoclonal Antibody (Clone 100B272.1)
 36 13772 MBD1 Polyclonal Antibody
 36 13777 MBD2 Binding Zinc Finger Polyclonal Antibody
 37 13773 MBD2/3 Monoclonal Antibody (Clone 106B691)
 37 13775 MeCP2 Polyclonal Antibody

Biochemicals

- 37 13968 UNC1215

Proteins

- 36 11235 HP1- α (human recombinant)
 36 11286 MBD2 (human recombinant; methyl binding domain aa 150-220)
 37 11287 MeCP2 (human recombinant; methyl binding domain aa 77-166)


HP1- α (human recombinant)

11235

Antigen p25, CBX5, Chromobox Homolog 5, Heterochromatin Protein 1- α **M_r**: 24.1 kDa **Purity**: $\geq 95\%$ **Stability**: ≥ 6 months at -80°C

Source: Recombinant N-terminal His-tagged protein consisting of amino acids 2-191 expressed in *E. coli* • HP1- α is involved in gene regulation and heterochromatin formation. The chromodomain of HP1- α has been shown to recognize di- and trimethylated histone H3K9. The methyltransferase SETDB1 can then be recruited which performs H3K9 trimethylation, leading to chromatin condensation and gene silencing. The chromoshadow domain facilitates protein-protein interactions and is responsible for recruitment to sites of DNA damage, where HP1- α helps to reorganize chromatin as part of the DNA damage response system.

25 μg
 50 μg
 100 μg




MBD1 Monoclonal Antibody (Clone 100B272.1)

13771

*Methyl-CpG-Binding Domain 1*Protein G-purified IgG **Stability**: ≥ 1 year at -20°C

Summary: Antigen: human MBD1 amino acids 391-405 • Host: mouse, clone 100B272.1 • Isotype: IgG₁ • Cross Reactivity: (+) human MBD1 • Application(s): WB • MBD1 contains an MBD that allows it to bind specifically to methylated DNA and to repress transcription from methylated gene promoters.

1 ea




MBD1 Polyclonal Antibody

13772

*Methyl-CpG-Binding Domain 1*Protein G-purified IgG **Stability**: ≥ 1 year at -20°C

Summary: Antigen: human MBD1 amino acids 98-113 and 391-405 • Host: rabbit • Cross Reactivity: (+) human MBD1 • Application(s): WB • MBD1 contains an MBD that allows it to bind specifically to methylated DNA and to repress transcription from methylated gene promoters.

1 ea




MBD2 (human recombinant; methyl binding domain aa 150-220)

11286

*Methyl-CpG Binding Domain 2, Methyl Cytosine Binding Domain 2***M_r**: 34.8 kDa **Purity**: $\geq 95\%$ **Stability**: ≥ 1 year at -80°C

Source: Recombinant N-terminal GST-tagged protein consisting of amino acids 150-220 expressed in *E. coli* • MBD2 specifically binds to methylated promoters on CpG islands. MBD2 binding to 5mC facilitates the recruitment of chromatin remodeling and transcriptional repressor complexes, which results in a repressive chromatin state.

25 μg
 50 μg
 100 μg




MBD2 Binding Zinc Finger Polyclonal Antibody

13777

*HINFP, Histone H4 Transcription Factor, Methyl-CpG-Binding Domain 2, MIZF*Protein G-purified IgG **Stability**: ≥ 1 year at -20°C

Summary: Antigen: human MIZF amino acids 180-194, 331-346, and 371-388 • Host: rabbit • Cross Reactivity: (+) human MIZF • Application(s): WB • MIZF protein represses transcription by associating with MBD2 in a histone deacetylase complex.

1 ea




MBD2/3 Molyclonal Antibody (Clone 106B691)

13773

*Methyl-CpG-Binding Domain 2/3*Protein G-purified IgG **Stability**: ≥ 1 year at -20°C

Summary: Antigen: human MBD3 amino acids 215-230 • Host: mouse, clone 106B691 • Isotype: IgG_{1k} • Cross-reactivity: (+) human MBD2/3 • Application(s): WB • MBD2 and MBD3 are members of a family of nuclear proteins related by the presence in each of an MBD. MBD2 is capable of binding specifically to methylated DNA, whereas MBD3 cannot either *in vitro* or *in vivo*.

1 ea




MeCP2 (human recombinant; methyl binding domain aa 77-166)

11287

*Methyl-CpG Binding Protein 2***M_r**: 37 kDa **Purity**: $\geq 95\%$ **Stability**: ≥ 1 year at -80°C

Source: Recombinant N-terminal GST-tagged protein consisting of amino acids 77-166 expressed in *E. coli* • MeCP2 specifically binds to methylated promoters on CpG islands and mediates gene silencing by recruiting corepressor complexes. *In vitro* work suggests high affinity binding of MeCP2 is facilitated by DNA fragments containing A/T bases adjacent to the MeCpG.

25 μg
 50 μg
 100 μg




MeCP2 Polyclonal Antibody

13775

*Methyl-CpG-Binding Protein 2*Protein G-purified IgG **Stability**: ≥ 1 year at -20°C

Summary: Antigen: human MeCP2 amino acids 11-25 and 181-195 • Host: rabbit • Cross Reactivity: (+) human MeCP2 • Application(s): WB • MeCP2 may function as a mediator of the biological consequences of the methylation signal. It is also reported that this protein functions as a demethylase to activate transcription, as DNA methylation causes gene silencing.

1 ea




UNC1215

13968

*[1415800-43-9]***MF**: C₃₂H₄₃N₅O₂ **FW**: 529.7 **Purity**: $\geq 98\%$ A crystalline solid **Stability**: ≥ 2 years at -20°C

Summary: A potent, selective chemical probe for the methyl lysine reading function of L3MBTL3 (K_d = 120 nM; IC₅₀ = 40 nM) that competitively displaces mono- or dimethyl-lysine containing peptides

1 mg
 5 mg
 10 mg
 50 mg



Methyltransferases

Antibodies

- 40 13536 DNA Methyltransferase 1-Associated Protein 1 Polyclonal Antibody
40 13479 DNA Methyltransferase 1 Monoclonal Antibody (Clone 60B1220.1)
40 13481 DNA Methyltransferase 2 Monoclonal Antibody (Clone 102B1259.2)
40 13480 DNA Methyltransferase 2 Polyclonal Antibody
40 13483 DNA Methyltransferase 3a Monoclonal Antibody - Biotinylated (Clone 64B814.1)
40 13484 DNA Methyltransferase 3a Monoclonal Antibody (Clone 64B1446)
40 13482 DNA Methyltransferase 3a Monoclonal Antibody (Clone 64B814.1)
40 13485 DNA Methyltransferase 3b Monoclonal Antibody (Clone 52A1018)
41 13487 EZH1 Polyclonal Antibody
42 13782 I2PP2A/SET Polyclonal Antibody
42 13727 Methylated Lysine Polyclonal Antibody
43 13728 Methylated Lysine Polyclonal Antibody-biotin
43 13729 Methylated Lysine Polyclonal Antibody HRP Conjugate
45 13552 PRMT4 Polyclonal Antibody
46 13559 PRMT5 Polyclonal Antibody
46 13558 PRMT6 Polyclonal Antibody
46 13551 PRMT7 Polyclonal Antibody
47 13731 SET7/9 Polyclonal Antibody
50 13780 SET7/9 (FL) Polyclonal Antibody

Biochemicals

- 38 13956 S-(5'-Adenosyl)-L-methionine chloride (hydrochloride)
38 13965 AMI-1 (sodium salt)
39 11164 5-Azacytidine
39 13373 2',3',5'-triacetyl-5-Azacytidine
39 13124 BIX01294 (hydrochloride hydrate)
39 13156 Chaetocin
39 13828 3-Deazaneplanocin A
39 11102 3-Deazaneplanocin A (hydrochloride)
39 11166 Decitabine
41 10569 Ellagic Acid
44 11620 MI-2 (hydrochloride)
44 11621 MI-nc (hydrochloride)
47 13302 RG-108
52 13631 UNC0224
52 10582 UNC0321 (trifluoroacetate salt)
52 10734 UNC0638

Kits

- 40 589324 DNA Methylation EIA Kit
42 700500 G9a Methyltransferase Inhibitor Screening Assay Kit
42 600570 GLP SAM-Screener™ Assay Kit
43 700140 Methyltransferase Colorimetric Assay Kit
43 700150 Methyltransferase Fluorometric Assay Kit
43 600580 MLL1 SAM-Screener™ Assay Kit
47 700270 SET7/9 Methyltransferase Inhibitor Screening Assay Kit
50 600490 SET7/9 SAM-Screener™ Assay Kit
51 700350 SET8 Methyltransferase Inhibitor Screening Assay Kit

Proteins

- 39 10946 Ash2L (human recombinant)
40 10770 DNA Methyltransferase 3L (human recombinant)
41 10354 Dot1L (human recombinant)
41 11178 DPY-30 (human recombinant)
41 10628 EED (human recombinant)
41 10353 G9a (human recombinant)
42 10755 G9a-like protein (human recombinant)
44 10658 MLL1 (human recombinant)
44 10756 MLL1/WAR complex (human recombinant)
44 10945 MLL1/WAR complex (human recombinant)
44 10757 NSD1 (human recombinant)
45 10758 NSD2 (human recombinant)
45 11209 PRDM9 (human recombinant)
45 10350 PRMT1 (human recombinant)
45 11642 PRMT3 (human recombinant)
45 10750 PRMT4 (human recombinant)
46 10752 PRMT6 (human recombinant)
46 13866 PRMT6 (human recombinant; baculovirus expressed)
46 11644 PRMT8 (human recombinant)
47 10947 RbBP5 (human recombinant)
47 10765 RIZ1 (human recombinant)
47 10320 SET7/9 (human recombinant)
50 10319 SET8 (human recombinant)
50 10767 SETD2 (human recombinant)
51 10761 SMYD1 (human recombinant)
51 10762 SMYD3 (human recombinant)
51 10763 SUV4-20H1 (human recombinant)
51 10764 SUV4-20H2 (human recombinant)
52 10228 TAF 10 Peptide
52 10944 WDR5 (human recombinant)


S-(5'-Adenosyl)-L-methionine chloride (hydrochloride) 13956

[86867-01-8] AdoMet, SAM

MF: C₁₅H₂₃ClN₆O₅S • 2HCl FW: 507.8 Purity: ≥95%

A lyophilized powder Stability: ≥1 year at -80°C

Summary: A ubiquitous methyl donor involved in a wide variety of biological reactions, including those mediated by DNA and protein methyltransferases


1 mg
2.5 mg
5 mg

AMI-1 (sodium salt) 13965

Arginine N-Methyltransferase Inhibitor-1


MF: C₂₁H₁₂N₂O₉S₂ • 4Na FW: 592.4 Purity: ≥99%

A crystalline solid Stability: ≥2 years at -20°C

Summary: A cell permeable inhibitor of PRMTs; inhibits both yeast Hmt1p and human PRMT1 (IC₅₀ = 3.0 and 8.8 μM, respectively); also effectively blocks the activity of PRMTs 3, 4, and 6 but not that of lysine methyltransferases; inhibits HIV-1 reverse transcriptase (IC₅₀ = 5.0 μM)5 mg
25 mg

Ash2L (human recombinant) 10946

Absent, small, or homeotic discs 2-like, Set1/Ash2 Histone Methyltransferase Complex Subunit Ash2 Isoform A


M_r: 60.1 kDa Purity: ≥90% Stability: ≥6 months at -80°CSource: Recombinant protein consisting of amino acids 96-628 expressed in *E. coli* with a SUMO tag • ASH2L is a component of various multisubunit protein complexes, including the large complex of proteins associated with the SET1 (MLL) family of lysine methyltransferases.25 μg
50 μg
100 μgLane 1: MW Ladders
Lane 2: ASH2L (1 μg)
Lane 3: ASH2L (2 μg)
Lane 4: ASH2L (4 μg)
Lane 5: ASH2L (6 μg)

5-Azacytidine 11164

[320-67-2] Antibiotic U 18496, 5-AzaC, Ladakamycin, Mylosar, NSC 102816, NSC 103-627, WR 183027

MF: C₈H₁₂N₄O₅ FW: 244.2 Purity: ≥95%

A crystalline solid Stability: ≥2 years at -20°C

Summary: An inhibitor of DNA methyltransferases that reduces hypermethylation associated with certain diseases, including myelodysplastic syndromes (IC₅₀s = 2.4 and 2.6 μM for *in vitro* anti-myeloma activity) and cancer (IC₅₀s = ~ 0.4 μM for inhibiting proliferation of various cancer cell lines)50 mg
100 mg
250 mg


2',3',5'-triacetyl-5-Azacytidine 13373

[10302-78-0]

MF: C₁₄H₁₈N₄O₈ FW: 370.3 Purity: ≥95%

A crystalline solid Stability: ≥2 years at -20°C


Summary: A prodrug form of 5-azacytidine, an inhibitor of DNA methyltransferases, that may reverse epigenetic changes

5 mg
10 mg
50 mg
100 mg

BIX01294 (hydrochloride hydrate) 13124

MF: C₂₈H₃₈N₆O₂ • 3HCl [XH₂O] FW: 600.0 Purity: ≥98%

A crystalline solid Stability: ≥2 years at -20°C


Summary: A selective inhibitor of G9a histone methyltransferase (IC₅₀ = 1.7 μM); less effectively inhibits G9a-like protein (IC₅₀ = 38 μM) and has no effect on other known histone methyltransferases1 mg
5 mg
10 mg
50 mg

Chaetocin 13156

[28097-03-2]

MF: C₃₀H₂₈N₆O₆S₄ FW: 696.8 Purity: ≥95%

A crystalline solid Stability: ≥2 years at -20°C

Summary: A fungal mycotoxin that inhibits the Lys9-specific histone methyltransferases SU(VAR)3-9 (IC₅₀ = 0.8 μM), G9a (IC₅₀ = 2.5 μM), and DIM5 (IC₅₀ = 3 μM)1 mg
5 mg
10 mg


3-Deazaneplanocin A 13828

[102052-95-9] DZNeP, NSC 617989

MF: C₁₂H₁₄N₄O₃ FW: 262.3 Purity: ≥98%

A crystalline solid Stability: ≥2 years at -20°C

Summary: An inhibitor of SAH hydrolase; depletes EZH2 levels and inhibits trimethylation of lysine 27 on histone H3 in acute myeloid leukemia cells in a dose-dependent manner (0.2-1 μM); increases expression of the cell-cycle regulators p21, p27, and FBXO32 leading to cell cycle arrest and apoptosis

500 μg
1 mg
5 mg

*Also Available: 3-Deazaneplanocin A (hydrochloride) (11102)


Decitabine 11166

[2353-33-5] 5-aza-2'-Deoxycytidine, DAC, Dacogen, NSC 127716

MF: C₈H₁₂N₄O₄ FW: 228.2 Purity: ≥98%

A crystalline solid Stability: ≥2 years at -20°C

Summary: A 2' deoxy analog of 5-azacytidine which is incorporated into DNA and causes hypomethylation by inhibiting DNA methyltransferases in a concentration-dependent manner; useful in conditions characterized by DNA hypermethylation, as is found in myelodysplastic syndromes

5 mg
10 mg
25 mg
50 mg

DNA Methylation EIA Kit

589324


Stability: ≥1 year at -20°C
Sensitivity: 50% B/B₀; 12 ng/ml • 80% B/B₀; 3 ng/ml
Summary: DNA methylation is an important epigenetic process regulating gene expression. Methylation occurs on carbon 5 of 2'-deoxy-cytidine yielding the modified base 5-methyl-2'-deoxy cytidine. The methylation pattern of cells is tightly regulated during development with the methylation profile being transmitted from parent to daughter cells during cell division. Methylation results in long-term silencing of genes, while unmethylated regions of DNA can be actively transcribed. Global changes in methylation can be quantified by measuring plasma or urinary levels of 5-methyl-2'-deoxy cytidine. These changes in methylation can provide valuable information about cancer status of an individual. For example, patients with leukemia excrete significantly elevated levels of 5-methyl-2'-deoxy cytidine compared to healthy individuals. Global methylation within tissues can be measured in a similar manner, allowing study of tissue-specific changes that occur as a result of differentiation, aging, or carcinogenesis. Cayman's DNA Methylation EIA Kit is a competitive assay that can be used for the quantification of 5-methyl-2'-deoxy cytidine in urine, culture supernatants, plasma, and other sample matrices.

Specificity:

5-Methyl-2'-deoxycytidine	100%
5-Methylcytidine	20%
2-Deoxycytidine	0.1%
Cytidine	0.1%

For a full specificity profile, please go to www.caymanchem.com


96 solid wells
96 strip wells
480 solid wells
480 strip wells



DNA Methyltransferase 1-Associated Protein 1 Polyclonal Antibody

13536


DMAP1
Protein G-purified IgG **Stability:** ≥1 year at -20°C
Summary: Antigen: peptide from human DMAP1 within the region of amino acids 250-300 • **Host:** rabbit • **Cross Reactivity:** (+) chimpanzee, bovine, canine, human, mouse, and rat DMAP1 • **Application(s):** IHC (paraffin-embedded sections) and WB • DMAP1 is involved in the repression or activation of transcription. It is a component of the NuA4 histone acetyltransferase complex and interacts with the transcriptional corepressor tumor susceptibility gene 101 and the pro-apoptotic death-associated protein 6.



DNA Methyltransferase 3L (human recombinant)

10770

DNMT3L, DNMT3-Like Protein
M_r: 53.2 kDa **Purity:** ≥95% **Stability:** ≥9 months at -80°C
Source: Recombinant N-terminal GST-tagged protein consisting of amino acids 160-387 expressed in *E. coli* • DNMT3L is required to stimulate the DNA methylation activity of DNMT3a and 3b through interactions with the catalytic domain of DNMT3a and 3b.




DNMT Antibodies			
Item No.	Product Name	Application(s)	Species Specificity
13479	DNA Methyltransferase 1 Monoclonal Antibody (Clone 60B1220.1)	WB, IHC (paraffin-embedded sections), ChIP, and IP	(+) Human, mouse, and zebrafish DNMT1
13481	DNA Methyltransferase 2 Monoclonal Antibody (Clone 102B1259.2)	WB	(+) Human and mouse DNMT2
13480	DNA Methyltransferase 2 Polyclonal Antibody	WB	(+) Human and mouse DNMT2
13483	DNA Methyltransferase 3a Monoclonal Antibody - Biotinylated (Clone 64B814.1)	ELISA	(+) Human and mouse DNMT3a
13484	DNA Methyltransferase 3a Monoclonal Antibody (Clone 64B1446)	WB, IF/ICC, ChIP, and IHC (paraffin)	(+) Human and mouse DNMT3a
13482	DNA Methyltransferase 3a Monoclonal Antibody (Clone 64B814.1)	WB, IF, and ICC	(+) Human and mouse DNMT3a
13485	DNA Methyltransferase 3b Monoclonal Antibody (Clone 52A1018)	WB, IP, IF, ICC, ChIP, and IHC (paraffin)	(+) Human and mouse DNMT3b

Dot1L (human recombinant)

10354


Disruptor of Telomeric Signaling 1-Like, Dot1-like
M_r: 74.1 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C
Source: Active recombinant N-terminal GST-tagged protein consisting of amino acids 2-416 expressed in *E. coli* • Dot1L is a non-SET domain containing methyltransferase that is the only enzyme known to methylate histone 3 at lysine 79, where it catalyzes mono-, di-, and trimethylation. Proper Dot1L function is necessary for transcriptional activation of many genes, DNA damage repair, and cell cycle regulation.



DPY-30 (human recombinant)

11178

DPY-30-Like protein, hDPY-30, SAF19
M_r: 11.2 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C
Source: Recombinant protein consisting of amino acids 2-99 expressed in *E. coli* • DPY-30 is a component of the MLL1 methylation complex, which interacts directly with Ash2L.




NOTE: SUMO Affinity tag used under license from LifeSensors, Inc.

EED (human recombinant)

10628

Embryonic Ectoderm Development
M_r: 53.5 kDa **Purity:** ≥95% **Stability:** ≥9 months at -80°C
Source: Recombinant N-terminal His-tagged protein consisting of amino acids 1-441 expressed in Sf21 cells • EED is a WD40 repeat-containing protein that forms part of the PRC2. The EED subunit does not contain methyltransferase activity. However, transcriptional repression by PRC2-mediated trimethylation of lysine 27 on Histone H3 has been shown to be dependent on EED binding to repressive histone marks.




Ellagic Acid

10569

[476-66-4] Gallogen, Lagistase, TBBD
MF: C₁₄H₆O₈ **FW:** 302.2 **Purity:** ≥95%
A crystalline solid **Stability:** ≥2 years at -20°C
Summary: A polyphenolic antioxidant that is abundant in many fruits, vegetables, plant bark, and peels; has anti-carcinogenic, anti-mutagenic, anti-inflammatory, and organ-preserving properties; blocks methylation of H3R17 by CARM1 without significantly altering histone acetylase or DNA methyltransferase activity


100 mg
500 mg
1 g



EZH1 Polyclonal Antibody

13487

Protein G-purified IgG **Stability:** ≥1 year at -20°C
Summary: Antigen: peptides of human EZH1 • **Host:** rabbit • **Cross Reactivity:** (+) human and mouse EZH1 • **Application(s):** WB • EZH1 is a human homolog of the *Drosophila* gene-enhancer of zeste, a member of the polycomb group of transcriptional repressors. It has a potential role in human development as a transcriptional regulator and a component of protein complexes that stably maintain heterochromatin.





G9a (human recombinant)

10353

EHMT2, Euchromatic Histone-Lysine N-Methyltransferase 2, KMT1C
M_r: 75.4 kDa **Purity:** ≥70% **Stability:** ≥6 months at -80°C
Source: Active recombinant N-terminal GST-tagged protein consisting of amino acids 785-1,210 expressed in *E. coli* • G9a is a SET domain-containing methyltransferase that specifically mono- and di-methylates Histone H3K9 at lysine 9.

25 µg
50 µg
100 µg




G9a-like protein (human recombinant) 10755*EHMT1, Euchromatic Histone-Lysine N-Methyltransferase 1, Eu HMTase 1, GLP, KMT1D***M_r:** 60.3 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal GST-tagged protein consisting of amino acids 1,004-1,298 expressed in *E. coli* • GLP is a SET domain-containing methyltransferase that specifically mono- and di-methylates H3K9.25 µg
50 µg
100 µgLane 1: MW Markers
Lane 2: GLP (1 µg)
Lane 3: GLP (2 µg)
Lane 4: GLP (5 µg)**G9a Methyltransferase Inhibitor Screening Assay Kit**


700500

*EHMT2***Stability:** ≥6 months at -80°C**Summary:** G9a is a SET domain-containing mammalian histone methyltransferase that can mono- or dimethylate lysine 9 and lysine 27 on histone H3. G9a is overexpressed in various cancers and is a potential inhibitory target for cancer treatment. In Cayman's G9a Methyltransferase Inhibitor Screening Assay, the transfer of the methyl group from SAM to the acceptor peptide by G9a generates SAH, which is rapidly converted to urate and H₂O₂ using an enzyme mixture provided in the kit. A subsequent reaction between H₂O₂ and ADHP produces the highly fluorescent compound resorufin.

96 wells

**GLP SAM-Screener™ Assay Kit**


600570

*G9a-like protein***Stability:** ≥6 months at -80°C **Z' Factor:** 0.56**Summary:** G9a-like protein is a SET domain-containing methyltransferase that specifically mono- and di-methylates H3K9. GLP and G9a function as major euchromatic H3K9me1 and H3K9me2 histone methyltransferases and also have been found to methylate several nonhistone substrates, including p53(K372). This fluorescence polarization assay is based upon a proprietary small molecule fluorescent probe that binds to the SAM binding pocket in GLP. Binding of the small molecule probe to GLP induces an increase in fluorescence polarization. Binding of the probe can be competed with the endogenous cofactor SAM or by the inhibitor sinefungin, but is unaffected by the histone H3 peptide substrate. The GLP SAM-Screener Assay is robust and exhibits a greater than 80 mP shift over a range of 0-250 nM GLP. The assay is suitable for high-throughput screening in the provided 384-well plate or can be scaled to higher density plate formats (e.g., 1,536-well) if desired.384 wells
1,920 wells**I2PP2A/SET Polyclonal Antibody**

13782

*Inhibitor of granzyme A-activated DNase, PHAPII, Phosphatase 2A Inhibitor, Template-Activating Factor 1*Protein G-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: synthetic peptides from human I2PP2A/SET amino acids 79-94 and 148-164 • Host: rabbit • Cross Reactivity: (+) human I2PP2A/SET • Application(s): WB • I2PP2A/SET is a multitasking protein, involved in apoptosis, transcription, nucleosome assembly, and histone binding.

1 ea

**Methylated Lysine Polyclonal Antibody**

13727

Affinity-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: methylated KLH • Host: rabbit • Cross Reactivity: (+) mono- and di-methylated lysine residues; (-) acetylated lysine • Application(s): ELISA, IHC, IP, and WB • Lysine can be methylated once, twice, or three times by lysine methyltransferases. The transfer of methyl groups from SAM to histones is catalyzed by histone methyltransferases.

400 µl

Methylated Lysine Polyclonal Antibody-biotin 13728Affinity-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: methylated KLH • Host: rabbit • Cross Reactivity: (+) methylated lysine residues • Application(s): ELISA, IP, and WB • Lysine can be methylated once, twice, or three times by lysine methyltransferases. The transfer of methyl groups from SAM to histones is catalyzed by histone methyltransferases.

400 µl

Methylated Lysine Polyclonal Antibody HRP Conjugate

13729

Affinity-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: methylated KLH • Host: rabbit • Cross Reactivity: (+) methylated lysine residues • Application(s): ELISA and WB • Lysine can be methylated once, twice, or three times by lysine methyltransferases. The transfer of methyl groups from SAM to histones is catalyzed by histone methyltransferases.


400 µl

Methyltransferase Colorimetric Assay Kit


700140

*MT***Stability:** ≥6 months at -80°C**Summary:** Cayman's Methyltransferase Colorimetric Assay Kit is a continuous enzyme-coupled assay that can continuously monitor SAM-dependent methyltransferases. The removal of the methyl group from SAM generates AdoHcy, which is rapidly converted to urate and H₂O₂ by an enzyme mixture provided in the kit. The rate of production of H₂O₂ is measured with the colorimetric reagent, 3,5-dichloro-2-hydroxybenzenesulfonic acid, by an increase in absorbance at 500-520 nm. The assay is supplied with AdoHcy as a positive control. The assay can be used with any purified SAM-dependent methyltransferase.


96 wells

**Methyltransferase Fluorometric Assay Kit** 700150*MT***Stability:** ≥6 months at -80°C**Summary:** Cayman's MT Fluorometric Assay is a continuous enzyme-coupled assay that can continuously monitor SAM-dependent MTs. The removal of the methyl group from SAM generates AdoHcy, which is rapidly converted to urate and H₂O₂ by an enzyme mixture provided in the kit. The reaction between H₂O₂ and ADHP produces the highly fluorescent compound resorufin, which is analyzed with an excitation wavelength of 530-540 nm and an emission wavelength of 585-595 nm. The assay is supplied with AdoHcy as a positive control. The assay can be used with any purified SAM-dependent MT.

96 wells


**MLL1 SAM-Screener™ Assay Kit**

600580

*Mixed-Lineage Leukemia-1***Stability:** ≥6 months at -80°C **Z' Factor:** 0.71**Summary:** MLL1 is a member of the trithorax group (trxG)/Set1-like family of gene activators that contains histone methyltransferase activity specific for lysine 4 of histone H3. This methylation plays an important role in gene activation at various developmentally regulated loci, such as the Hox gene loci. This fluorescence polarization assay is based upon a proprietary small molecule fluorescent probe that binds to the SAM binding pocket in MLL1. Binding of the small molecule probe to MLL1 induces an increase in fluorescence polarization. Binding of the probe can be competed with the endogenous cofactor SAM or by the inhibitor sinefungin, but is unaffected by the histone H3 peptide substrate. The MLL1 SAM-Screener Assay is robust and exhibits a greater than 100 mP shift over a range of 0-500 nM MLL1. The assay is suitable for high-throughput screening in the provided 384-well plate or can be scaled to higher density plate formats (e.g., 1,536-well) if desired.384 wells
1,920 wells

MI-2 (hydrochloride)


11620

MF: C₁₈H₂₅N₅S₂ • 2HCl **FW:** 448.5 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** Potently binds menin, blocks the menin-MLL fusion protein interaction (IC₅₀ = 0.45 μM), and induces apoptosis in cells expressing MLL fusion proteins1 mg
5 mg
10 mg
25 mg

• Also Available: MI-nc (hydrochloride) (11621)

MLL1 (human recombinant)


10658

*ALL1, HRX, KMT2A, Lysine Methyltransferase 2A, Mixed Lineage Leukemia-1***M_r:** 24.1 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C**Source:** Recombinant protein consisting of amino acids 3,762-3,969 expressed in *E. coli* • MLL1 plays a major role in epigenetic regulation through methylation of H3K4 to activate gene transcription.50 μg
100 μg
250 μgLane 1: MW Markers
Lane 2: MLL1 (4 μg)
Lane 3: MLL1 (2 μg)
Lane 4: MLL1 (1 μg)

NOTE: SUMO Affinity tag used under license from LifeSensors, Inc.

MLL1/WAR complex (human recombinant)


10756

*Mixed-Lineage Leukemia-1, MLL Core Complex***M_r:** 180 kDa **Purity:** ≥90% **Stability:** ≥6 months at -80°C**Source:** Active recombinant proteins expressed in *E. coli*. MLL1, WDR5, Ash2L, and RbBP5 are expressed with N-terminal His-SUMOpro affinity tags. The His-SUMO affinity tags were removed using recombinant SUMO Protease 1 • The MLL complex methylates H3K4 to upregulate transcription. MLL1 is the catalytic subunit which exhibits a low basal methyltransferase activity. Addition of WDR5, Ash2L, and RbBP5 (WAR) enhances Histone H3K4 methylation 300-fold.100 μg
250 μg
500 μgAutoradiograph of Core Histone (Cayman Item No. 11010) separated on SDS-PAGE following reaction with [³H]-SAM and the MLL1/WAR complex.Lane 1: MLL1/WAR complex
(1 μM each of MLL1, Ash2L,
RbBP5, WDR5)
Lane 2: MLL1/WAR complex
(3 μM each of MLL1, Ash2L,
RbBP5, WDR5)

NOTE: SUMO Affinity tag used under license from LifeSensors, Inc.

MLL1/WAR complex (human recombinant)


10945

*KMT2A, Mixed-Lineage Leukemia-1, MLL Core Complex***M_r:** ~200.1 kDa **Purity:** ≥90% **Stability:** ≥6 months at -80°C**Source:** Active recombinant proteins expressed in *E. coli*. MLL1, WDR5, Ash2L, RbBP5, and DPY-30 are expressed with N-terminal His-SUMOpro affinity tags. The His-SUMO affinity tags were removed using recombinant SUMO Protease 1 • The MLL complex methylates H3K4 to upregulate transcription. MLL1 is the catalytic subunit and contains the core components WDR5, Ash2L, and RbBP5 (MLL/WAR complex). MLL1 alone exhibits a low basal methyltransferase activity, but is enhanced 300-fold by the addition of the WAR components. A further 2-fold activation has been reported with the addition of DPY-30 (MLL/WAR complex).100 μg
250 μg
500 μgAutoradiograph of Core Histone (Cayman Item No. 11010) separated on SDS-PAGE following reaction with [³H]-SAM and the MLL1/WAR complex or MLL1/WAR complex.Lane 1: MLL1/WAR complex (1 μM each of MLL1,
Ash2L, RbBP5, WDR5, DPY-30)
Lane 2: MLL1/WAR complex (3 μM each of MLL1,
Ash2L, RbBP5, WDR5, DPY-30)
Lane 3: MLL1/WAR complex (1 μM each of MLL1,
Ash2L, RbBP5, WDR5)
Lane 4: MLL1/WAR complex (3 μM each of MLL1,
Ash2L, RbBP5, WDR5)

NOTE: SUMO Affinity tag used under license from LifeSensors, Inc.


NSD1 (human recombinant)

10757

*AR267, KMT3B, Nuclear Receptor SET domain-containing protein 1, WHSC1***M_r:** 60.2 kDa **Purity:** ≥60% **Stability:** ≥1 year at -80°C**Source:** Active recombinant N-terminal GST-tagged protein consisting of amino acids 1,700-1,986 expressed in *E. coli* • NSD1 methylates H3K36 and H4K20 and is important for maintaining chromatin integrity. It contains a catalytic lysine methyltransferase SET domain and four zinc-binding PHD fingers.25 μg
50 μg
100 μgLane 1: MW Markers
Lane 2: NSD1 (1 μg)
Lane 3: NSD1 (2 μg)
Lane 4: NSD1 (5 μg)


NSD2 (human recombinant)

10758

*Histone-lysine N-methyltransferase NSD2 isoform 1, KIAA1090, MMSET, Nuclear Receptor SET domain-containing protein 2, Protein trithorax-5, TRX5, WHSC1***M_r:** 61.2 kDa **Purity:** ≥90% **Stability:** ≥1 year at -80°C**Source:** Active recombinant N-terminal GST-tagged protein consisting of amino acids 941-1,240 expressed in *E. coli* • NSD2 has been implicated in multiple myeloma, an incurable malignancy in mature plasma cells, being involved in recurrent t(4;14) translocations with the immunoglobulin promotor/enhancer. NSD2 is responsible for the post-translational modification of histones H3 and H4. The methylation target of NSD2 is dependent on the nature of the substrate.25 μg
50 μg
100 μgLane 1: MW Markers
Lane 2: NSD2 (1 μg)
Lane 3: NSD2 (5 μg)
Lane 4: NSD2 (10 μg)


PRDM9 (human recombinant)

11209

*Meisetz, MSBP3, PFM6, PR Domain-containing Protein 9, ZNF899***M_r:** 60.5 kDa **Purity:** ≥65% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal His-tagged protein consisting of amino acids 2-511 expressed in *E. coli* • PRDM9 is a histone methyltransferase that binds specifically to recombination hotspots catalyzing trimethylation of H3K4 in nucleosomes near its binding site, thus initiating genetic recombination by recruiting recombination initiation machinery.10 μg
25 μg
50 μgLane 1: MW Markers
Lane 2: PRDM9 (4 μg)
Lane 3: PRDM9 (2 μg)


PRMT1 (human recombinant)

10350

*Protein Arginine Methyltransferase 1***M_r:** 69.2 kDa **Purity:** ≥80% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal GST-tagged protein consisting of amino acids 2-371 expressed in *E. coli* • PRMT1 is a class I arginine methyltransferase that methylates arginine residues at a number of glycine and arginine rich regions including histone H4 at arginine 3.25 μg
50 μg
100 μgLane 1: MW ladder
Lane 2: PRMT1 (1 μg)
Lane 3: PRMT1 (2 μg)
Lane 4: PRMT1 (4 μg)


PRMT3 (human recombinant)

11642

*Heterogeneous Nuclear Ribonucleoprotein Methyltransferase-like Protein 3, Protein Arginine methyltransferase 3***M_r:** 86.6 kDa **Purity:** ≥90% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal GST-tagged protein consisting of amino acids 2-531 expressed in *E. coli* • PRMT3 is a type-1 PRMT, catalyzing the formation of MMA and ADMA.25 μg
50 μg
100 μgLane 1: MW Markers
Lane 2: PRMT3 (2.5 μg)
Lane 3: PRMT3 (5 μg)

PRMT4 (human recombinant)

10750


*CARM1, Protein Arginine Methyltransferase 4***M_r:** 92.6 kDa **Purity:** ≥70% **Stability:** ≥1 year at -80°C**Source:** Active recombinant N-terminal GST-tagged protein consisting of amino acids 2-608 expressed in *E. coli* • PRMT4, also known as CARM1, is a type-1 PRMT, catalyzing MMA and ADMA on histone H3, Arg-17, and Arg-26.25 μg
50 μg
100 μgLane 1: MW Markers
Lane 2: PRMT4 (5 μg)
Lane 3: PRMT4 (2.5 μg)
Lane 4: PRMT4 (1.25 μg)

PRMT4 Polyclonal Antibody

13552

*CARM1*Protein G-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: human PRMT4 amino acid sequences 45-69 and 595-608 • Host: rabbit • Cross Reactivity: (+) human PRMT4 • Application(s): WB • PRMT4, also known as CARM1, belongs to a family of proteins that catalyzes the methylation of arginine residues.

1 ea



Lane 1: HeLa cell lysate


PRMT5 Polyclonal Antibody

13559

JBP1, Skb1HS

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human PRMT5 • Host: rabbit • Cross Reactivity:(+) human PRMT5 • Application(s): WB • PRMT5, also known as JBP1 and human homolog of Skb1 of fission yeast (Skb1HS), can catalyze the formation of MMA and symmetric dimethylarginine in a variety of proteins. Recombinant PRMT5 can mono- and dimethylate histone 2A and myelin basic protein.




PRMT6 (human recombinant)

10752

Histone Arginine N-methyltransferase PRMT6, HRMT1L6, Protein Arginine Methyltransferase 6

M_r: 68.7 kDa **Purity:** ≥80% **Stability:** ≥1 year at -80°C

Source: Active recombinant N-terminal GST-tagged protein consisting of amino acids 2-375 expressed in *E. coli* • PRMT6 is a nuclear type-1 PRMT, catalyzing the formation of MMA and ADMA on both histone and non-histone targets.




PRMT6 (human recombinant; baculovirus expressed)

13866

Histone Arginine N-methyltransferase PRMT6, HRMT1L6, Protein Arginine Methyltransferase 6

M_r: 43.7 kDa **Purity:** ≥60% **Stability:** ≥1 year at -80°C

Source: Active recombinant N-terminal His-tagged protein consisting of amino acids 2-375 expressed in Sf21 cells • PRMT6 is a nuclear type-1 PRMT, catalyzing the formation of MMA and ADMA on both histone and non-histone targets. Non-histone targets of PRMT6 include the nuclear high-mobility group (HMG) protein HMGA1a, a protein important in several processes relating to the maintenance of DNA integrity.




PRMT6 Polyclonal Antibody

13558

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human PRMT6 amino acids 23-43 • Host: rabbit • Cross Reactivity:(+) human and mouse PRMT6 • Application(s): WB • PRMT6 is a protein with an approximate molecular weight of 42 kDa consisting of a catalytic core sequence common to other PRMT enzymes. PRMT6 demonstrates type I PRMT activity, capable of forming both MMA and ADMA derivatives on recombinant glycine- and arginine-rich substrates.




PRMT7 Polyclonal Antibody

13551

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human PRMT7 amino acids 346-360 • Host: rabbit • Cross Reactivity:(+) human and mouse PRMT7 • Application(s): WB • PRMT7 can catalyze the formation of MMA in peptides.




PRMT8 (human recombinant)

11644

HRMT1L, Protein Arginine Methyltransferase 8

M_r: 65.7 kDa **Purity:** ≥90% **Stability:** ≥9 months at -80°C

Source: Active recombinant N-terminal GST-tagged protein consisting of amino acids 61-394 expressed in *E. coli* • PRMT8 is a type I methyltransferase shown to methylate the glycine-arginine rich region in fibrillarin. Its expression is limited to the brain.




RbBP5 (human recombinant)

10947

RBQ-3, Retinoblastoma-binding Protein 5, SWD1/Set1c WD40 repeat protein homolog

M_r: 72.9 kDa **Purity:** ≥90% **Stability:** ≥6 months at -80°C

Source: Recombinant protein consisting of amino acids 2-538 expressed in *E. coli* • RbBP5 is a ubiquitously expressed nuclear protein that contains WD40 repeat-like domains. RbBP5 binds directly to tumor suppressor retinoblastoma protein and regulates cell proliferation. RbBP5 is also an important component of the multi-subunit SET1 lysine methyltransferase protein complex, which includes MLL1.




RG-108

13302

[48208-26-0] N-Phthalyl-L-Tryptophan

MF: C₁₉H₁₄N₂O₄ **FW:** 334.3 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C

Summary: A non-nucleoside DNA methyltransferase inhibitor (IC₅₀ = 115 nM *in vitro*) that significantly reduces the methylation of genomic DNA in cells at 10 µM without detectable toxicity

5 mg
10 mg
50 mg
100 mg


RIZ1 (human recombinant)

10765

KMT8, PRDM2, PR Domain Zinc Finger Protein 2, Retinoblastoma Protein-interacting Zinc Finger Protein

M_r: 49.5 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C

Source: Recombinant N-terminal GST-tagged protein consisting of amino acids 2-200 expressed in *E. coli* • RIZ1 is a SAM-dependent histone methyltransferase that specifically methylates H3K9. RIZ1 is a tumor suppressor that can arrest the cell cycle and induce apoptosis.




SET7/9 (human recombinant)

10320

KMT7, SETD7/9, SET Domain-Containing Protein 7/9

M_r: 43.3 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C

Source: Active recombinant N-terminal His-tagged SET7/9 amino acids 1-366 expressed in *E. coli* • SET7/9 is exclusively a mono-methylase that methylates histone H3, tumor suppressor p53, and transcription factor TAF10. SET7/9 methylates p53 in response to DNA damage thereby activating p53 for subsequent acetylation.




SET7/9 Methyltransferase Inhibitor Screening Assay Kit

700270

Stability: ≥6 months at -80°C

Summary: SET7/9 is an MT that acts on various substrates including H3K4, p53, and the transcription factor TAF 10. In Cayman's SET7/9 MT Inhibitor Screening Assay the transfer of the methyl group from SAM by SET7/9 to the acceptor peptide (TAF 10) generates SAH, which is rapidly converted to urate and H₂O₂ using an enzyme mixture provided in the kit. A subsequent reaction between H₂O₂ and ADHP produces the highly fluorescent compound resorufin.

96 wells


Inhibition of SET7/9 by 5'-Deoxy-5'-(methylthio)adenosine (IC₅₀ ~2.5 mM)

SET7/9 Polyclonal Antibody

13731

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human SET7/9 amino acids 131-145 and 336-352 • Host: rabbit • Cross Reactivity: (+) mouse and human SET7/9 • Application(s): WB • SET7/9 is a histone specific HMTase that methylates histone H3K4.



Histone Methylation and the Language of Epigenetics

by [Thomas G. Brock, Ph.D.]

The basic processes of chromatin modification which influence gene expression are familiar to many of us: the addition and removal of acetyl and methyl marks, phosphorylation, ubiquitination, and so on. However, the exciting part, the complexity of how the marks are used, is just emerging. Some of the new terminology related to this, such as 'epigenetic editing', 'reading', and 'writing', begins to hint that chromatin modification is much like a language.^{1,2} This means that, perhaps, we have just learned some of the basic letters or words, with the more challenging issues of grammar, literal meaning, innuendo, and misleading information still to be contemplated.³ Consider the discovery of the four 'letters' of DNA, which was found to encode messages that were divided into segments, complete with phrases for starting and stopping and much, much more. We are still learning how to decipher the more complex information generated by these four letters. By analogy, a methylated lysine mark is one of the letters of epigenetic signaling. This article considers histone methylation and suggests directions for understanding its meaning.

Lysine Methyltransferases

Braille, the writing system developed for the visually impaired, uses a simple raised mark as its basic unit. Importantly, it is the organization of raised marks which creates a letter and, of course, you need a series of these to produce words and sentences. In the same way, a methyl mark is a basic unit of epigenetic 'reading'. This mark can be attached to lysine (K) or arginine (R)

Common Name	KMT Name	Mark
SUV39H1	KMT1A	H3K9me3
SUV39H2	KMT1B	H3K9me3
EHMT2, G9A	KMT1C	H3K9me1,2, H3K27me, H3K56me1
EHMT1, GLP	KMT1D	H3K9me1,2, H3K27me
SETDB1	KMT1E	H3K9me3
SETCB2	KMT1F	H3K9me3
MLL	KMT2A	H3K4me
MLL4	KMT2B	H3K4me
MLL3	KMT2C	H3K4me
MLL2	KMT2D	H3K4me
MLL5	KMT2E	H3K4me1,2
SETD1A, SET1A	KMT2F	H3K4me
SETD1B, SET1B	KMT2G	H3K4me
ASH1L	KMT2H	H3K4me, H3K36me
SETD2, SET2	KMT3A	H3K36me
NSD1	KMT3B	H3K36me, H4K20me
SMYD2	KMT3C	H3K4me, H3K36me2
DOT1L	KMT4	H3K79me
SET8, SETD8	KMT5A	H4K20me1
SUV420H1	KMT5B	H4K20me3
SUV420H2	KMT5C	H4K20me3
EZH2	KMT6	H3K9me, H3K27me1,2,3
SET7/9, SETD7	KMT7	H3K4me

Table 1. Human lysine methyltransferases




Figure 1. Structural organization of representative SET domain-containing proteins

on proteins as well as to cytosine on DNA. Up to three methyl groups can be added to K while R can be monomethylated as well as symmetrically or asymmetrically dimethylated. Should each of these be viewed, like Braille letters with different arrangements of raised dots, as being distinct letters of the epigenetic language? Certainly monomethylation at the three types of targets (K, R, cytosine) produces three distinct structures that are written, erased, and, most importantly, read by binding domains that are unique to each.

Focusing on lysine methyltransferases (KMTs), a list of those that have been assigned KMT names (Table 1) provides more information. There clearly is duplication of KMTs, with two similar enzymes making the same mark. This may indicate that the function of each pair is important enough to require a backup should one of the pair fail. Alternatively, the pair afford distinct regulatory pathways or tissue specific expression. In addition, methylation at H3K4 can be performed by a wide number of KMTs. Possibly, this mark is essential to chromatin housekeeping and plays only a general role in the epigenetic language. Finally, some enzymes specialize in monomethylating a particular residue (*e.g.*, SET8 at H4K20) while others specifically di- or trimethylate (*e.g.*, SUV420H1 and 2 at H4K20). This suggests that, at least in some cases, the accumulation of methyl groups at a given site is regulated, indicating that being trimethylated can have a different meaning from being monomethylated.

SET Domain Proteins

Additional insights may be derived from the methyl writers. To make the mark, many KMTs contain a SET domain, named after regions shared by three *Drosophila* proteins recognized as being involved in epigenetic processes: Su(var)3-9, Enhancer of zeste, and Trithorax. The SET domain includes conserved N- and C- terminal regions (SET-N, SET-C) and an intervening insert region (SET-I). Flanking pre- and post-SET regions are typically also required for full KMT activity.

Some KMTs have multiple functions. For example, the SMYD proteins are short KMTs that contain SET and MYND-type zinc finger domains (Figure 1). Like other zinc finger domains, MYND domains are involved in protein-

protein interactions, commonly binding a co-repressor protein, like NCoR or SMRT. SMYD1 acts as a transcriptional repressor, is essential for cardiomyocyte differentiation, and interacts with HDACs. SMYD3 specifically methylates H3K4, inducing di- and tri-methylation, but not monomethylation. The related SMYD2 methylates p53 and RB1, as well as H3K4 and H3K36.

The human SUV proteins are homologs of the *Drosophila* Su(var) proteins. There are two homologs, SUV39H1 and SUV39H2, which specifically trimethylate H3K9 after it has already been monomethylated. Both proteins contain N-terminal chromatin organization modifier (chromo) domains, which facilitate the condensation of heterochromatin.⁴ They function mainly in these condensed heterochromatin regions, suppressing gene expression. Trimethylation on H3K9 facilitates DNA methylation in this context. Two additional Su(var) homologs, SUV420H1 and SUV420H2, specifically trimethylate H4K20. Like the SUV39 homologs, these proteins are targeted to heterochromatin and are involved in epigenetic transcriptional repression.

Another structurally-defined family, the PRDM series, contains a PR domain, an evolutionarily conserved region of about 100 amino acids that is involved in protein-protein interactions. PRDM proteins also contain classical C2H2-type zinc finger domains which mediate DNA binding. PRDM1, also known as BLIMP1, acts as a transcriptional repressor, binding to the promoter of β -interferon, and in this way regulates B cell maturation. PRDM2, also known as RIZ, is another important family member. It methylates H3K9, binds the retinoblastoma protein, and is highly expressed in brain tumors.

The Enhancer of zeste homologs, EZH1 and EZH2, are polycomb group (PcG) proteins that can mono-, di- and trimethylate H3K27. The EZH proteins contain WD repeat binding domains, which mediate interaction with EED (embryonic ectoderm development) protein to form, with

KATs	Target	KMTs
	H3K4	ASH1L, MLL, MLL2, MLL3, MLL4, MLL5, SET1A, SET1B, SET5, SET7/9, SETMAR, SMYD1, SMYD2, SMYD3, WHSC1L1
ELP3, GCN5, MYST3, PCAF	H3K9	EHMT1, EHMT2, EZH2, SETDB1, SETDB2, SUV39H1, SUV39H2
CBP, CLOCK, P300, GCN5, GTF3C4, MYST2, MYST3, NCOA1, TAF1	H3K14	
CBP, ELP3, P300	H3K18	
CBP, P300	H3K27	EHMT2, EZH1, EZH2, WHSC1, WHSC1L1
	H3K36	ASH1L, MLL5, NSD1, SET2, SETMAR, SMYD2, WHSC1
CBP, P300, GCN5	H3K56	EHMT2
	H3K79	DOT1L
CBP, P300, GCN5, HAT1, MYST2, TIP60	H4K5	
CBP, P300, GCN5, MYST2, TIP60	H4K8	
CBP, P300, GCN5, HAT1, MYST2, TIP60	H4K12	
CBP, P300, GCN5, MYST1, TIP60	H4K16	
	H4K20	SET8, SUV420H1, SUV420H2, WHSC1
GCN5	H4K61	H3K36me

Table 2. Known HATs and KMTs of select H3, H4 lysine residues




Figure 2. Readers of nucleosomal marks may work alone (A), in combination on a given site (*cis*, B), in pairs on distinct intranucleosomal sites (*trans*, C), or in pairs on different internucleosomal targets (*trans*, D).

SUZ12 (suppressor of zeste 12 homolog), the polycomb repressor complex 2 (PRC2).⁵ Both EZH complexes play important roles in embryonic stem cell function.

A diverse group of SET domain-containing proteins is denoted as SETD. Two key members, SETD1A and SETD1B, methylate H3K4, but not if H3K9 is already methylated. Both proteins, which function as components of multimeric complexes, contain RNA binding domains (RBD). SETD2, unlike the SET1 proteins, methylates H3K36, binds DNA at promoters, and directly binds hyperphosphorylated RNA polymerase II large subunit. This latter interaction is mediated by a Set2Rpb1 interacting (SRI) domain and serves to couple H3K36 methylation with transcript elongation.

Thus, KMTs do not simply generate a methyl mark. All are involved in reading the nucleosomal setting, some requiring methylated marks in order to function. Others facilitate multiple changes, including adding or subtracting marks in addition to the methyl groups that it adds.

The Language of Lysine Methylation

In languages, meaning comes from combinations of letters or words. In the language of epigenetics, the message produced by chromatin remodeling enzymes may also be read as groups of marks. Inspection of the N-terminal tails of histones reveals that many of the residues that may be modified with negatively-charged phosphorylation marks are adjacent to lysines, whose charge may have been altered by acetylation. Remarkably, a recent report found that, in histones from *Plasmodium*, the majority of phosphorylation marks were found adjacent to acetylated lysines.⁶ This suggests a language 'rule' that pairs these marks together. Another observation comes from organizing histone targets according to the KATs and KMTs which modify them (Table 2). Certain targets (H3K9, H3K27, and H3K56) can be modified by either type of enzyme, indicating these sites may be controlling the message. This becomes more interesting if, as commonly thought, the rate of acetylation/deacetylation is fast while methylation/demethylation is slow.^{7,8} Moreover, it is possible that di- and trimethylation differ in stability compared to monomethylation and, as a result, may be functionally distinct, perhaps acting as punctuation (*e.g.*, stop) marks. It should be noted that the information in Table 2 reflects currently known enzyme-target assignments and dramatically underestimates the marks detected by proteomic analysis of histones. Such analyses, compiled at sites like PhosphoSitePlus, reveal that the majority of lysine residues on histones which can be targeted by KATs or KMTs are, in fact, modified by both.


Of course, the message of the marks depends on the readers. Just as nucleosomal marks should not be expected to function in isolation, binding proteins (readers) may assemble together.⁹ Readers of nucleosomal marks may work alone, in combination on a given site (*cis*), or in pairs on distinct intranucleosomal sites or on different internucleosomal targets (*trans*), as shown in Figure 2. Since readers serve to recruit or stabilize other proteins at the nucleosome, the interactions between different readers may be indirect, involving associated proteins. In this way, clusters of marks at distant sites may act together to create a complex message.

SET7/9 (FL) Polyclonal Antibody

13780

KMT7, SETD7/9, SET Domain-Containing Protein
Protein A-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: human recombinant SET7/9 (amino acids 1-366) • Host: rabbit
• Cross Reactivity: (+) human and mouse SET7/9 • Application(s): WB • SET7/9 is a histone specific HMTase that methylates histone H3 lysine 4.

500 µl



Lane 1: Standard
Lane 2: SET7/9 Recombinant Protein (0.001 µg)
Lane 3: SET7/9 Recombinant Protein (0.005 µg)
Lane 4: SET7/9 Recombinant Protein (0.01 µg)
Lane 5: SET7/9 Recombinant Protein (0.1 µg)
Lane 6: K562 cell lysate (15 µg)
Lane 7: K562 cell lysate (30 µg)
Lane 8: K562 cell lysate (50 µg)
Lane 9: HeLa cell lysate (50 µg)

SET7/9 SAM-Screener™ Assay Kit


600490

*SET7/9 SAM-Binding Site Inhibitor Screening***Stability:** ≥6 months at -80°C **Z' Factor:** 0.71**Summary:** SET7/9 (KMT7) is a SET domain-containing mono-methyltransferase that acts on a large number of histone and non-histone targets including histone H3, TAF10, p53, viral Tat, and estrogen receptor α . This fluorescence polarization assay is based upon a proprietary small molecule fluorescent probe that binds to the SAM binding pocket in SET7/9. Binding of the small molecule probe to SET7/9 induces an increase in fluorescence polarization. Binding of the probe can be competed with the endogenous cofactor SAM or by the inhibitor sinefungin, but is unaffected by the histone H3 peptide substrate. The SET7/9 SAM-Screener Assay is robust and exhibits a greater than 100 mP shift over a range of 0-250 nM SET7/9. The assay is suitable for high-throughput screening in the provided 384-well plate or can be scaled to higher density plate formats (e.g., 1,536-well) if desired.

384 wells


1,920 wells

SAM-Binding Site Probe Displacement



SET8 (human recombinant)

10319

KMT5a, PR-Set7, SETD8, SET domain-containing (lysine methyltransferase) 8
M_r: 21.1 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal His-tagged protein amino acids 190-352 expressed in *E. coli* • SET8 selectively mono-methylates histone H4 at lysine 20, an event proven to have an important role in chromatin structure and transcriptional activation. SET8 is also a novel regulator of p53, mono-methylating lysine 382 of the tumor suppressor. SET8's ability to suppress p53 transcriptional activity implies that it may play a significant role in tumorigenesis.25 µg
50 µg
100 µg


Lane 1: MW Markers
Lane 2: SET8 (4 µg)
Lane 3: SET8 (2 µg)
Lane 4: SET8 (1 µg)

SET8 Methyltransferase
Inhibitor Screening Assay Kit

700350


Stability: ≥6 months at -80°C**Summary:** SET8 is a methyltransferase that selectively mono-methylates H4K20, an event proven to have an important role in chromatin structure and transcriptional activation. Cayman's SET8 Methyltransferase Inhibitor Screening Assay provides a convenient method for screening human SET8 inhibitors. The transfer of the methyl group from SAM by SET8 (provided in the kit) to the acceptor peptide (H4K20) generates SAH, which is rapidly converted to urate and H₂O₂ using an enzyme mixture provided in the kit. The H₂O₂ formed is quantified using ADHP to produce the highly fluorescent compound resorufin (excitation 530-540 nm; emission 585-595 nm).

96 wells



SETD2 (human recombinant)


10767

Huntington Interacting Protein B, HYPB, KMT3A, p231HBP
M_r: 33.7 kDa **Purity:** ≥95% **Stability:** ≥9 months at -80°C**Source:** Active recombinant N-terminal His-tagged protein consisting of amino acids 1,435-1,711 expressed in *E. coli* • SETD2 is a histone methyltransferase that catalyzes the trimethylation of H3K36. The WW domain of SETD2 has been shown to interact with hyperphosphorylated RNA polymerase II, indicating a broad role in transcriptional regulation. Evidence has also been found indicating SETD2 can automethylate itself, possibly influencing its methyltransferase activity.25 µg
50 µg
100 µg

Lane 1: MW Markers
Lane 2: SETD2 (5 µg)
Lane 3: SETD2 (2.5 µg)
Lane 4: SETD2 (1 µg)

SMYD1 (human recombinant)

10761


*SET and MYND domain-containing protein 1, smk-BOP***M_r:** 69.6 kDa **Purity:** ≥50% **Stability:** ≥9 months at -80°C**Source:** Active recombinant N-terminal His- and SUMOpro-tagged protein consisting of amino acids 2-490 expressed in *E. coli* • SMYD1 is a cardiac and muscle-specific histone methyltransferase, specifically methylating H3K4, and is crucial for cardiomyocyte differentiation and maturation.10 µg
25 µg
50 µg

Lane 1: MW Markers
Lane 2: SMYD1 (4 µg)
Lane 3: SMYD1 (2 µg)

NOTE: SUMO Affinity tag used under license from LifeSensors, Inc.

SMYD3 (human recombinant)

10762


KMT3E, SET and MYND domain-containing protein 3
M_r: 58.5 kDa **Purity:** ≥60% **Stability:** ≥9 months at -80°C**Source:** Active recombinant N-terminal His- and SUMOpro-tagged protein consisting of amino acids 35-428 expressed in *E. coli* • SMYD3 is a histone methyltransferase that is overexpressed in several different cancers, such as breast, colorectal, and liver cancer. The N-terminus of full-length SMYD3 interacts with Hsp90 α and increases its activity, suggesting that the N-terminal region is involved in regulating methyltransferase activity.10 µg
25 µg
50 µg

Lane 1: MW Markers
Lane 2: SMYD3 (1 µg)
Lane 3: SMYD3 (2 µg)
Lane 4: SMYD3 (5 µg)

NOTE: SUMO Affinity tag used under license from LifeSensors, Inc.

SUV4-20H1 (human recombinant)


10763

*KMT5B, Lysine N-methyltransferase 5B, Suppressor of Variegation 4-20 Homolog 1***M_r:** 70.7 kDa **Purity:** ≥60% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal GST-tagged protein consisting of amino acids 2-387 expressed in *E. coli* • SUV4-20H1 is a SET domain containing methyltransferase responsible for di- and trimethylation of histone H4 lysine 20 (H4K20me2 and H4K20me3) at pericentric heterochromatin.25 µg
50 µg
100 µg

Lane 1: MW Ladders
Lane 2: KMT5B (0.5 µg)
Lane 3: KMT5B (1 µg)
Lane 4: KMT5B (2 µg)
Lane 5: KMT5B (4 µg)

SUV4-20H2 (human recombinant)

10764

*KMT5C, Suppressor of Variegation 4-20 Homolog 2***M_r:** 58.6 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal GST-tagged protein consisting of amino acids 2-280 expressed in *E. coli* • SUV4-20H2 is a SET domain containing methyltransferase responsible for di- and trimethylation of histone H4 lysine 20 (H4K20me2 and H4K20me3) at pericentric heterochromatin.25 µg
50 µg
100 µg

Lane 1: MW Ladder
Lane 2: KMT5C (1 µg)
Lane 3: KMT5C (2 µg)
Lane 4: KMT5C (4 µg)
Lane 5: KMT5C (6 µg)

TAF10 Peptide

10228

TAB-Associated Factor 10, TAFI130, TAF10 RNA polymerase II, TATA Box Binding Protein (TBP)-Associated Factor

FW: 1,267.0 **Supplied as:** 1 mg peptide lyophilized peptide from bicarbonate buffer **Stability:** ≥1 year at -20°C

Summary: TAF10 is one of many protein factors or coactivators associated with RNA polymerase II activity. One vial of this peptide may be used as a methyltransferase acceptor peptide for more than 200 reactions at 15 μM.

1 ea

WDR5 (human recombinant)


10944

BIG3, BMP2-induced 3-kb Gene Protein, SSET1c WD40 Repeat Protein, SWD3, WD-Repeat Protein 5

M_r: 34.4 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C

Source: Recombinant protein consisting of amino acids 23-334 expressed in *E. coli*
• WDR5 has been demonstrated to bind histone H3 by recognizing the first three amino acids of the N-terminal tail. Binding of WDR5 to a conserved arginine-containing motif in MLL-1, the so-called WDR5 interaction ("Win") motif, promotes the assembly and activity of the MLL core complex.

50 μg
100 μg
250 μg



NOTE: SUMO Affinity tag used under license from LifeSensors, Inc.

UNC0224

13631


[1197196-48-7]

MF: C₂₆H₄₃N₇O₂ **FW:** 485.7 **Purity:** ≥98%

A crystalline solid **Stability:** ≥2 years at -20°C

Summary: A potent and selective G9a HMTase inhibitor (IC₅₀ = 15 nM; K_d = 23 mM); more than 1,000-fold selective for G9a over SET7/9 and SET8

1 mg
5 mg
10 mg
50 mg



UNC0321 (trifluoroacetate salt)


10582

MF: C₂₇H₄₅N₇O₃ • 3(CF₃COOH) **FW:** 857.8 **Purity:** ≥95%

A solution in methyl acetate **Stability:** ≥1 year at -20°C

Summary: A potent and selective G9a HMTase inhibitor (IC₅₀ = 6 nM; K_i = 63 pM); more than 40,000-fold selective for G9a over SET7/9, SET8, PRMT3, and JMJD2E

1 mg
5 mg
10 mg
25 mg



UNC0638

10734


[1255580-76-7]

MF: C₃₀H₄₇N₅O₂ **FW:** 509.7 **Purity:** ≥98%

A crystalline solid **Stability:** ≥2 years at -20°C

Summary: A potent, selective G9a and GLP HMTase inhibitor (IC₅₀s = <15 and 19 nM, respectively); inhibits H3K9 dimethylation in MDA-MB231 cells (IC₅₀ = 81 nM) and demonstrates favorable separation of functional and toxic effects

1 mg
5 mg
10 mg
25 mg



Demethylases

Antibodies

- 55 10382** JMJD2A Polyclonal Antibody
57 13787 JMJD6 Peptide Affinity-Purified Polyclonal Antibody
57 13554 LSD1 Polyclonal Antibody (aa 100-150)
57 13553 LSD1 Polyclonal Antibody (aa 400-450)
57 13486 LSD1 Polyclonal Antibody (aa 450-500)
60 13555 LSD1 Polyclonal Antibody (aa 800-850)

Biochemicals

- 53 10599** Cl-Amidine
53 10610 F-Amidine (trifluoroacetate salt)
53 12033 Daminozide
54 11690 Gemcitabine
54 12054 GSK-J1 (sodium salt)
54 12073 GSK-J4 (hydrochloride)
54 12074 GSK-J5 (hydrochloride)
54 11572 IOX1
60 13944 N-Oxalylglycine
60 10010494 2-PCPA (hydrochloride)

Kits

- 53 700390** Demethylase (Jumonji-type) Activity Assay Kit
54 700400 Demethylase (LSD-type) Activity Assay Kit
55 700360 JMJD2A Inhibitor Screening Assay Kit
56 700370 JMJD2D Inhibitor Screening Assay Kit
57 700120 LSD1 Inhibitor Screening Assay Kit
60 700560 PAD4 Inhibitor Screening Assay Kit

Proteins

- 54 10336** JMJD2A (human recombinant)
55 11299 JMJD2A-Strep tagged (human recombinant)
55 10776 JMJD2C-Strep tagged (human recombinant)
56 10335 JMJD2D (human recombinant)
56 11300 JMJD2D-Strep tagged (human recombinant)
56 11237 JMJD2E-Strep tagged (human recombinant)
56 10773 JMJD6 (human recombinant)
60 10500 PAD4 (human recombinant)
60 10774 UTX (human recombinant)

Cl-Amidine

10599


[913723-61-2]

MF: C₁₄H₁₉ClN₄O₂ • CF₃CO₂H **FW:** 424.8 **Purity:** ≥95%

A crystalline solid **Stability:** ≥1 year at -20°C

Summary: An inhibitor of PAD4 deamination activity (IC₅₀ = 5.9 μM) that also inhibits PAD1 and PAD3 (IC₅₀ = 0.8 and 6.2 μM, respectively); dose dependently decreases the citrulline content in serum and joints and reduces the development of IgG autoantibodies in a CIA mouse model of inflammatory arthritis

1 mg
5 mg
10 mg
50 mg



NOTE: Sold under license from University of South Carolina under U.S. Patent No. 7,964,363

F-Amidine (trifluoroacetate salt)

10610


[877617-46-4]

MF: C₁₄H₁₉FN₄O₂ • CF₃COOH **FW:** 408.4 **Purity:** ≥95%

A solution in methanol **Stability:** ≥1 year at -20°C

Summary: Inhibits PAD4 activity (IC₅₀ = 21.6 μM) as well as PAD1 and PAD3 activity (IC₅₀s = 29.5 and 350 μM, respectively); cytotoxic to HL-60, MCF7, and HT-29 cancer cell lines (IC₅₀s = 0.5, 0.5 and 1 μM, respectively)

100 μg
250 μg
500 μg
1 mg



NOTE: Sold under license from University of South Carolina under U.S. Patent No. 7,964,363

Daminozide

12033


[1596-84-5] Alar, Aminozide, B 995, DIMG, DMASA, Kylar, SADH, Succinic Acid

MF: C₆H₁₂N₂O₃ **FW:** 160.2 **Purity:** ≥95%

A crystalline solid **Stability:** ≥2 years at -20°C

Summary: A selective inhibitor of the human 2-oxoglutarate (JmjC) histone demethylases KDM2A, PHF8, and KDM7A (IC₅₀s = 1.5, 0.55, and 2.1 μM, respectively)

5 g
10 g
25 g




Demethylase (Jumonji-type) Activity Assay Kit

700390

Stability: ≥6 months at -80°C

Summary: Lysine demethylases containing Jumonji C (JmjC) domains produce formaldehyde following 2-oxoglutarate-dependent demethylation. Cayman's Demethylase (Jumonji-type) Activity Assay provides a convenient fluorescence-based method for assaying JmjC-mediated demethylase activity from cell lysates or purified enzyme preparations. The assay is based on the production of formaldehyde during the demethylation of a methylated peptide substrate. Cyclization of formaldehyde and acetoacetanilide in the presence of ammonia gives a fluorescent product for quantitation.


96 wells



Demethylase (LSD-type) Activity Assay Kit 700400


Stability: ≥6 months at -80°C**Summary:** Cayman's Demethylase (LSD-type) Activity Assay provides a fluorescence-based method for assaying LSD-type demethylase activity. In this assay, formaldehyde is measured directly, eliminating the need for a coupled-enzyme reaction system. Formaldehyde, produced during demethylation of lysine 4 on a histone H3 peptide, reacts with the detection reagents provided in the kit to give a brightly fluorescent product. The assay is easy to use and can be completed in under two hours.

96 wells




Gemcitabine 11690


[95058-81-4] DDFC, Folfugem, Gemcel, GemLip, Gemzar, LY 188011, NSC 613327, Zefei

MF: C₉H₁₁N₃O₄F₂ **FW:** 263.2 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A nucleoside analog that arrests tumor growth and induces apoptosis by inhibiting DNA replication and repair; inhibits repair-mediated DNA demethylation inducing epigenetic gene silencing and has broad antiretroviral activity10 mg
25 mg
50 mg
100 mg


GSK-J1 (sodium salt) 12054

MF: C₂₂H₂₂N₅O₂ • Na **FW:** 411.4 **Purity:** ≥95%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A potent, cell impermeable H3K27 histone demethylase inhibitor highly selective for human JMJD3 (IC₅₀ = 60 nM *in vitro*)1 mg
5 mg
10 mg
50 mg

GSK-J4 (hydrochloride) 12073


MF: C₂₄H₂₇N₅O₂ • HCl **FW:** 454.0 **Purity:** ≥95%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** An ethyl ester prodrug of the JMJD3 selective histone demethylase inhibitor GSK-J1; reduces LPS-induced proinflammatory cytokine production, including that of TNFα (IC₅₀ = 9 µM) in human primary macrophages1 mg
5 mg
10 mg
50 mg

GSK-J5 (hydrochloride) 12074


MF: C₂₄H₂₇N₅O₂ • HCl **FW:** 454.0 **Purity:** ≥95%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A pyridine regio-isomer of the JMJD3 inhibitor GSK-J4; cell-permeable and hydrolyzed to a free base, which is a weak inhibitor of JMJD3 (IC₅₀ = 100 µM), making it an ideal negative control molecule1 mg
5 mg
10 mg

IOX1 11572

[5852-78-8]

MF: C₁₀H₇NO₃ **FW:** 189.2 **Purity:** ≥97%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A broad-spectrum inhibitor of 2OG oxygenases that inhibits JMJD2 demethylase activity (IC₅₀ = 87 µM); inhibits JMJD2A, JMJD2E and the 2OG oxygenases PHF8, PHD2, and FIH (IC₅₀s = 1.7, 2.4, 13.3, 14.3, and 20.5 µM, respectively)1 mg
5 mg


JMJD2A (human recombinant) 10336

*JHDM3A, Jumonji Domain Containing 2A, KDM4A***M_r:** 42.7 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal His-tagged protein consisting of amino acids 1-350 expressed in *E. coli* • JMJD2A catalyzes the demethylation of tri- and di-methylated forms of histone H3 at lysine residues 9 and 36.25 µg
50 µg
100 µgLane 1: MW Markers
Lane 2: JMJD2A (1 µg)
Lane 3: JMJD2A (2 µg)
Lane 4: JMJD2A (4 µg)

JMJD2A Inhibitor Screening Assay Kit 700360

*Jumonji Domain Containing 2A, KDM4A***Stability:** ≥6 months at -80°C**Summary:** JMJD2A is a JmjC histone demethylase that catalyzes the demethylation of trimethylated lysine 9 and lysine 36 of histone H3. Cayman's JMJD2A Inhibitor Screening Assay Kit provides a convenient fluorescence-based method for screening JMJD2A inhibitors. The assay is based on the multistep reaction in which JMJD2A first produces formaldehyde during the demethylation of the trimethylated peptide substrate, histone H3 trimethyl lys9, with the concomitant oxidation/decarboxylation of 2-oxoglutarate. The detection reaction involves the cyclization between formaldehyde and acetoacetanilide in the presence of ammonia. The resulting fluorescent product is analyzed using an excitation wavelength between 365-375 nm and an emission wavelength between 465-475 nm.


96 wells




JMJD2A Polyclonal Antibody 10382

*Jumonji Domain Containing 2A, KDM4A, Lysine-Specific Demethylase 4A*Antigen affinity-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: recombinant human JMJD2A amino acids 1-350 • Cross Reactivity: (+) human JMJD2A • Application(s): WB • JMJD2A is a lysine specific demethylase with emerging roles in histone modification or epigenetic remodeling. This JMJD2A polyclonal antibody was raised against an N-terminal recombinant fragment of JMJD2A. This fragment (amino acids 1-350) includes the JMJD and JMJC domains but not the two LAP/PHD zinc finger or Tudor domains of the 1,064 amino acid protein.


1 ea

Lane 1: DLD1 cell lysate (30 µg)
Lane 2: DLD1 cell lysate (60 µg)


JMJD2A-Strep tagged (human recombinant) 11299

*JHDM3A, Jumonji Domain Containing 2A, KDM4A***M_r:** 42.7 kDa **Purity:** ≥95% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal Strep II-tagged protein consisting of amino acids 2-350 expressed in *E. coli* • JMJD2A catalyzes the demethylation of trimethylated forms of histone at lysine residues 9 and 36. JMJD2A is an α-ketoglutarate-dependent Fe (II) oxygenase. Purification of Fe-dependent JmjC family members by IMAC can result in displacement of the catalytic iron and decreased activity. Therefore this protein is purified by Strep-Tactin affinity chromatography.25 µg
50 µg
100 µgLane 1: MW Markers
Lane 2: JMJD2A (2 µg)
Lane 3: JMJD2A (4 µg)
Lane 4: JMJD2A (6 µg)

JMJD2C-Strep tagged (human recombinant) 10776

*JMJD3C, Jumonji Domain Containing 2C, KDM4C, KIAA0780, Lysine-specific Demethylase 4C***M_r:** 45.3 kDa **Purity:** ≥90% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal Strep II-tagged protein consisting of amino acids 2-372 expressed in *E. coli* • JMJD2C catalyzes the demethylation of trimethylated histone H3 at lysine residues 9 or 36 (me 2/3), leading to transcriptional changes. JMJD2C is an α-ketoglutarate-dependent Fe (II) oxygenase. Purification of Fe-dependent JmjC family members by IMAC can result in displacement of the catalytic iron and decreased activity. Therefore this protein is purified by Strep-Tactin affinity chromatography.25 µg
50 µg
100 µgLane 1: MW Markers
Lane 2: JMJD2C (1 µg)
Lane 3: JMJD2C (2 µg)
Lane 4: JMJD2C (4 µg)

JMJD2D (human recombinant) 10335

*Jumonji Domain Containing 2D, KDM4D, Lysine-Specific Demethylase 4D***M_r:** 42.7 kDa **Purity:** ≥75% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal His-tagged protein consisting of amino acids 1-354 expressed in *E. coli* • JMJD2D catalyzes the demethylation of di- and trimethylated forms of histone H3 at lysine residue 9 (me2/3), leading to transcriptional repression and activation, respectively.25 µg
50 µg
100 µgLane 1: MW ladder
Lane 2: JMJD2D (1 µg)
Lane 3: JMJD2D (2 µg)
Lane 4: JMJD2D (4 µg)
Lane 5: JMJD2D (6 µg)


JMJD2D Inhibitor Screening Assay Kit

700370

*Jumonji Domain Containing 2D***Stability:** ≥6 months at -80°C

Summary: JMJD2D is a JmjC histone demethylase that catalyzes the demethylation of di- and trimethylated lysine 9 of histone H3. Cayman's JMJD2D Inhibitor Screening Assay Kit provides a convenient fluorescence-based method for screening JMJD2D inhibitors. The assay is based on the multistep reaction in which JMJD2D first produces formaldehyde during the demethylation of the trimethylated peptide substrate, histone H3 trimethyl lysine 9, with the concomitant oxidation/decarboxylation of 2-oxoglutarate. The detection reaction involves the cyclization between formaldehyde and acetoacetanilide in the presence of ammonia. The resulting fluorescent product is analyzed using an excitation wavelength between 365-375 nm and an emission wavelength between 465-475 nm.

96 wells




JMJD2D-Strep tagged (human recombinant)

11300

*Jumonji Domain Containing 2D, KDM4D, Lysine-specific Demethylase 4D***M_r:** 42.6 kDa **Purity:** ≥70% **Stability:** ≥6 months at -80°C

Source: Active recombinant N-terminal Strep II-tagged protein consisting of amino acids 4-354 expressed in *E. coli*. JMJD2D catalyzes the demethylation of di- and trimethylated forms of histone H3 at lysine residue 9 (me2/3), leading to transcriptional repression and activation, respectively. JMJD2D is an α-ketoglutarate-dependent Fe (II) oxygenase. Purification of Fe-dependent JmjC family members by IMAC can result in displacement of the catalytic iron and decreased activity. Therefore this protein is purified by Strep-Tactin affinity chromatography.


25 μg
50 μg
100 μg

JMJD2E-Strep tagged (human recombinant)

11237

*Jumonji Domain Containing 2E, KDM4D-Like***M_r:** 40.8 kDa **Purity:** ≥75% **Stability:** ≥6 months at -80°C

Source: Active recombinant N-terminal Strep II-tagged protein consisting of amino acids 2-337 expressed in *E. coli*. JMJD2E catalyzes the demethylation of histone H3 at lysine residue 9. JMJD2E is an α-ketoglutarate-dependent Fe (II) oxygenase. Purification of Fe-dependent JmjC family members by IMAC can result in displacement of the catalytic iron and decreased activity. Therefore this protein is purified by Strep-Tactin affinity chromatography.


25 μg
50 μg
100 μg

JMJD6 (human recombinant)

10773

*Jumonji Domain Containing 6, KIAA0585, PTDSR1***M_r:** 73.2 kDa **Purity:** ~50% **Stability:** ≥6 months at -80°C

Source: Recombinant N-terminal GST-tagged protein consisting of amino acids 2-403 expressed in *E. coli*. JMJD6 was initially proposed to function as a histone arginine demethylase. Further studies suggest that JMJD6 has hydroxylase activity towards histone tails, the splicing regulatory protein LUC7-Like2, or single-stranded RNA.

25 μg
50 μg
100 μg

JMJD6 Peptide Affinity-Purified


Polyclonal Antibody

13787

*Jumonji Domain Containing 6, PTDSR*Antigen affinity-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: human JMJD6 amino acids 127-144 • Host: rabbit • Cross Reactivity: (+) chimpanzee, ovine, canine, equine, human, mouse, and opossum JMJD6 • Application(s): WB • JMJD6 is a 403 amino acid nuclear protein lysyl-hydroxylase that has been reported to have arginine demethylase activity for histone H3 at 'Arg-2' and histone H4 at 'Arg-3'. JMJD6 has been suggested to function in the differentiation of multiple organs during embryogenesis and regulate hematopoietic differentiation and macrophage cytokine responses.

1 ea




LSD1 Inhibitor Screening Assay Kit

700120

*Lysine-Specific Demethylase 1***Stability:** ≥6 months at -80°C

Summary: LSD1 is a histone demethylase whose actions on specific lysine residues repress transcription of chromosomal DNA. LSD1 also inhibits the tumor suppressor activity of p53 by demethylating a specific lysine residue. Cayman's LSD1 Inhibitor Screening Assay Kit provides a convenient fluorescence-based method for screening LSD1-specific inhibitors. The assay is based on the multistep enzymatic reaction in which LSD1 first produces H₂O₂ during the demethylation of lysine 4 on a peptide corresponding to the first 21 amino acids of the N-terminal tail of histone H3. In the presence of horseradish peroxidase, H₂O₂ reacts with ADHP to produce the highly fluorescent compound resorufin that can be analyzed with an excitation wavelength of 530-540 nm and an emission wavelength of 585-595 nm.

96 wells




LSD1 Polyclonal Antibody (aa 100-150)

13554

*Amine Oxidase (flavin containing) Domain 2*Protein affinity-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: peptide corresponding to a portion of human LSD1 amino acids 100-150 • Host: rabbit • Cross Reactivity: (+) canine, human, mouse, rat, Rhesus monkey, and zebrafish LSD1 • Application(s): WB • LSD1 functions as a transcriptional corepressor and catalyzes the flavin-dependent demethylation of Lys4 of histone 3 resulting in the formation of methyl-free lysine and release of formaldehyde. It is typically associated with CoREST and HDACs 1 and 4 and participates in the silencing of endogenous neuron-specific genes.

1 ea




LSD1 Polyclonal Antibody (aa 400-450)

13553

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: peptide from a portion of human LSD1 amino acids 400-450 • Host: rabbit • Cross Reactivity: (+) chimpanzee, bovine, canine, human, monkey, and mouse LSD1 • Application(s): WB

1 ea




LSD1 Polyclonal Antibody (aa 450-500)

13486

Protein G-purified IgG **Stability:** ≥1 year at -20°C

Summary: Antigen: peptide within the region of human LSD1 amino acids 450-500 • Host: rabbit • Cross Reactivity: (+) chimpanzee, bovine, canine, equine, human, mouse, orangutan, and porcine LSD1 • Application(s): WB

1 ea



Histone Demethylation

by [Daniel A. Bochar, Ph.D.]

The alteration of chromatin structure provides a key regulatory step for all processes that act upon DNA. One mechanism for inducing these changes is through the posttranslational modification of histones. Of the multitude of covalent modifications that occur on histones, lysine methylation plays a central role in the epigenetic regulation of the genome.¹ The methylation state of histone residues is dynamically regulated through the opposing activities of histone methyltransferases and demethylases. Histones are methylated on numerous lysine or arginine residues, and to add further complexity, lysine residues may be mono-, di-, or tri-methylated, and arginine residues may be mono-, symmetrically, or asymmetrically di-methylated. Unlike the other common histone modifications of acetylation or phosphorylation, methylation does not alter the charge of the modified residue. Charge alteration of histone can directly impact histone-DNA or histone-histone interactions, suggesting that histone methylation elicits its effect by recruiting or blocking association of chromatin associated factors that can direct relevant activities to appropriate regions of the genome.

The histone modifications of acetylation and phosphorylation are dynamically regulated through opposing enzymatic actions; i.e. acetyltransferases/deacetylases or kinases/phosphatases. Until recently, histone methylation

was thought to be a static mark lacking this dynamic relationship. This assumption was based on the high stability of the carbon-nitrogen bond, and therefore removal of the methyl mark was thought to occur only through histone turnover or proteolytic removal of the histone tails.² The discovery of enzymes capable of removing this mark has solidified the dynamic nature of histone methylation, and has set the stage for a wide range of biological and pharmacological discoveries.

In 2004, LSD1 was identified as the first histone demethylase.³ This enzyme, which primarily demethylates H3K4, belongs to a larger family of FAD-dependent amine oxidases. LSD1 catalyzed demethylation occurs through the two electron oxidation of the amine by FAD to produce an iminium ion. This intermediate spontaneously hydrolyzes to produce formaldehyde and a product lysine less one methyl group (Figure 1). It is important to note that formation of the imine intermediate requires free electrons on the epsilon nitrogen and therefore only mono- or di-methylated lysine can serve as a substrate. LSD1 was first identified in the BRAF-HDAC complex along with the transcriptional corepressor, CoREST.⁴ This association suggests LSD1-catalyzed demethylation is important for transcriptional repression. Subsequently, LSD1 was identified as an androgen receptor-interacting





Figure 1. Lysine demethylation by LSD1 and JmjC containing proteins.
 A. FAD-dependent lysine demethylation catalyzed by LSD1 proceeds via an amine oxidase mechanism.
 B. Fe(II)-dependent lysine demethylation catalyzed by JmjC containing proteins proceeds via a dioxygenase mechanism.



protein and in this context LSD1 changed substrate specificity to H3K9,⁵ suggesting that LSD1 can act as a corepressor or coactivator based on its associated factors. In humans, there exists only one LSD1 paralog, LSD2, and this protein has also been shown to possess histone demethylase activity for H3K4 mono- and di-methyl.⁶

The LSD1 family lacks the ability to demethylate tri-methylated lysine residues and has only been shown to demethylate at H3K4 or H3K9. This is in contrast to the diversity of histone methyltransferases and resulting multitude of methyl marks that can occur on histones. For histone methylation to be truly a dynamic mark, a larger, more diverse family of histone demethylases would be needed to match the complexity of methyltransferases. A new family of histone demethylases was soon discovered with the characterization of FBXL11 (JHDM1A).⁷ This protein was identified following biochemical purification of an H3K36 demethylase activity from HeLa nuclear extracts. FBXL11 was previously uncharacterized, but contained an interesting conserved domain; the Jumonji-C domain (JmjC). The JmjC domain is a member of a larger domain family found in α -ketoglutarate dependent Fe(II) dioxygenases that catalyze a variety of cellular reactions. JmjC-catalyzed demethylation proceeds through the oxidative decarboxylation of α -ketoglutarate coupled to hydroxylation of the methyl group. This hydroxylated intermediate spontaneously decomposes to produce formaldehyde and a product lysine less one methyl group (Figure 1). Unlike LSD1, this mechanism does not require free electrons on the epsilon nitrogen to create an imine intermediate, allowing JmjC demethylases to utilize tri-methylated substrates. Bioinformatic searches of the human genome identify at least 31 JmjC domain-containing proteins, many of which possess histone demethylase activity (Figure 2). The

identification of a larger, more diverse family of histone demethylases adds the complexity needed to match the histone methyltransferases; thus creating a truly dynamic mark.

As shown in Figure 2, not all JmjC domain-containing proteins have an identified substrate. The search for JmjC domain substrates is complicated by the fact that these enzymes may not function as histone demethylases, as illustrated by JMJD6. The enzyme JMJD6 has been shown to be a lysyl hydroxylase catalyzing lysine C-5 hydroxylation of arginine-serine-rich regions of splicing regulatory proteins.⁸ The crystal structure of the related JMJD5 suggests that this protein may also function as a lysyl hydroxylase.⁹ A third member of this phylogenetic clade, FIH, also functions as a hydroxylase, this time as an asparaginyl hydroxylase targeting the transcription factor HIF1 α . Therefore, it is clear that not all JmjC domain-containing proteins will function as histone demethylases and identifying the JmjC domain-containing proteins with epigenetic activity is an area of active research.


Initial comparisons of the catalytic activity of JmjC domain-containing proteins to LSD1 suggest that JmjC domain-containing proteins are over a magnitude lower in activity than that of LSD1.¹⁰ This had initially raised speculation as to whether or not these enzymes are biologically relevant demethylases. Interestingly, JmjC domain-containing proteins are inhibited by divalent transition metals (e.g. Ni(II)) that can displace the catalytically required Fe(II).¹¹ Common purification schemes for many of the JmjC domain-containing proteins have included metal affinity chromatography, raising the possibility of inhibition of the purified target enzyme. Indeed, metal analysis of metal affinity purified JMJD2A and JMJD2D revealed a Ni(II) content of approximately 70%.¹⁰ In order to avoid

Ni(II) contamination, Krishnan, *et al.* devised a purification strategy utilizing a strep-tag to minimize the exposure to inhibitory transition metals.¹⁰ Using this approach, they demonstrate that the activity of JMJD2A and JMJD2D are approximately equal to that of the LSD1. These results validate the role of JmjC domain-containing proteins in the demethylation of histone substrates.

JmjC domain-containing demethylases work in concert with histone methyltransferases to control the methylation patterns of chromatin. These patterning enzymes are important for maintaining normal gene transcription and genomic stability, and disruption of this patterning is seen in many human diseases. One important example is the EZH2/UTX regulation of H3K27 methylation. Overexpression or mutation in EZH2, the histone methyltransferase responsible for this methyl mark, is seen in many prostate and breast cancers.¹² Mutations in UTX, the histone demethylase that removes this mark, have also been described in human cancers.¹³ Together, these findings highlight the importance of these balancing enzymes. While there has been a great deal of interest in the development of histone methyltransferase inhibitors as targeted cancer chemotherapeutics, the opposing action of histone demethylases suggests that these enzymes may also represent important drug targets. Indeed, the first selective inhibitor of JmjC domain-containing demethylases was recently described.¹⁴ The discovery of this JMJD3/UTX specific inhibitor demonstrates that JmjC domain-containing demethylases can be selectively inhibited and that this class of enzymes is a tractable area for epigenetic drug discovery. This discovery will hopefully expedite the development of inhibitors for other JmjC domain-containing family members that have been shown to malfunction in human diseases.

LSD1 Polyclonal Antibody (aa 800-850)


13555

Protein G-purified IgG **Stability:** ≥1 year at -20°C**Summary:** Antigen: peptide from human LSD1 within the range of amino acids 800-850 • Host: rabbit • Cross Reactivity: (+) canine, human, mouse, rat, and Rhesus monkey LSD1 • Application(s): IHC (paraffin-embedded sections) and WB

N-Oxalylglycine


13944

[5262-39-5] NOG

MF: C₄H₅NO₅ **FW:** 147.1 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** A cell permeable inhibitor of α-ketoglutarate-dependent enzymes, including JMJD2A, JMJD2C, and JMJD2E (IC₅₀s = 250, 500, and 24 μM, respectively); inhibits the prolyl hydroxylase domain-containing proteins PHD1 and PHD2 with IC₅₀ values of 2.1 and 5.6 μM, respectively5 mg
10 mg
50 mg
100 mg

PAD4 (human recombinant)

10500


*Peptidylarginine Deiminase 4, Protein Arginine Deiminase 4***M_r:** 75.8 kDa **Purity:** ≥95% **Stability:** ≥9 months at -80°C**Source:** Active recombinant N-terminal His-tagged protein consisting of amino acids 2-663 expressed in *E. coli* • PAD4 is a homodimer that functions as a transcriptional coregulator to catalyze the conversion of specific arginine residues to citrulline in a calcium-dependent manner. PAD4 substrates include histones H2A, H3, and H4, whose post-translational modifications play a large role in gene regulation.50 μg
100 μg
250 μg

PAD4 Inhibitor Screening Assay Kit

700560


*Peptidylarginine Deiminase 4, Protein Arginine Deiminase 4***Stability:** ≥6 months at -80°C**Summary:** PAD4 is a guanidino-modifying enzyme that functions as a transcriptional coregulator catalyzing the conversion of specific arginine residues to citrulline. Substrates for PAD4 include histones H2A, H3, and H4. PAD4 autocitrullinates itself at several sites, inhibiting its enzymatic activity. PAD4 activity is increased in rheumatoid arthritis, producing an abundance of citrulline-containing proteins that generate an immune response resulting in production of autoantibodies that ultimately attack the host tissues. PAD4 has also been implicated in several other diseases including multiple sclerosis, Alzheimer's disease, glaucoma, and cancer. Cayman's PAD4 Inhibitor Screening Assay provides a convenient, fluorescence-based method for screening human PAD4 inhibitors.

96 wells




2-PCPA (hydrochloride)

10010494

*trans-2-Phenylcyclopropylamine, Translycpropamine***MF:** C₉H₁₁N • HCl **FW:** 169.7 **Purity:** ≥98%A crystalline solid **Stability:** ≥2 years at -20°C**Summary:** An irreversible, mechanism-based inhibitor of LSD1 with an IC₅₀ value of 20.7 μM and a K_i value of 242 μM that effectively inhibits histone demethylation *in vivo*; irreversibly inhibits monoamine oxidases (MAO) A and MAO B with IC₅₀ values of 2.3 and 0.95 μM and K_i values of 101.9 and 16 μM, respectively10 mg
50 mg
100 mg
250 mg

UTX (human recombinant)

10774

*KDM6A, Ubiquitously Transcribed Tetratricopeptide Repeat X***M_r:** 70.5 kDa **Purity:** ≥85% **Stability:** ≥6 months at -80°C**Source:** Active recombinant N-terminal proprietary tagged protein consisting of amino acids 930-1,410 expressed in *E. coli* • UTX plays a crucial role in epigenetic regulation of gene expression by catalyzing the demethylation of tri-methylated lysine 27 on histone H3.25 μg
50 μg
100 μg

Article References

DNA Methylation
establishing & editing the methylome
page 4

- Heyn, H., Li, N., Ferreira, H.J., et al. *Proc. Natl. Acad. Sci.* **109**, 10522-10527 (2012).
- Klose, R.J. and Bird, A.P. *Trends Biochem. Sci.* **31**(2), 89-97 (2006).
- Clouaire, T., de las Heras, J.I., Merusi, C. et al. *Nucleic Acids Res.* **38**(14), 4620-4634 (2010).
- Chédin, F., Lieber, M.R., and Hsieh, C.L. *Proc. Natl. Acad. Sci.* **99**(26), 16916-16921 (2002).
- Kelsey, G. *Cell Res.* **21**, 1155-1156 (2011).
- Kelsey, G. and Feil, R. *Phil. Trans. R. Soc. B* **368**, 20110336 (2012).
- Bird, A. *Genes Dev.* **16**, 6-21 (2002).
- Smallwood, S.A. and Kelsey, G. *Trends Genet.* **28**(1), 33-42 (2012).
- Pawlak, M. and Jaenisch, R. *Genes Dev.* **25**, 1035-1040 (2011).
- Ooi, S.K.T., and Bestor, T.H. *Cell* **133**, 1145-1148 (2008).
- Bhutani, N., Burns, D.M., and Blau, H.M. *Cell* **146**, 866-872 (2011).
- Yang, X., Lay, F., Han, H. et al. *Trends Pharmacol. Sci.* **31**(11), 536-546 (2010).
- McCabe, M.T., Brandes, J.C., and Vertino, P.M. *Clin. Cancer Res.* **15**, 3927-3937 (2009).
- Handy, D.E., Castro, R., and Loscalzo, J. *Circulation* **123**, 2145-2156 (2011).

Histone Acetylation, Deacetylation, and
Genomic Bistability
page 16

- Ptashne, M. *FASEB J.* **27**(1), 1-2 (2013).
- Walsh, C.A., Qin, L., Tien, J.C.Y., et al. *Int. J. Biol. Sci.* **8**(4), 470-485 (2012).
- Martin, D.G.E., Grimes, D.E., Baetz, K., et al. *Mol. Cell Biol.* **26**(8), 3018-3028 (2006).
- Katan-Khaykovich, Y. and Struhl, K. *Genes Dev.* **16**, 743-752 (2002).
- Dodd, I.B., Micheelsen, M.A., Sneppen, K., et al. *Cell* **129**, 813-822 (2007).
- Danzer, J.R. and Wallrath, L.L. *Development* **131**(15), 3571-3580 (2004).
- Johnson, A., Li, G., Sikorski, T.W., et al. *Mol. Cell* **35**(6), 769-781 (2009).
- Margueron, R., Justin, N., Ohno, K., et al. *Nature* **461**, 762-767 (2009).
- Wysocka, J., Swigut, T., Milne, T.A., et al. *Cell* **121**, 859-872 (2013).
- Lai, K., Robertson, M.J., and Schaffer, D.V. *Biophys. J.* **86**(5), 2748-2757 (2004).
- Lopes, F.J.P., Vieira, F.M.C., Holloway, D.M., et al. *PLoS Comput. Biol.* **4**(9), (2008).
- Sneppen, K. and Dodd, I.B. *PLoS Comput. Biol.* **8**(8), (2012).

Bromodomains, ChromoHub, & STRING
page 26

- Filippakopoulos, P., Picaud, S., Mangos, M., et al. *Cell* **149**(1), 214-231 (2012).
- Denis, G.V. *Discov. Med.* **10**(55), 489-499 (2010).
- Belkina, A.C. and Denis, G.V. *Nat. Rev. Cancer* **12**, 465-477 (2012).
- Cliffe, S.T., Bloch, D.B., Suryani, S., et al. *J. Allergy Clin. Immunol.* **130**(3), 735-742 (2012).
- Field, M., Tarpey, P.S., Smith, R., et al. *Am. J. Hum. Genet.* **81**(2), 367-374 (2007).
- Bueno, C., Montes, R., Catalina, P., et al. *Leukemia* **25**, 400-410 (2011).
- Tamkun, J.W., Deuring, R., Scott, M.P., et al. *Cell* **68**, 561-572 (1992).
- Zipppo, A., Serafini, R., Rocchigiani, M., et al. *Cell* **138**, 1122-1136 (2009).
- Liu, L., Zhen, X.T., Denton, E., et al. *Bioinformatics* **28**(16), 2205-2206 (2012).
- Szklarczyk, D., Franceschini, A., Kuhn, M., et al. *Nucleic Acids Res.* **39**, D561-D568 (2011).

Histone Methylation and the Language
of Epigenetics
page 48

- de Groot, M.L., Verschure, P.J., and Rots, M.G. *Nucleic Acids Res.* **40**(21), (2012).
- Janzen, W.P., Wigle, T.J., Jin, J., et al. *Drug Discov. Today Technol.* **7**(1), e59-e65 (2010).
- Baker, L.A., Allis, C.D., and Wang, G.G. *Mutat. Res.* **647**(1-2), 3-12 (2008).
- Maurer-Stroh, S., Dickens, N.J., Hughes-Davies, L., et al. *Trends Biochem. Sci.* **28**(2), 69-74 (2003).
- Migliori, V., Mapelli, M., and Guccione, E. *Epigenetics* **7**(8), 815-822 (2012).
- Dastidar, E.G., Dzeyk, K., Krijgsveld, J., et al. *PLoS One* **8**(1), (2013).
- Katan-Khaykovich, Y. and Struhl, K. *Genes Dev.* **16**, 743-752 (2002).
- Klose, R.J. and Zhang, Y. *Nat. Rev. Mol. Cell Biol.* **8**(4), 307-318 (2007).
- Musselman, C.A., Lalonde, M.E., Côté, J., et al. *Nat. Struct. Mol. Biol.* **19**, 1218-1227 (2012).

Histone Demethylation
page 58

- Bannister, A.J. and Kouzarides, T. *Cell Res.* **21**, 381-395 (2011).
- Bannister, A.J., Schneider, R., and Kouzarides, T. *Cell* **109**, 801-806 (2002).
- Shi, Y., Lan, F., Matson, C., et al. *Cell* **119**, 941-953 (2004).
- Hakimi, M.-A., Bochar, D.A., Chenoweth, J., et al. *Proc. Natl. Acad. Sci. USA* **99**(11), 7420-7425 (2002).
- Metzger, E., Wissmann, M., Yin, N., et al. *Nature* **437**, 436-439 (2005).
- Ciccone, D.N., Su, H., Hevi, S., et al. *Nature* **461**, 415-418 (2009).
- Tsukada, Y., Fang, J., Erdjument-Bromage, H., et al. *Nature* **439**(16), 811-816 (2006).
- Webby, C.J., Wolf, A., Gromak, N., et al. *Science* **325**, 90-93 (2009).
- Del Rizzo, P.A., Krishnan, S., and Trievel, R.C. *Mol. Cell Biol.* **32**(19), 4044-4052 (2012).
- Krishnan, S., Collazo, E., Ortiz-Tello, P.A., et al. *Anal. Biochem.* **420**, 48-53 (2012).
- Sekirnik, R., Rose, N.R., Mecinovic, J., et al. *Metalomics* **2**(6), 397-399 (2010).
- Chase, A. and Cross, N.C.P. *Clin. Cancer Res.* **17**(9), 2613-2618 (2011).
- Grasso, C.S., Wu, Y.-M., Robinson, D.R., et al. *Nature* **487**, 239-243 (2012).
- Kruidenier, L., Chung, C.-W., Cheng, Z., et al. *Nature* **488**, 404-408 (2012).

Alphabetical Index

A

Absent, small, or homeotic discs 2-like (Ash2L (human recombinant))	39
N-Acetyldinaline (CI-994)	13
Acetyl Lysine Monoclonal Antibody (Clone 7F8)	9
Acetyl Lysine Polyclonal Antibody-biotin	9
Acetyl Lysine Polyclonal Antibody HRP Conjugate	9
S-(5'-Adenosyl)-L-methionine chloride (hydrochloride)	38
AdoMet (S-(5'-Adenosyl)-L-methionine chloride (hydrochloride))	38
AGK2	12
AK-7	12
Alar (Daminozide)	53
ALL1 (MLL1 (human recombinant))	44
AMI-1 (sodium salt)	38
CI-Amidine	53
F-Amidine (trifluoroacetate salt)	53
Amine Oxidase (flavin containing) Domain 2 (LSD1 Polyclonal Antibody (aa 100-150))	57
N-phenyl-N'-(2-Aminophenyl)hexamethylenediamide (CAY10433)	13
Aminoazide (Daminozide)	53
Anacardic Acid	9
Antibiotic U 18496 (5-Azacytidine)	39
Antigen p25 (HP1- α (human recombinant))	36
Apicidin	12
AR42 ((S)-HDAC-42)	13
AR267 (NSD1 (human recombinant))	44
Arginine N-Methyltransferase Inhibitor-1 (AMI-1 (sodium salt))	38
Ash2L (human recombinant)	39
ATP-dependent helicase SMARCA4 (BRG1 bromodomain (human recombinant))	31
5-AzaC (5-Azacytidine)	39
5-aza-2'-Deoxycytidine (Decitabine)	39
5-Azacytidine	39
2',3',5'-triacetyl-5-Azacytidine	39

B

B 995 (Daminozide)	53
BAF180 (Polybromo-1D bromodomain 1 (human recombinant))	33
BAF190A Mitotic growth and transcription activator (BRG1 bromodomain (human recombinant))	31
BA2R (TAF1 bromodomain 1 (human recombinant))	33
BAZ2B bromodomain (human recombinant)	24
BAZ2B bromodomain TR-FRET Assay Kit	24
BIG3 (WDR5 (human recombinant))	52
BIX01294 (hydrochloride hydrate)	39
BML-210 (CAY10433)	13
BMP2-induced 3-kb Gene Protein (WDR5 (human recombinant))	52
BPTF bromodomain (human recombinant)	25
BR140-like protein (BRD1 bromodomain (human recombinant))	25
BRD and WD repeat-containing protein 1 (WDR9 bromodomain 2 (human recombinant))	33
BRD-containing protein testis specific	See BRDT
BRD1 bromodomain (human recombinant)	25
BRD2 bromodomain 1 (human recombinant; GST-tagged)	25
BRD2 bromodomain 1 TR-FRET Assay Kit	25
BRD2 bromodomains 1 and 2 (human recombinant)	25
BRD2 bromodomain 2 (human recombinant; GST-tagged)	25
BRD2 bromodomain 2 TR-FRET Assay Kit	28
BRD3 bromodomain 1 (human recombinant)	28
BRD3 bromodomain 1 TR-FRET Assay Kit	28
BRD4 bromodomain 1 (human recombinant; GST-tagged)	28
BRD4 bromodomain 1 (human recombinant; His-tagged)	29
BRD4 bromodomain 1 TR-FRET Assay Kit	29
BRD4 bromodomains 1 and 2 (human recombinant)	29
BRD4 bromodomain 2 (human recombinant; GST-tagged)	29
BRD4 bromodomain 2 (human recombinant; His-tagged)	29
BRD4 bromodomain 2 TR-FRET Assay Kit	30
BRD4/HUNK1 Polyclonal Antibody	30
BRD6	See BRDT
BRD9 bromodomain (human recombinant)	30
BRDT bromodomain 1 (human recombinant)	30
BRDT bromodomain 1 TR-FRET Assay Kit	30
BRDT bromodomain 2 (human recombinant)	31
BRG1-associated factor 180 (Polybromo-1D bromodomain 1 (human recombinant))	33

BRG1-associated factor 190A (BRG1 bromodomain (human recombinant))	31
BRG1 bromodomain (human recombinant)	31
BRL (BRD1 bromodomain (human recombinant))	25
BRM bromodomain (human recombinant)	31
BRM bromodomain TR-FRET Assay Kit	31
Bromodomain Adjacent to Zinc Finger Domain 2B	See BAZ2B bromodomain
Bromodomain and PHD Finger Containing 3 (BRPF3 bromodomain 1 (human recombinant))	31
Bromodomain and PHD finger-containing protein 2 (BRD1 bromodomain (human recombinant))	25
Bromodomain containing	See BRD
Bromodomain PHD Finger Transcription Factor (BPTF bromodomain (human recombinant))	25
BRPF1 (BRD1 bromodomain (human recombinant))	25
BRPF2 (BRD1 bromodomain (human recombinant))	25
BRPF3 bromodomain 1 (human recombinant)	31
Butyric Acid (Sodium Butyrate)	23
Butyrolactone 3	9

C

Camboginol (Garcinol)	10
cAMP-responsive element-binding protein 1 CREB-1 (CREB-binding protein bromodomain (human recombinant))	32
Cancer/testis antigen 9	See BRDT
CAP-C	See Chromosome Associated Protein-C Polyclonal Antibody
m-Carboxycinnamic Acid bis-Hydroxamide (CBHA)	13
CARM1	See PRMT4
Cat Eye Syndrome Critical Region Protein 2 (CECR2 bromodomain (human recombinant))	32
CAY10398	13
CAY10433	13
CAY10591	13
CAY10603	13
CBHA	13
CBP (CREB-binding protein bromodomain (human recombinant))	32
CBX5 (HP1- α (human recombinant))	36
CCG1 (TAF1 bromodomain 1 (human recombinant))	33
CCGS (TAF1 bromodomain 1 (human recombinant))	33
CECR2 bromodomain (human recombinant)	32
Chaetocin	39
Chidamide	13
Chromobox Homolog 5 (HP1- α (human recombinant))	36
Chromosome Associated Protein-C Polyclonal Antibody (aa 47-61)	6
Chromosome Associated Protein-C Polyclonal Antibody (aa 281-297)	6
CI-994	13
Click Tag™ 4-pentynoyl-CoA (4-pentynoyl-Coenzyme A (trifluoroacetate salt))	9
4-pentynoyl-Coenzyme A (trifluoroacetate salt)	9
Core Histones (human)	6
coumarin-SAHA	19
coumarin-Suberoylanilide Hydroxamic Acid (coumarin-SAHA)	19
cPLA(2) Interacting Protein (TIP60 (human recombinant))	10
CPTH2	9
CREB-binding protein bromodomain (human recombinant)	32
CREBBP (CREB-binding protein bromodomain (human recombinant))	32
CT9	See BRDT
CTG Repeat Protein 26 (NCOR2/SMRT (human recombinant))	19

D

D237 (M 344)	18
DAC (Decitabine)	39
Dacogen (Decitabine)	39
Daminozide	53
DDFC (Gemcitabine)	54
3-Deazaneplanocin A	39
3-Deazaneplanocin A (hydrochloride)	39
Decitabine	39
Delphinidin chloride	10
Demethylase (Jumonji-type) Activity Assay Kit	53
Demethylase (LSD-type) Activity Assay Kit	54
5-aza-2'-Deoxycytidine (Decitabine)	39

DIMG (Daminozide)	53
Disruptor of Telomeric Signaling 1-Like (Dot1L (human recombinant))	41
DKFZp686F0325 (BRD1 bromodomain (human recombinant))	25
DMAP1 (DNA Methyltransferase 1-Associated Protein 1 Polyclonal Antibody)	40
DMASA (Daminozide)	53
DNA Methylation EIA Kit	40
DNA Methyltransferase 1-Associated Protein 1 Polyclonal Antibody	40
DNA Methyltransferase 1 Monoclonal Antibody (Clone 60B1220.1)	40
DNA Methyltransferase 2 Monoclonal Antibody (Clone 102B1259.2)	40
DNA Methyltransferase 2 Polyclonal Antibody	40
DNA Methyltransferase 3a Monoclonal Antibody - Biotinylated (Clone 64B814.1)	40
DNA Methyltransferase 3a Monoclonal Antibody (Clone 64B1446)	40
DNA Methyltransferase 3a Monoclonal Antibody (Clone 64B814.1)	40
DNA Methyltransferase 3b Monoclonal Antibody (Clone 52A1018)	40
DNA Methyltransferase 3L (human recombinant)	40
DNMT3L (DNA Methyltransferase 3L (human recombinant))	40
DNMT3-Like Protein (DNA Methyltransferase 3L (human recombinant))	40
Dot1L (human recombinant)	41
Dot1-like (Dot1L (human recombinant))	41
DPY-30 (human recombinant)	41
DPY-30-like protein (DPY-30 (human recombinant))	41
DYT3 (TAF1 bromodomain 1 (human recombinant))	33
DZNep (3-Deazaneplanocin A)	39

E

E3 Ubiquitin-protein Ligase TRIM24 (TRIM24 bromodomain (human recombinant))	33
EED (human recombinant)	41
EHMT	See G9a
Ellagic Acid	41
Embryonic Ectoderm Development (EED (human recombinant))	41
Entinostat (MS-275)	18
Ephdine (Delphinidin chloride)	10
Epigenetics Screening Library (96-Well)	34
Esa1 (TIP60 (human recombinant))	10
Euchromatic Histone-Lysine N-Methyltransferase (G9a-like protein (human recombinant))	42
Eu HMTase 1 (G9a-like protein (human recombinant))	42
EX-527	13
EZH1 Polyclonal Antibody	41

F

Folfugem (Gemcitabine)	54
------------------------	----

G

G9a (human recombinant)	41
G9a-like protein (GLP SAM-Screener™ Assay Kit)	42
G9a-like protein (human recombinant)	42
G9a Methyltransferase Inhibitor Screening Assay Kit	42
Gallogen (Ellagic Acid)	41
Garcinol	10
GATA Zinc Finger Domain-Containing Protein 2A (p66 α Polyclonal Antibody)	19
GCK 1026 (Scriptaid)	20
Gcn5 (human recombinant)	10
Gemcel (Gemcitabine)	54
Gemcitabine	54
GemLip (Gemcitabine)	54
Gemzar (Gemcitabine)	54
General control of amino acid synthesis protein 5-like 2 (Gcn5 (human recombinant))	10
GLP (G9a-like protein (human recombinant))	42
GLP SAM-Screener™ Assay Kit	42
Goe 5549 (CI-994)	13
GSK-J1 (sodium salt)	54
GSK-J4 (hydrochloride)	54
GSK-J5 (hydrochloride)	54

H

H3.3	See Histone H3.3
H3F3	See Histone H3.3
H3K9me3 Peptide (Histone H3 Trimethyl Lys9 Peptide)	8
H3 Peptide (Histone H3 Peptide Substrate (1-21))	8
H4 Peptide (Histone H4 Peptide Substrate (1-21))	8
HAT (pCAF Histone Acetyltransferase)	10
HAT Inhibitor Screening Assay Kit	10
HC Toxin	13
(S)-HDAC-42	13
HDAC Activity Assay Kit	14
HDAC Cell-Based Activity Assay Kit	14
HDAC Inhibitor II (CBHA)	13
HDAC1 (human recombinant)	14
HDAC1 Inhibitor Screening Assay Kit	14
HDAC2 (human recombinant)	14
HDAC3/NCOR2 (human recombinant)	14
HDAC3 Polyclonal Antibody	15
HDAC4 (human recombinant)	15
HDAC4 Polyclonal Antibody	15
HDAC5 (human recombinant)	15
HDAC6 (human recombinant)	15
HDAC6 Polyclonal Antibody	15
HDAC7 (Phospho-Ser ¹⁵⁵) Polyclonal Antibody	15
HDAC8 (human recombinant)	15
HDAC8 Inhibitor Screening Assay Kit	18
HDAC9 (human recombinant)	18
HDAC11 Polyclonal Antibody	18
hDPY-30 (DPY-30 (human recombinant))	41
HeLa Core Histones (Core Histones (human))	6
Heterochromatin Protein 1- α (HP1- α (human recombinant))	36
Heterogeneous Nuclear Ribonucleoprotein Methyltransferase-like Protein 3 (PRMT3 (human recombinant))	45
HINFP (MBD2 Binding Zinc Finger Polyclonal Antibody)	36
Histone Acetyltransferase (HAT Inhibitor Screening Assay Kit)	10
Histone Arginine N-methyltransferase PRMT6	See PRMT6
Histone Deacetylase Inhibitor III (M 344)	18
Histone Deacetylase Inhibitor VI (HNHA)	18
Histone Deacetylase Inhibitor VII (Pimelic Diphenylamide 106)	19
Histone H2A Polyclonal Antibody	7
Histone H2A (<i>Xenopus</i> recombinant)	6
Histone H2B (C-Term) Polyclonal Antibody	7
Histone H2B (N-Term) Polyclonal Antibody	7
Histone H2B (<i>Xenopus</i> recombinant)	7
Histone H3 (human recombinant)	7
Histone H3 (Phospho-Ser ²⁹) Monoclonal Antibody (Clone 117C826)	7
Histone H3 Peptide Substrate (1-21)	8
Histone H3 Trimethyl Lys9 Peptide	8
Histone H3.3 Polyclonal Antibody	8
Histone H4 (human recombinant)	8
Histone H4K20 Peptide (Histone H4 Peptide Substrate (15-24))	8
Histone H4 Peptide Substrate (1-21)	8
Histone H4 Peptide Substrate (15-24)	8
Histone H4 Polyclonal Antibody	8
Histone H4 Transcription Factor (MBD2 Binding Zinc Finger Polyclonal Antibody)	36
Histone-lysine N-methyltransferase NSD2 isoform 1 (NSD2 (human recombinant))	45
HIV-1 Tat Interacting Protein 60 kDa (TIP60 (human recombinant))	10
HNHA	18
HP1- α (human recombinant)	36
hPB1 (Polybromo-1D bromodomain 1 (human recombinant))	33
HRMT1L (PRMT8 (human recombinant))	46
HRMT1L6	See PRMT6
HRX (MLL1 (human recombinant))	44
Hs.6364 (TIP60 (human recombinant))	10
HTATIP (TIP60 (human recombinant))	10
hTIF1 (TRIM24 bromodomain (human recombinant))	33
HUNK1	See BRD4
Huntington Interacting Protein B (SETD2 (human recombinant))	51
HYPB (SETD2 (human recombinant))	51

I2PP2A/SET Polyclonal Antibody	42
ILS-JGB-1741 (JGB1741)	18
Inhibitor of granzyme A-activated DNase (I2PP2A/SET Polyclonal Antibody)	42
4-iodo-SAHA	20
4-iodo-Suberoylanilide Hydroxamic Acid (4-iodo-SAHA)	20
IOX1	54

J

JBP1 (PRMT5 Polyclonal Antibody)	46
JGB1741	18
JHDM3A	See <i>JMJD2A</i>
JMJCD3C (JMJD2C-Strep tagged (human recombinant))	55
JMJD2A (human recombinant)	54
JMJD2A Inhibitor Screening Assay Kit	55
JMJD2A Polyclonal Antibody	55
JMJD2A-Strep tagged (human recombinant)	55
JMJD2C-Strep tagged (human recombinant)	55
JMJD2D (human recombinant)	55
JMJD2D Inhibitor Screening Assay Kit	56
JMJD2D-Strep tagged (human recombinant)	56
JMJD2E-Strep tagged (human recombinant)	56
JMJD6 (human recombinant)	56
JMJD6 Peptide Affinity-Purified Polyclonal Antibody	57
(+)-JQ1	32
(-)-JQ1	32
Jumonji Domain Containing	See <i>JMJD</i>

K

KAT2A (Gcn5 (human recombinant))	10
KAT2B (PCAF bromodomain (human recombinant))	32
KAT4 (TAF1 bromodomain 1 (human recombinant))	33
KAT5	See <i>TIP60</i>
KDM4A	See <i>JMJD2A</i>
KDM4C (JMJD2C-Strep tagged (human recombinant))	55
KDM4D	See <i>JMJD2D</i>
KDM6A (UTX (human recombinant))	60
KIAA0585 (JMJD6 (human recombinant))	56
KIAA0780 (JMJD2C-Strep tagged (human recombinant))	55
KIAA1090 (NSD2 (human recombinant))	45
KIAA1286 (BRPF3 bromodomain 1 (human recombinant))	31
KIAA1476	See <i>BAZ2B bromodomain</i>
KMT1C (G9a-like protein (human recombinant))	42
KMT2A	See <i>MLL1</i>
KMT3A (SETD2 (human recombinant))	51
KMT3B (NSD1 (human recombinant))	44
KMT3E (SMYD3 (human recombinant))	51
KMT5a (SET8 (human recombinant))	50
KMT5B (SUV4-20H1 (human recombinant))	51
KMT5C (SUV4-20H2 (human recombinant))	51
KMT7	See <i>SET7/9</i>
KMT8 (RIZ1 (human recombinant))	47
KRS2 (Mammalian STE-20-Like Kinase 1 Polyclonal Antibody)	8
Kylar (Daminozide)	53

L

Ladakamycin (5-Azacytidine)	39
Lagistase (Ellagic Acid)	41
LSD1 Inhibitor Screening Assay Kit	57
LSD1 Polyclonal Antibody (aa 100-150)	57
LSD1 Polyclonal Antibody (aa 400-450)	57
LSD1 Polyclonal Antibody (aa 450-500)	57
LSD1 Polyclonal Antibody (aa 800-850)	60
LY 188011 (Gemicitabine)	54
Lysine acetyltransferase 2A (Gcn5 (human recombinant))	10
Lysine Acetyltransferase 5 (TIP60 Polyclonal Antibody)	11
Lysine Methyltransferase 2A (MLL1 (human recombinant))	44
Lysine N-methyltransferase 5B (SUV4-20H1 (human recombinant))	51
Lysine-Specific Demethylase 1 (LSD1 Inhibitor Screening Assay Kit)	57
Lysine-Specific Demethylase 4	See <i>JMJD2</i>

M

M 344	18
Mammalian STE-20-Like Kinase 1 Polyclonal Antibody	8
MBD1 Monoclonal Antibody (Clone 100B272.1)	36
MBD1 Polyclonal Antibody	36
MBD2 (human recombinant; methyl binding domain aa 150-220)	36
MBD2 Binding Zinc Finger Polyclonal Antibody	36
MBD2/3 Monoclonal Antibody (Clone 106B691)	37
MCAP	See <i>BRD4</i>
MD 85 (CAY10398)	13
MeCP2 (human recombinant; methyl binding domain aa 77-166)	37
MeCP2 Polyclonal Antibody	37
Meisetz (PRDM9 (human recombinant))	45
Metacept 3 (Oxamflatin)	19
Metastasis Associated 1 Family Member 2 Polyclonal Antibody	18
Methylated Lysine Polyclonal Antibody	42
Methylated Lysine Polyclonal Antibody-biotin	43
Methylated Lysine Polyclonal Antibody HRP Conjugate	43
Methyl-CpG-Binding Domain	See <i>MBD</i>
Methyl-CpG Binding Protein 2	See <i>MeCP2</i>
Methyl Cytosine Binding Domain 2 (MBD2 (human recombinant; methyl binding domain aa 150-220))	36
Methyltransferase Colorimetric Assay Kit	43
Methyltransferase Fluorometric Assay Kit	43
MGC58603 (BRPF3 bromodomain 1 (human recombinant))	31
MGC165155 (Polybromo-1D bromodomain 1 (human recombinant))	33
MGC165156 (Polybromo-1D bromodomain 1 (human recombinant))	33
Ml-2 (hydrochloride)	44
Ml-nc (hydrochloride)	44
Mitochondrial Nicotinamide Adenine Dinuclear-dependent Deacetylase (SIRT3 (human recombinant))	22
MIZF (MBD2 Binding Zinc Finger Polyclonal Antibody)	36
Mixed-Lineage Leukemia-1	See <i>MLL1</i>
MLL Core Complex	See <i>MLL1</i>
MLL1 (human recombinant)	44
MLL1 SAM-Screener™ Assay Kit	43
MLL1/WAR complex (human recombinant)	44
MLL1/WARD complex (human recombinant)	44
MMSSET (NSD2 (human recombinant))	45
MS-275	18
MS 344 (M 344)	18
MSBP3 (PRDM9 (human recombinant))	45
MST-1 (Mammalian STE-20-Like Kinase 1 Polyclonal Antibody)	8
MT	See <i>Methyltransferase</i>
MTA-L1 Protein (Metastasis Associated 1 Family Member 2 Polyclonal Antibody)	18
MTA2 (Metastasis Associated 1 Family Member 2 Polyclonal Antibody)	18
Mylosar (5-Azacytidine)	39

N

NAD-dependent ADP-ribosyltransferase Sirtuin 4 (SIRT4 (human recombinant))	22
NAD-dependent Deacetylase	See <i>SIRT</i>
1-Naphthalenepropanoic Acid (Splitomicin)	23
NCOR2/SMRT (human recombinant)	19
NOG (N-Oxalylglycine)	60
NSC 103-627 (5-Azacytidine)	39
NSC 617989 (3-Deazaneplanocin A)	39
NSC 102816 (5-Azacytidine)	39
NSC 127716 (Decitabine)	39
NSC 613327 (Gemicitabine)	54
NSCL2 (TAF1 bromodomain 1 (human recombinant))	33
NSD1 (human recombinant)	44
NSD2 (human recombinant)	45
Nuclear Receptor Corepressor 2 (NCOR2/SMRT (human recombinant))	19
Nuclear Receptor SET domain-containing protein	See <i>NSD</i>

O

ORFX	See <i>BRD3</i>
OSI 2040 (Apicidin)	12
N-Oxalylglycine	60
Oxamflatin	19

P

p53 Target Protein in Deacetylase Complex (Metastasis Associated 1 Family Member 2 Polyclonal Antibody)	18
p66α Polyclonal Antibody	19
p231HBP (SETD2 (human recombinant))	51
p300/CBP-associated factor (PCAF bromodomain (human recombinant))	32
p300/(CREB binding protein) Associated Factor (pCAF Histone Acetyltransferase)	10
PAD4 (human recombinant)	60
PAD4 Inhibitor Screening Assay Kit	60
PB1 (Polybromo-1D bromodomain 1 (human recombinant))	33
PBRM1 (Polybromo-1D bromodomain 1 (human recombinant))	33
PCAF bromodomain (human recombinant)	32
pCAF Histone Acetyltransferase	10
PCI 34051	19
PD-PCPA (hydrochloride)	60
PD 123654 (CI-994)	13
Peptidylarginine Deiminase 4	See <i>PAD4</i>
6-pentadecyl Salicylic Acid (Anacardic Acid)	9
4-pentynoyl-Coenzyme A (trifluoroacetate salt)	9
PFI-1	32
PFM6 (PRDM9 (human recombinant))	45
PHAPII (I2PP2A/SET Polyclonal Antibody)	42
trans-2-Phenylcyclopropylamine (2-PCPA (hydrochloride))	60
N-phenyl-N'-(2-Aminophenyl)hexamethylenediamide (CAY10433)	13
PHH3 (Histone H3 (Phospho-Ser ²⁸) Monoclonal Antibody (Clone 117C826))	7
Phosphatase 2A Inhibitor (I2PP2A/SET Polyclonal Antibody)	42
N-Phthalyl-L-Tryptophan (RG-108)	47
PID (Metastasis Associated 1 Family Member 2 Polyclonal Antibody)	18
Pimelic Diphenylamide 106	19
PLIP (TIP60 (human recombinant))	10
Polybromo-1D bromodomain 1 (human recombinant)	33
Pracinostat (SB 939)	20
PRDM2 (RIZ1 (human recombinant))	47
PRDM9 (human recombinant)	45
PR Domain-containing Protein 9 (PRDM9 (human recombinant))	45
PR Domain Zinc Finger Protein 2 (RIZ1 (human recombinant))	47
PRMT1 (human recombinant)	45
PRMT3 (human recombinant)	45
PRMT4 (human recombinant)	45
PRMT4 Polyclonal Antibody	45
PRMT5 Polyclonal Antibody	46
PRMT6 (human recombinant)	46
PRMT6 (human recombinant; baculovirus expressed)	46
PRMT6 Polyclonal Antibody	46
PRMT7 Polyclonal Antibody	46
PRMT8 (human recombinant)	46
2-Propylvaleric Acid (Valproic Acid (sodium salt))	23
Protein Arginine Deiminase 4	See <i>PAD4</i>
Protein Arginine Methyltransferase	See <i>PRMT</i>
Protein BRG-1 (BRG1 bromodomain (human recombinant))	31
Protein Polybromo-1 (Polybromo-1D bromodomain 1 (human recombinant))	33
Protein trithorax-5 (NSD2 (human recombinant))	45
PR-Set7 (SET8 (human recombinant))	50
PTC6 (TRIM24 bromodomain (human recombinant))	33
PTDSR	See <i>JMJD6</i>
PX 089274 (CAY10398)	13
Pyroxamide	19

Q

R

RbBP5 (human recombinant)	47
RBQ-3 (RbBP5 (human recombinant))	47
Retinoblastoma-binding Protein 5 (RbBP5 (human recombinant))	47
Retinoblastoma Protein-interacting Zinc Finger Protein (RIZ1 (human recombinant))	47
RG-108	47
RING finger protein 82 (TRIM24 bromodomain (human recombinant))	33

RING3	See <i>BRD and BRDT</i>
RIZ1 (human recombinant)	47
RNF3	See <i>BRD2</i>
RNF82 (TRIM24 bromodomain (human recombinant))	33

S

SADH (Daminozide)	53
SAF19 (DPY-30 (human recombinant))	41
SAHA	19
SAHA-BPyne	19
coumarin-SAHA	19
4-iodo-SAHA	20
Salemide	20
6-pentadecyl Salicylic Acid (Anacardic Acid)	9
SAM (S-(5'-Adenosyl)-L-methionine chloride (hydrochloride))	38
SB 939	20
SBHA (Suberohydroxamic Acid)	23
Scriptaid	20
SET and MYND domain-containing protein	See <i>SMYD</i>
Set1/Ash2 Histone Methyltransferase Complex Subunit Ash2 Isoform A (Ash2L (human recombinant))	39
SET7/9 (human recombinant)	47
SET7/9 Methyltransferase Inhibitor Screening Assay Kit	47
SET7/9 Polyclonal Antibody	47
SET7/9 (FL) Polyclonal Antibody	50
SET7/9 SAM-Binding Site Inhibitor Screening (SET7/9 SAM-Screener™ Assay Kit)	50
SET7/9 SAM-Screener™ Assay Kit	50
SET8 (human recombinant)	50
SET8 Methyltransferase Acceptor Peptide (Histone H4 Peptide Substrate (15-24))	8
SET8 Methyltransferase Inhibitor Screening Assay Kit	50
SETD2 (human recombinant)	51
SETD7/9	See <i>SET7/9</i>
SETD8 (SET8 (human recombinant))	50
SET Domain-Containing Protein 7/9	See <i>SET7/9</i>
SET domain-containing (lysine methyltransferase) 8 (SET8 (human recombinant))	50
Silencing Mediator of Retinoic Acid and Thyroid Hormone Receptor (NCOR2/SMRT (human recombinant))	19
Silent Information Regulator	See <i>SIRT</i>
SIR2L	See <i>SIRT</i>
SIR2-like Protein	See <i>SIRT</i>
SIRT1 Activator 3 (CAY10591)	13
SIRT1 Direct Fluorescent Screening Assay Kit	20
SIRT1 FRET-Based Screening Assay Kit	21
SIRT1 (human recombinant)	21
SIRT2 Direct Fluorescent Screening Assay Kit	20
SIRT2 (human recombinant)	21
SIRT3 Direct Fluorescent Screening Assay Kit	21
SIRT3 (human recombinant)	22
SIRT4 (human recombinant)	22
SIRT5 (human recombinant)	22
SIRT6 Direct Fluorescent Screening Assay Kit	21
SIRT6 (human recombinant)	22
SIRT7 (human recombinant)	22
SIRT7 Polyclonal Antibody	22
Sirtinol	23
Sirtuin	See <i>SIRT</i>
Sirtuin 1 Activator 3 (CAY10591)	13
Sir Two Inhibitor Naphthol (Sirtinol)	23
Skb1HS (PRMT5 Polyclonal Antibody)	46
SMAP270 (NCOR2/SMRT (human recombinant))	19
SMARCA2A/B (BRM bromodomain (human recombinant))	31
smk-BOP (SMYD1 (human recombinant))	51
SMRT (NCOR2/SMRT (human recombinant))	19
SMYD1 (human recombinant)	51
SMYD3 (human recombinant)	51
SNDX 275 (MS-275)	18
Sodium Butyrate	23
Sodium Valproate (Valproic Acid (sodium salt))	23
Splitomicin	23
SSET1c WD40 Repeat Protein (WDR5 (human recombinant))	52
STAF97 (Gcn5 (human recombinant))	10

STK4 (Mammalian STE-20-Like Kinase 1 Polyclonal Antibody)	8
Suberic bis-Hydroxamic Acid (Suberohydroxamic Acid)	23
Suberohydroxamic Acid	23
Suberoylanilide Hydroxamic Acid	See SAHA
Succinic Acid (Daminozide)	53
Suppressor of Variegation 4-20 Homolog	See SUV-20H
SUV4-20H1 (human recombinant)	51
SUV4-20H2 (human recombinant)	51
SWD1/Set1c WD40 repeat protein homolog (RbBP5 (human recombinant))	47
SWD3 (WDR5 (human recombinant))	52
SWI/SNF ATPase (BRM bromodomain (human recombinant))	31

T

T3 Receptor-Associating Factor (NCOR2/SMRT (human recombinant))	19
TAB-Associated Factor 10 (TAF 10 Peptide)	52
Tacedinaline (CI-994)	13
TAF1 bromodomain 1 (human recombinant)	33
TAF1 RNA Polymerase II (TAF1 bromodomain 1 (human recombinant))	33
TAF2A (TAF1 bromodomain 1 (human recombinant))	33
TAF 10 Peptide	52
TAF10 RNA polymerase II (TAF 10 Peptide)	52
TAFII30 (TAF 10 Peptide)	52
TAFII250 (TAF1 bromodomain 1 (human recombinant))	33
TATA Box Binding Protein (TBP)-Associated Factor	See TAF
TBBD (Ellagic Acid)	41
TC-H 106 (Pimelic Diphenylamide 106)	19
Template-Activating Factor I (I2PP2A/SET Polyclonal Antibody)	42
Tenovin-1	23
Tenovin-6	23
Thyroid-Retinoic-Acid-Receptor-Associated Corepressor (NCOR2/SMRT (human recombinant))	19
TIP60 (human recombinant)	10
TIP60 Polyclonal Antibody	11
Toxin I (<i>Helminthosporium carbonum</i>) (HC Toxin)	13
TRAC (NCOR2/SMRT (human recombinant))	19
Transcriptional Intermediary Factor 1- α (TRIM24 bromodomain (human recombinant))	33
Transcriptional Repressor p66 α (p66 α Polyclonal Antibody)	19
Tranylcypromine (2-PCPA (hydrochloride))	60
2',3',5'-triacetyl-5-Azacytidine	39
Trichostatin A	23
TRIM24 bromodomain (human recombinant)	33
Tripartite Motif Containing 24 (TRIM24 bromodomain (human recombinant))	33
TRX5 (NSD2 (human recombinant))	45
TSA (Trichostatin A)	23
Tubastatin A (trifluoroacetate salt)	23

U

Ubiquitously Transcribed Tetratricopeptide Repeat X (UTX (human recombinant))	60
UNC0224	52
UNC0321 (trifluoroacetate salt)	52
UNC0638	52
UNC1215	37
UTX (human recombinant)	60

V

Valproic Acid (sodium salt)	23
Vorinostat (SAHA)	19

W

WALp4	See BAZ2B bromodomain
WDR5 (human recombinant)	52
WDR9 bromodomain 2 (human recombinant)	33
WD repeat-containing protein 9 (WDR9 bromodomain 2 (human recombinant))	33
WD-Repeat Protein 5 (WDR5 (human recombinant))	52
WHSC1	See NSD
WR 183027 (5-Azacytidine)	39

X

XDP (TAF1 bromodomain 1 (human recombinant))	33
--	----

Y

Z

Zefei (Gemcitabine)	54
ZNF899 (PRDM9 (human recombinant))	45
Zolinza™ (SAHA)	19

Item No. Index

10228	52	11285	28	13559	46
10261	6	11286	36	13631	52
10262	7	11287	37	13686	13
10263	7	11288	32	13725	9
10264	8	11289	31	13726	9
10315	22	11299	55	13727	42
10316	22	11300	56	13728	43
10317	22	11507	25	13729	43
10318	22	11509	30	13731	47
10319	50	11548	30	13771	36
10320	47	11549	33	13772	36
10335	55	11572	54	13773	37
10336	54	11620	44	13775	37
10350	45	11621	44	13776	8
10353	41	11633	19	13777	36
10354	41	11642	45	13778	18
10380	8	11644	46	13780	50
10382	55	11649	31	13782	42
10443	20	11650	31	13784	8
10444	19	11652	33	13785	19
10495	20	11653	33	13787	57
10500	60	11690	54	13789	11
10523	23	11720	29	13828	39
10530	8	11721	29	13866	46
10547	9	11917	24	13870	19
10559	23	11918	25	13944	60
10566	10	11920	32	13956	38
10569	41	11922	33	13965	38
10572	20	12033	53	13968	37
10574	23	12054	54	14004	12
10575	12	12073	54	14133	32
10576	13	12074	54	19380	15
10582	52	12084	13	89730	23
10599	53	12086	9	89740	13
10610	53	12095	9	589324	40
10628	41	13033	23	600150	14
10641	18	13085	23	600490	50
10658	44	13086	23	600500	25
10671	19	13121	23	600510	28
10675	19	13124	39	600520	29
10734	52	13144	9	600530	30
10750	45	13145	12	600570	42
10752	46	13146	13	600580	43
10755	42	13156	39	600630	28
10756	44	13168	23	600650	30
10757	44	13172	13	600710	24
10758	45	13174	18	600730	31
10761	51	13176	19	700120	57
10762	51	13178	20	700140	43
10763	51	13212	19	700150	43
10764	51	13277	13	700230	18
10765	47	13284	18	700270	47
10767	51	13295	18	700280	20
10770	40	13302	47	700290	21
10773	56	13373	39	700350	50
10774	60	13477	22	700360	55
10776	55	13479	40	700370	56
10782	10	13480	40	700390	53
10783	10	13481	40	700400	54
10854	8	13482	40	700500	42
10877	8	13483	40	700560	60
10944	52	13484	40	10005019	12
10945	44	13485	40	10006515	10
10946	39	13486	57	10009115	10
10947	47	13487	41	10009231	14
11010	6	13493	15	10009232	14
11012	10	13494	15	10009377	14
11052	29	13497	30	10009379	15
11066	29	13499	15	10009465	15
11068	28	13500	15	10009466	18
11069	25	13501	6	10009652	15
11070	25	13503	6	10009797	13
11071	25	13504	18	10009798	13
11076	34	13535	7	10009929	19
11102	39	13536	40	10010401	20
11155	32	13538	7	10010494	60
11164	39	13539	7	10010567	9
11166	39	13540	7	10010991	21
11178	41	13543	8	10011563	14
11187	32	13551	46	10011564	14
11209	45	13552	45	10011566	21
11232	32	13553	57	10011190	21
11235	36	13554	57	10011191	21
11237	56	13555	60	10011194	22
11284	31	13558	46		



(800) 364 9897
www.caymanchem.com