 CAPP APS  
LANDBRUGSVEJ 10  
DK-5260 ODENSE  
DENMARK

TEL. (+45) 6613 6140  
FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK  
WEB: WWW.CAPP.DK

**CAPP**  
D E N M A R K

# Capp *Tempo*

PIPETTE CONTROLLER

The Capp Tempo pipette controller has revolutionary features that separate it from other pipette controllers. The Capp Tempo has an LED display for ultimate user friendliness and the patented Electronic Break System (EBS).

The LED display gives the user relevant information about aspiration speed, dispense speed and battery life. The EBS is unique to Capp Tempo and allows users to pipette at maximum speed without compromising accuracy.

The Capp Tempo is loaded with features in a stylish Danish design and is available in four color choices.

The charging stand and wall holder are included in the box with the Capp Tempo pipette controller.

Capp Tempo gives you everything you need to accomplish your goals.

#### LED Display

Display shows aspiration and dispense speeds

Displays Remaining battery life and flashes when power is low

Display turns-off after 20 seconds of non-use to conserve battery

#### Adjustable speeds

Powerful motor fills 25ml in under 12 seconds

Separate speed controls for aspiration and dispense modes

Six different speed settings for both aspiration and dispense modes

#### Electronic Break System (EBS)

EBS maximizes speed and precision

Apply EBS to slow down at the last second to dial in the most accurate volume


Press EBS a second time to return to max speed

#### Power Supply

Lithium ion battery, full charge provides 8 hours of continuous use

Battery can be charged while pipette is being used





# Ordering Info



## CappTempo

| Cat No. | Description      |
|---------|------------------|
| T100-B  | CappTempo Blue   |
| T100-G  | CappTempo Green  |
| T100-R  | CappTempo Red    |
| T100-V  | CappTempo Violet |