



---

## CERTIFICATE OF ANALYSIS

Product **AGAROSE PHASEOLUS VULGARIS LEUCOAGGLUTININ (PHA-L)**

Catalog No. AL-1113

Amount 2 ml of settled gel (1:1, v/v slurry)

Lot No. W0728

Suspension Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.1 mM Ca<sup>++</sup>,  
0.08% sodium azide, 0.01 mM Mn<sup>++</sup>

Protein Concentration 3.0 mg/ml of settled gel

Binding Capacity >4.0 mg thyroglobulin/ml of gel

Storage Conditions Refrigerate - DO NOT FREEZE

### Other Remarks

Recommended product for eluting glycoconjugates bound to this agarose-lectin:  
Glycoprotein Eluting Solution, Cat. No. ES-2100.  
Alternatively, 100 mM acetic acid or 50 mM glycine, pH 3.0 can be used.

After use, wash the gel with several column volumes of buffered saline, then  
resuspend gel in buffered saline containing 0.08% sodium azide for storage.