



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED CONCANAVALIN A**

Catalog **B-1005**

Amount **5 mg**

Concentration **5 mg/ml**

Lot No. **Z1203**

Active Conjugate **98%**

Lectin Biotinylated **99%**

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.01 mM Mn⁺⁺, 0.1 mM Ca⁺⁺,
0.08% sodium azide.

Storage Conditions **Refrigerate; for long term storage, aliquots may be stored frozen**

For most applications we recommend a freshly prepared working solution of 5-20 $\mu\text{g/ml}$ in the above buffer