



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED SOYBEAN AGGLUTININ**

Catalog B-1015

Amount 5 mg

Concentration 5 mg/ml

Lot No. ZA0512

Active Conjugate 99%

Lectin Biotinylated 99%

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺

Storage Conditions Refrigerate; for long term storage, aliquot and store frozen

For most applications we recommend a freshly prepared working solution of 5-20 μ g/ml in the above buffer