



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *ULEX EUROPAEUS* AGGLUTININ I**

Catalog B-1065

Amount 2 mg active conjugate

Concentration 2 mg/ml

Lot No. ZC0331

Active Conjugate 98%

Lectin Biotinylated 97%

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.1 mM Ca⁺⁺, 0.08% sodium azide.

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen.

Remarks:

For most applications we recommend a freshly prepared working solution of 5-20 µg/ml in the above buffer.