



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *GRIFFONIA (BANDEIRAEA)***
***SIMPLICIFOLIA* LECTIN I**
(GSL I, BSL I)

Catalog **B-1105**

Amount **2 mg**

Concentration **2 mg/ml**

Lot No. **Z1022**

Active Conjugate **99%**

Lectin Biotinylated **99%**

Solution **10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide,**
0.1 mM Ca⁺⁺.

Storage Conditions **Refrigerate; for long term storage, aliquots may be stored frozen.**

Remarks

For most applications we recommend a freshly prepared working solution of 5-20 μ g/ml in the above buffer.