



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED PHASEOLUS VULGARIS ERYTHROAGGLUTININ (PHA-E)**

Catalog B-1125

Amount 2 mg

Concentration 2 mg/ml

Lot No. ZA0214

Active Conjugate 98%

Lectin Biotinylated 98%

Storage Conditions 4 °C

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.1 mM Ca⁺⁺, 0.08% sodium azide.

For most applications we recommend a freshly prepared working solution of 10 μ g/ml in the above buffer