



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *LYCOPERSICON ESCULENTUM***
(TOMATO) LECTIN

Catalog B-1175

Amount 1 mg active conjugate

Lot No. W0424

Active Conjugate 99%

Lectin Biotinylated 99%

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen.

Other Remarks:

Reconstitute by adding 0.5 ml of water. The resulting solution will have the following composition: 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 0.1 mM Ca⁺⁺.

For most applications, we recommend a freshly prepared working solution of 5-20 $\mu\text{g/ml}$ in the above buffer.