



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED *HIPPEASTRUM* HYBRID (AMARYLLIS) LECTIN**

Catalog B-1385

Amount 2 mg

Concentration 2 mg/ml

Lot No. V0408

Active Conjugate 99%

Lectin Biotinylated 99%

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide

Storage Refrigerate

Other Remarks:

For most applications, we recommend a freshly prepared working solution of 5-20 $\mu\text{g/ml}$ in the above buffer.