



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED ALEURIA AURANTIA LECTIN**

Catalog No. B-1395

Amount 1 mg active conjugate

Concentration 2 mg/ml

Lot No. Z0415

Active Conjugate 99%

Lectin Biotinylated 99%

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08%

Storage Conditions Refrigerate; for long term storage, aliquots may be stored frozen

Other Remarks:

Binding inhibitors: 100-400 mM L-fucose