



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED ANTI-FLUORESCIN**
AFFINITY PURIFIED

Catalog BA-0601

Amount 0.5 mg

Concentration 1 mg/ml

Lot No. W1203

Produced in goat

Solution 10 mM phosphate, 0.15 M NaCl, pH 7.8, 0.08% sodium azide.
Immunohistochemical grade bovine serum albumin has been added for stabilization.

Storage Refrigerate

Biotinylation 98%

Antibody activity after biotinylation (by solid-phase binding assay) 96%

Other Remarks:

The recommended concentration range for use is 2 - 10 µg/ml.