



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED ANTI-DINITROPHENYL (DNP)**
AFFINITY PURIFIED

Catalog **BA-0603**

Amount **0.5 mg**

Concentration **0.5mg/ml**

Lot No. **ZC0708**

Produced in **rabbit**

Solution 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide. Sucrose has been added for stabilization.

Storage **Refrigerate**

Antibody biotinylated **97%**

Other Remarks

The recommended concentration range for use is 2-10 µg/ml.