



CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED ANTI-RHODAMINE**
AFFINITY PURIFIED

Catalog BA-0605

Amount 0.5 mg

Lot No. V0317

Produced in goat

Solution Reconstitute by adding 0.5 ml water. The resulting solution will have the following composition: 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide. Immunohistochemical grade bovine serum albumin has been added for stabilization.

Storage Refrigerate

Other Remarks:

This antibody conjugate has a broad specificity toward various rhodamines and recognizes most, if not all rhodamine related fluorochromes, including Texas Red[®], tetramethylrhodamine and rhodamine₆₀₀ (x-rhodamine). The recommended concentration range for most applications is 2-10 µg/ml.

Texas Red[®] is a trademark of Molecular Probes