



CERTIFICATE OF ANALYSIS

Product **DyLight® 488 ANTI-DIGOXIGENIN/DIGOXIN**
****AFFINITY PURIFIED****

Catalog No. DI-7488

Amount 0.5 mg active conjugate

Lot No. ZA0304

Produced in goat

Concentration 1.0 mg/ml active conjugate

Solution 10 mM Phosphate, 0.15 M NaCl, pH 7.8, 0.08% sodium azide.


Maximum Excitation @ 493 nm
Maximum Emission @ 518 nm

F/P (molar) 5.0

Storage Conditions Refrigerate in the dark.

The recommended concentration range for use is 5-20 μ g/ml.

DyLight® 488 Conjugate Spectral Properties



DyLight® is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.