



CERTIFICATE OF ANALYSIS

Product **FLUORESCCEIN *DOLICHOS BIFLORUS* AGGLUTININ**

Catalog No. FL-1031

Amount 2 mg active conjugate

Lot No. ZA0722

Concentration 2 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺

Free Fluorescein None

F/P (molar) 2.2

Active Conjugate 97%

Storage Conditions Refrigerate in the dark

Other Remarks:

A precipitate may form during storage. This does not have a significant adverse effect on the product. If a precipitate forms upon long-term storage at 4 °C, warm to 37 °C and centrifuge before use.

The recommended concentration range for use is 5-20 µg/ml.