



CERTIFICATE OF ANALYSIS

Product FLUORESCEIN *LENS CULINARIS* AGGLUTININ

Catalog No. EL-1041

Amount 5 mg

Lot No. ZD0104

Concentration 5 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺, 0.01 mM Mn⁺⁺.

Free Fluorescein None

F/P (molar) 3.5

Active Conjugate 95%

Storage Conditions Refrigerate in the dark

Other Remarks:

A precipitate may form during storage. This does not have a significant adverse effect on the product. If a precipitate forms upon long-term storage at 4 °C, warm to 37 °C and centrifuge before use.

The recommended concentration range for use is 5-20 µg/ml.