



CERTIFICATE OF ANALYSIS

Product **FLUORESCHEIN ULEX EUROPAEUS AGGLUTININ I**

Catalog No. EL-1061

Amount 2 mg active conjugate

Lot No. ZC0426

Concentration 2 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺

Free Fluorescein None

F/P (molar) 2.6

Active Conjugate 98%

Storage Conditions Refrigerate in the dark

Other Remarks:

The recommended concentration range for use is 5-20 µg/ml.