



CERTIFICATE OF ANALYSIS

Product **FLUORESCIN *PHASEOLUS VULGARIS***
LEUCOAGGLUTININ (PHA-L)

Catalog No. EL-1111

Amount 2 mg active conjugate

Lot No. ZC0302

Concentration 2 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺

Free Fluorescein None

F/P (molar) 4.9

Maximum Excitation @ 495-500 nm
Maximum Emission @ 514-521 nm

Storage Conditions Refrigerate in the dark

Other Remarks:

The recommended concentration range for use is 5-20 µg/ml.