



CERTIFICATE OF ANALYSIS

Product **FLUORESCHEIN *PHASEOLUS VULGARIS***
ERYTHROAGGLUTININ (PHA-E)

Catalog No. EL-1121

Amount 2 mg protein

Lot No. ZD0118

Concentration 2 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide, 0.1 mM Ca⁺⁺,
5mg/ml Beta-Cyclodextrin

Free Fluorescein None

F/P (molar) 3.4

Storage Conditions Refrigerate in the dark

Other Remarks:

The recommended concentration range for use is 5-20 µg/ml.