



CERTIFICATE OF ANALYSIS

Product FLUORESCEIN SOLANUM TUBEROSUM (POTATO) LECTIN

Catalog No. EL-1161

Amount 2 mg

Lot No. Zc0304

Concentration 2 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺

Free Fluorescein None

F/P (molar) 3.6

Active Conjugate 99%

Storage Conditions Refrigerate in the dark

Other Remarks:

If a precipitate forms upon long-term storage, warm to 37 °C.