



CERTIFICATE OF ANALYSIS

Product FLUORESCHEIN LYCOPERSICON ESCULENTUM
(TOMATO) LECTIN

Catalog No. EL-1171

Amount 1 mg active conjugate

Lot No. ZA0812

Concentration 2 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺, 5mg/ml β cyclodextrin

Free Fluorescein None

F/P (molar) 2.1

Active Conjugate 99%

Storage Conditions Refrigerate in the dark

Other Remarks:

If a precipitate forms upon long-term storage, warm to 37 °C.

The recommended concentration range for use is 5-20 μ g/ml.