



CERTIFICATE OF ANALYSIS

Product FLUORESCIN *GRIFFONIA (BANDEIRAEA) SIMPLICIFOLIA*
LECTIN II

Catalog No. EL-1211

Amount 2 mg active conjugate

Lot No. ZC0601

Concentration 2.0 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺

Free Fluorescein None

F/P (molar) 3.6

Storage Conditions Refrigerate in the dark

Other Remarks:

The recommended concentration range for use is 5-20 µg/ml.