



CERTIFICATE OF ANALYSIS

Product FLUORESCHEIN *VICIA VILLOSA* LECTIN

Catalog No. EL-1231

Amount 2 mg

Lot No. ZB0501

Concentration 2 mg/ml

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide,
0.1 mM Ca⁺⁺

Free Fluorescein None

F/P (molar) 3.5

Maximum Excitation @ 495-500 nm
Maximum Emission @ 514-521 nm

Active Conjugate 94%

Storage Conditions Refrigerate in the dark

Other Remarks:

If a precipitate forms upon long-term storage, warm to 37 °C.

The recommended concentration range for use is 5-20 µg/ml.