



CERTIFICATE OF ANALYSIS

Product FLUORESCHEIN *ALEURIA AURANTIA* LECTIN
AFFINITY PURIFIED

Catalog No. EL-1391

Amount 1 mg active conjugate

Concentration 1 mg/ml

Lot No. ZC1115

Solution 10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide.

Active Conjugate 99%

Free Fluorescein None

F/P (molar) 2.2

Storage Conditions Refrigerate in the dark

Other Remarks:

The recommended concentration range for use is 5-20 µg/ml.

The recommended mounting medium is VECTASHIELD® Hard+Set™ (Cat. No. H-1400) or VECTASHIELD® Hard+Set™ with DAPI (Cat. No. H-1500).

Binding is inhibited by 100 mM-400 mM L-fucose.