



CERTIFICATE OF ANALYSIS

Product **SOYBEAN AGGLUTININ (SBA)**

Catalog No. L-1010

Amount 10 mg of lyophilized lectin

Lot No. ZC0916

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

Although many buffers can be employed for reconstituting and diluting this lectin, we recommend 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca⁺⁺. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.