



## CERTIFICATE OF ANALYSIS

Product **WHEAT GERM AGGLUTININ (WGA)**

Catalog No. L-1020

Amount 25 mg of lyophilized lectin

Lot No. X1205

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

Salt-free.

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca<sup>++</sup> is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.