



CERTIFICATE OF ANALYSIS

Product **SUCCINYLATED WHEAT GERM AGGLUTININ (WGA)**

Catalog No. L-1020S

Amount 10 mg salt-free

Lot No. W0110

Storage Refrigerate

Other remarks:

This derivative of WGA was prepared from highly purified WGA which tested as homogeneous by SDS-PAGE.

Monsigny, et al, have reported (Eur. J. Biochem., 98, 39-45, 1979) that, in contrast to the navit lectin, succinylated-WGA does not recognize sialic acid residues. However, succinylated WGA has the same affinity for GlcNAc as the native lectin.