



## CERTIFICATE OF ANALYSIS

Product **ULEX EUROPAEUS AGGLUTININ I**

Catalog No. L-1060

Amount 2 mg of lyophilized lectin, salt-free

Lot No. ZA0218

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

This lot will agglutinate a 3% suspension of fresh human type O red blood cells at a lectin concentration of 2.0  $\mu\text{g/ml}$ .

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 7.5 - 8.5, 0.1 mM  $\text{Ca}^{++}$  is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.