



CERTIFICATE OF ANALYSIS

Product *ULEX EUROPAEUS* AGGLUTININ I

Catalog No. L-1060

Amount 5 mg of lyophilized lectin, salt-free

Lot No. U0920

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

This lot will agglutinate a 3% suspension of fresh human type O red blood cells at a lectin concentration of 2.0 µg/ml.

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 7.5 - 8.5, 0.1 mM Ca⁺⁺ is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.