



## CERTIFICATE OF ANALYSIS

Product **PEANUT AGGLUTININ**

Catalog No. L-1070

Amount 10 mg of lyophilized lectin. In addition to the 10 mg of lectin, the lyophilized powder contains approximately 37% salts.

Lot No. Y0207

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca<sup>++</sup> is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.