



---

## CERTIFICATE OF ANALYSIS

Product **RICINUS COMMUNIS AGGLUTININ I (RCA 120)**

Catalog No. L-1080

Amount 10mg

Lot No. ZB1022

Protein Concentration 5.0 mg/ml

Solution 10 mM phosphate, 0.15 M NaCl, pH 7.8, 0.08% sodium azide.

Storage 4°C

Purity Homogeneous by SDS-PAGE

Other Remarks:

**EXERCISE CAUTIONS IN ITS USE!** It has been found to have moderate toxicity in some animal models, but its toxicity to humans has not been fully determined, and a related protein, RCA<sub>60</sub>, is highly toxic to humans. Never pipet solutions of this lectin by mouth. **DO NOT LYOPHILIZE!**