



---

## CERTIFICATE OF ANALYSIS

Product **RICINUS COMMUNIS AGGLUTININ II (RCA 60, Ricin)**

Catalog No. L-1090

Amount 10 mg

Lot No. ZB1218

Protein Concentration 5.0 mg/ml

Solution 10 mM phosphate, 0.15 M NaCl, pH 7.8, 0.08% sodium azide

Storage 4°C

Toxicity: (50% inhibition)

Cellular toxicity: 3.5ng/ml

NOTE:

**THIS IS A TOXIC LECTIN! EXERCISE CAUTION IN ITS USE!**  
**NEVER PIPET SOLUTIONS OF THIS LECTIN BY MOUTH.**  
**DO NOT LYOPHILIZE!**

**\*\*\* BEFORE USE, READ ATTACHED CAUTION SHEET \*\*\***