

CERTIFICATE OF ANALYSIS

Product GRIFFONIA (BANDEIRAEA) SIMPLICIFOLIA LECTIN I

Catalog No. L-1100

Amount 5 mg of lyophilized lectin

Lot No. ZC0614

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

Although many buffers can be employed for reconstituting and diluting this lectin, we recommend 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca++. For preserving solutions stored at 4 °C, add 0.08% sodium azide to the buffer.

For dissolving the lectin, we recommend adding the lectin slowly to buffer rather than adding the buffer to the lectin. Allow the lectin to dissolve without vigorous agitation.