



CERTIFICATE OF ANALYSIS

Product **SOLANUM TUBEROSUM (POTATO) LECTIN**

Catalog No. L-1160

Amount 5 mg of lyophilized lectin

Lot No. Z0402

Storage Refrigerate; for long term storage, aliquots may be stored frozen or preserved with 0.08% sodium azide in the recommended buffer and stored at 4 °C

Purity Homogeneous by SDS-PAGE

Other remarks:

Essentially salt-free.

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca⁺⁺ is recommended.