

## CERTIFICATE OF ANALYSIS

## Product LYCOPERSICON ESCULENTUM (TOMATO) LECTIN

- Catalog No. L-1170
- Amount <u>2 mg of lyophilized lectin</u>
- Lot No. ZA0916
- Storage Refrigerate; for long term storage, aliquots may be stored frozen or preserved with 0.08% sodium azide and stored at 4 °C
- Purity Homogeneous by SDS-PAGE

Other remarks:

Essentially salt-free.

Although many buffers can be employed for reconstituting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca++ is recommended.