



## CERTIFICATE OF ANALYSIS

Product **DATURA STRAMONIUM LECTIN**  
Catalog No. **L-1180**  
Amount **5 mg of lyophilized lectin**  
Lot No. **ZB1116**  
Storage **Refrigerate; aliquots of reconstituted lectin can be stored frozen or preserved with 0.08% sodium azide and stored at 4 °C**

Purity **Homogeneous by SDS-PAGE**

Other remarks:

Essentially salt-free.

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca<sup>++</sup> is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.

Aggregation may occur with time if stored at concentrations greater than 2 mg/ml.