



CERTIFICATE OF ANALYSIS

Product **GRIFFONIA (BANDEIRAEA) SIMPLICIFOLIA LECTIN II (GSL II, BSL II)**

Catalog No. L-1210

Amount 5 mg of lyophilized lectin*

Lot No. Y0619

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

*In addition to the 5 mg of lectin, the lyophilized powder contains approximately 10% salts.

Reconstitute with 1 ml of purified water.
DO NOT VORTEX.

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca⁺⁺ is recommended. For preserving solutions stored at 4 °C, 0.8% sodium azide can be used.