



CERTIFICATE OF ANALYSIS

Product **VICIA VILLOSA LECTIN**

Catalog No. L-1230

Amount 5 mg of lyophilized lectin

Lot No. ZB0706

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

Essentially salt-free.

Although many buffers can be employed for reconstituting this lectin, we recommend 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca⁺⁺. For preserving solutions stored at 4 °C, 0.8% sodium azide can be used.