



CERTIFICATE OF ANALYSIS

Product **GALANTHUS NIVALIS (SNOWDROP) LECTIN**

Catalog No. L-1240

Amount 5 mg of lyophilized lectin*

Lot No. Y1016

Storage Refrigerate at 4 °C

Purity Homogeneous by SDS-PAGE

Other remarks:

*In addition to the 5 mg of lectin, the lyophilized powder contains approximately 15% salts.

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca⁺⁺ is recommended. Lectin may be preserved with 0.04% sodium azide.