



CERTIFICATE OF ANALYSIS

Product **BAUHINIA PURPUREA LECTIN (BPL)**

Catalog No. L-1280

Amount 5 mg of lyophilized lectin*

Lot No. ZC0318

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

*Essentially salt-free. This lot of BPL contains 18% w/w sugar added for stabilization.

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca⁺⁺ is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.