



CERTIFICATE OF ANALYSIS

Product RICIN B CHAIN

Catalog No. L-1290

Amount 1 mg

Lot No. ZC1215

Protein Concentration 1 mg/ml

Solution 10 mM sodium phosphate, 0.15 M NaCl, pH 7.5
0.08% sodium azide, 50 mM 2-mercaptoethanol

Storage 4 °C DO NOT FREEZE

Activity by solid phase binding assay 99%

NOTE:

Exercise caution when using Ricin B Chain.

DO NOT LYOPHILIZE!

If it is important in the application that the sulfydryl group on the Ricin B Chain is available, add fresh reducing agent to the Ricin B Chain, followed by dialysis or gel filtration at pH 6.5.