



## CERTIFICATE OF ANALYSIS

Product **MAACKIA AMURENSIS LECTIN I (MALI)**

Catalog No. L-1310

Amount 5 mg of lyophilized lectin

Lot No. ZB0612

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

### Other remarks:

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca<sup>++</sup> is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.

High resolution SDS-PAGE pattern shows two subunits of 36,000 and 38,000 daltons.