



CERTIFICATE OF ANALYSIS

Product **LOTUS TETRAGONOLOBUS LECTIN**

Catalog No. L-1320

Amount 5 mg of lyophilized lectin

Lot No. Z0621

Storage Refrigerate at 4 °C

Purity 2 bands by SDS-PAGE corresponding to the subunits of the lectin

Other remarks:

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 7.5 - 8.5, 0.1 mM Ca⁺⁺ is recommended. Lectin may be preserved with 0.04% sodium azide.