



CERTIFICATE OF ANALYSIS

Product **MACLURA POMIFERA LECTIN (MPL)**

Catalog No. L-1340

Amount 5 mg of lyophilized lectin*

Lot No. Z0918

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

*Essentially salt-free. This lot of MPL contains 26% w/w sugar added for stabilization.

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca⁺⁺ is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.