



CERTIFICATE OF ANALYSIS

Product **NARCISSUS PSEUDONARCISSUS (DAFFODIL) LECTIN**

Catalog No. L-1370

Amount 5 mg of lyophilized lectin*

Lot No. Y0620

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

*In addition to the 5 mg of lectin, the lyophilized powder contains approximately 15% salts.

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca⁺⁺ is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.