



## CERTIFICATE OF ANALYSIS

Product **HIPPEASTRUM HYBRID (AMARYLLIS) LECTIN**

Catalog No. L-1380

Amount 5 mg of lyophilized lectin\*

Lot No. ZB1027

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

\*Essentially salt-free. This lot of *Hippeastrum* Hybrid lectin contains 14% w/w sugar added for stabilization.

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca<sup>++</sup> is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.