



CERTIFICATE OF ANALYSIS

Product *ALEURIA AURANTIA LECTIN*

Catalog No. L-1390

Amount 2 mg of lyophilized lectin. In addition to the 2 mg of lectin, the lyophilized powder contains approximately 39% salts.

Lot No. Z0403

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM Ca⁺⁺ is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.

The recommended concentration range for use is 1-10 µg/ml.