



CERTIFICATE OF ANALYSIS

Product **MUSA PARADISIACA (BANANA) LECTIN**
(BanLec) **AFFINITY PURIFIED**

Catalog No. L-1410

Amount 5 mg of lyophilized lectin

Lot No. ZA1120

Storage Refrigerate

Purity Homogeneous by SDS-PAGE

Other remarks:

Reconstitute lectin with 10 mM HEPES buffered saline, pH 7.5, 0.08% sodium azide. Subsequent storage should be at 4 °C.

Inhibiting Sugar: 200 mM α -methylmannoside