



CERTIFICATE OF ANALYSIS

Product **ALKALINE PHOSPHATASE ANTI-DIGOXIGENIN/DIGOXIN (DIG)**
****AFFINITY-PURIFIED****

Catalog No. MB-7100

Amount 300 μ l

Lot No. ZB0928

Produced in goat

Concentration 0.2 mg/ml

Solution phosphate buffered saline, pH 7.4, 0.2 mg/ml BSA,
1 mM MgCl₂, 0.08% sodium azide

Storage 2-8 °C

Recommended dilution range and applications:

1:1000 - 1:2000 for blot applications (Southern, Northern, dot blot)

1:50 - 1:100 for *in situ* hybridization